


ACUMED[®]

THE ACUTRAK[®] FAMILY



A Heritage of Healing

Acutrak[®] Technology
for Headless Compression Screws



Introduction to Acutrak® Technology for Headless Compression Screws

Acumed® is a global leader of innovative orthopaedic and medical solutions.

We are dedicated to developing products, service methods and approaches that improve patient care.




Figure 1: Acutrak® & Acutrak 2® Screws

Acumed® developed the Acutrak® screw technology to provide a headless compression-holding solution to fix fractures using the first fully threaded bone screw with a continuously variable thread pitch along the entire length of the screw.

This white paper will discuss:

- the design rationale behind Acutrak® technology
- the differences between the first and second generation implementations of this technology (Acutrak® and Acutrak 2® respectively)
- examples of indications where each generation has application advantages

This white paper concludes with a comprehensive bibliography of published articles relating to Acutrak® technology biomechanics and clinical results.

Contents

Introducing	2
Design Rationale	2
Fundamental Design Features	2
Crossing the Fracture Site	3
Acutrak® Technology's "Window of Compression"	3
Enhanced Fracture Fixation Biomechanics	4
Acutrak 2 Screws	5
Simplified Surgical Technique	5
Self-Centering Versus Self-Cutting	7
Acutrak® Versus Acutrak 2® Screws	7
Selecting Acutrak® Versus Acutrak 2® Screws	8
Bibliography	9

Design Rationale

1. Minimal soft tissue irritation through Headless Fixation
2. Enhanced fracture fixation through a Fully Threaded Construct
3. Enhance window of compression through a Continuously Variable Screw Pitch
4. Versatility using a Cancellous-Based Thread Design

Incorporating the benefits and features above, Acumed® created a new category of bone screw fixation that goes beyond headed and differential pitch screw options by offering surgeons enhanced biomechanical performance in multiple clinical applications.

Fundamental Design Features: (See in Figure 2)

Acutrak® design features include a unique, patented thread pitch that varies continuously from tip to tail. This ensures each screw rotation engages threads into new bone along the screw's entire length. As each successive individual thread advances faster than the trailing thread counterpart, the conical shape becomes seated into bone. This radial expansion of the screw threads, combined with their axial advancement, creates the ability to reduce & compress bone fragments without a traditional screw head.




Figure 2: The three fundamental characteristics of Acutrak® technology

Crossing the Fracture Site

Conventional wisdom is that threads cannot cross the fracture site and create compression. Acutrak® technology changed conventional orthopaedic wisdom by creating and maintaining compression as the screw threads travel across the fracture site. The Acutrak® screw's continuously variable screw pitch overcame the inherent limitations of the constant pitch and differential pitch found in more traditional bone screws; these traditional approaches and technologies cannot create compression at the fracture site because the screw threads advance at the same rate within their threaded regions. Compression using these devices had to be created prior to screw installation (pre-compression) or through sacrificing thread purchase somewhere along the length of the bone screw. Acutrak® screw technology eliminates this compromise by utilizing continuously variable screw pitch. By allowing each thread along the entire length of the screw to aid in the reduction and compression, the thread location relative to the fracture site no longer became a limitation.




Figure 3: Acutrak® Technology passes threads across a fracture site & creates compression

Acutrak® Technology's Larger "Window of Compression"

All bone screw technologies have a "window of compression" that determines the number of screw rotations needed to reach a maximum compressive force (beyond which further rotations decrease this value). Traditional bone screws have a narrow window of compression. This narrow window characterizes a fixation construct that becomes very sensitive to loss of compression due to over-rotation and the stripping of thread purchase. Acutrak® technology has a wide window of compression, which is less sensitive to stripping the bone and more flexible in its placement for reaching a maximum amount of compression.

Figure 4 below illustrates the window of compression for Acutrak® technology, traditional bone screws (AO) and differential pitch screws (Herbert)




Figure 5: Headed bone screw has the lowest window to reach maximum

Figure 6: Acutrak® screw has the highest window to accomplish maximum compression

Figure 7: Differential pitch falls under both headed and Acutrak® screws in performance

Figure 4: Graph showing the window of screw rotations during which each screw delivers maximum compressive force.

The Acutrak® Screw has the largest window of compression due to the additive property of each variable thread pitch providing compression as a result of being fully threaded.

Introduction to Acutrak® Technology for Headless Compression Screws

Enhanced Fracture Fixation Biomechanics

Pullout strength, resistance to cyclic and torsional loading are the key measurable elements of bone screw fixation performance. Acumed® compared the performance of Acutrak® technology in each of these elements versus traditional (AO) and differential (Herbert) bone screws. A summary of the results is shown in the figures below:

1. Greater Pullout Strength





Figure 8: Acutrak® screws have the highest push-out force when compared to AO & Herbert Bone Screws (AC)

2. Greater Resistance to Cyclic Loading





Figure 9: Acutrak® screws have the highest amount of retained compression after cyclic loading when compared to AO & Herbert Bone Screws (AC)

3. Greater Resistance to Torsional Loading





Figure 10: Acutrak® Screws have the highest resistance to torsional loading when compared to AO & Herbert Bone Screws (AC)

Acutrak 2® Screws (Second Generation Acutrak® Technology)

The original Acutrak® screw families have a long history of success and are used in a wide variety of orthopaedic indications ranging from the tip of the finger to the tip of the toes. In 2005, Acumed® released a second generation of Acutrak® technology in the form of the Acutrak 2® screw families. This second generation of Acutrak® technology was based on the experience gained from the first generation of Acutrak® screw families.

The Acutrak 2® design objectives were:

1. Simplify the surgical technique:
 - a. Reduce the sensitivity to drill depth as related to screw length
 - b. Incorporate self-drilling & self-cutting features
2. Increase the strength of the screw and hex driver interface
3. Increase the stiffness of the guide wires
4. Reduce the radial stress transmitted to the bone without reducing compression or fixation

While retaining the wide “window of compression” and fracture fixation advantages of the original Acutrak® screw design, the Acutrak 2® also minimizes drill depth sensitivity and simplifies the surgical technique.

Simplified Surgical Technique

Drill depth sensitivity quantifies how the Acutrak® technology taper locks into a prepared bone profile. The Acutrak® bone screw taper locks along the entire screw length while the Acutrak 2® screw taper locks along the trailing 1/3 of the screw length. Depth sensitivity emerges when insertion torque rises rapidly upon the taper lock of the screw & prepared bone profile.




Figure 11: Graphs comparing rise in insertion torque at taper lock between Acutrak® & Acutrak 2® Screws

Introduction to Acutrak® Technology for Headless Compression Screws



Figure 12: Fully tapered Acutrak® screw root profile



Figure 13: Partially tapered Acutrak 2® screw root profile



Figure 14: Fully tapered Acutrak® screw profile




Figure 15: Diameter zones for Acutrak 2® screw

Simplified Surgical Technique (Continued)

The “Rules of Acutrak®” were created to manage a rapid increase in insertion torque as an effect of depth sensitivity in the original Acutrak® screws. To obtain an optimal outcome, a precise relationship between drill depth & screw length must be maintained. For example, a measured depth of 55mm, drilled to 55mm, required a screw of length of 50mm to be inserted. This Acutrak® screw rule ensured a smaller taper (downsized length) fit into the prepared bone profile and taper locked when the Acutrak® screw was buried beneath the bone surface. However, the mismatch was counterintuitive and created confusion for surgeons accustomed to the traditional technique of measuring, drilling and inserting to the depth of the matching screw length.

To simplify the surgical technique and enhance insertion flexibility, the Acutrak 2® screw was broken into zones. The outer diameter has three zones - tip, middle and tail - with the tail zone being identical for each screw length. The root diameter has two zones: a cylindrical section with a tapered tail zone, which is also identical for each screw length. The tail zones are approximately 1/3 of the screw length.




Figure 16: Acutrak® screw, conical drill

The Acutrak 2® Screw minimizes depth sensitivity by pairing a cylindrical drill with a tapered profile drill, which eliminates the requirement of “downsizing” as described for the Acutrak® screw

Please refer to the Acutrak® and Acutrak 2® surgical procedures for more specific details on appropriate drilling technique.




Figure 17: Acutrak 2® screw, taper and conical drill pair

Self Centering versus Self-Cutting: Increased Initial Distraction of Distal Bone Fragments

Acutrak® screws are self-centering and Acutrak 2® screws are self-cutting.

The fully tapered profile of an Acutrak® screw initially self-centers with the fully prepared bone profile by crossing over the fracture, fusion or osteotomy site.

The tip of an Acutrak 2® screw must cut into the opposing bony fragment and may distract slightly until the self-cutting tip engages. This momentary distraction is then eliminated as the screw is seated and does not impact the final compression.




Figure 18: Acutrak® screw, self-centering




Figure 19: Acutrak 2® screw, distraction and final compression

Performance differences between Acutrak® and Acutrak 2® screws

Both the Acutrak® and Acutrak 2® bone screw families have been characterized for compression and pullout strength.

Acutrak® and Acutrak 2® products provide equivalent levels of compression.

The Acutrak 2® screws have higher levels of pullout strength as shown to the right.

- The increase in pullout strength is largely due to the similarity of Acutrak 2® tip & tail diameters




Figure 20: Acutrak® screw compression levels




Figure 21: Acutrak® screw pullout strength

Introduction to Acutrak® Technology for Headless Compression Screws




Figure 22: Illustration of a Proximal Pole Fracture in a Scaphoid



Figure 23: 5th Metatarsal Fracture with Acutrak®




Figure 24: Illustration of bone graft for non-union of the scaphoid & a lateral column opening wedge osteotomy

Selecting Acutrak® vs Acutrak 2® bone screws

Acumed believes that both the Acutrak® and Acutrak 2® families of bone screws work well in a broad range of indications. However, surgeons may want to select one versus the other based on the unique aspects of a particular case.

Factors to consider when selecting Acutrak® vs Acutrak 2® are:

1. If significantly less than 50% of the distal bone fragment is present and initial distraction is a concern, Acutrak® may be a better choice. An example is proximal pole fractures of the scaphoid.
2. If fracture fixation requires the use of a screw in the intramedullary canal of a bone, and the anatomy has a conical shape, Acutrak® may be a better choice. Proximal 5th metatarsal fractures are a prime example of this.
3. If fracture fixation requires attempting to join more than two bone fragments with a single screw, Acutrak® screws may be a better choice due to the self-centering feature. Bone grafting for indications such as non-unions of the scaphoid or lateral column opening wedge osteotomies are examples of indications where this situation may exist.

In all other situations, the Acutrak 2® bone screws may provide a more “user friendly” surgical technique because of:

1. Reduced sensitivity to drill profile alignment with screw profile
2. Larger guide wire diameter, resulting in less sensitivity to bone density during guide wire insertion
3. Larger hex drivers, which increase the amount of torque at the driver & screw interface, making screw insertion easier
4. Micro screw option for small fractures, fusions or osteotomies when cannulation is desired & limited bone stock is available

Acutrak® Technology Articles

1. McQueen M.M., Gelbke M.K., Wakefield A., Will E.M., Gaebler A. Percutaneous Screw Fixation Versus Conservative Treatment for Fractures of the Waist of the Scaphoid, *The Journal of Bone and Joint Surgery*, 2007: 66-71.
2. Gholson J.J., Bae D.S., Zurakowski D., Waters P.M. Scaphoid Fractures in Children and Adolescents: Contemporary Injury Patterns and Factors Influencing Time to Union, *The Journal of Bone and Joint Surgery*, 2011: 1210-9.
3. Suh N., Benson E.C., Faber K.J., Macdermid J., Grewal R. Treatment of Acute Scaphoid Fractures: A Systematic Review and Meta- Analysis, *American Association for Hand Surgery*, 2010: 5: 343-353.
4. Richards A.A., Afifi A.M., Moneim M.S. Four-Corner Fusion and Scaphoid Excision Using Headless Compression Screws for SLAC and SNAC Wrist Deformities, *Techniques in Hand and Upper Extremity Surgery*, 2011: 15: 99-103.
5. Kobza P.E., Budoff J.E., Yeh M., Lou Z. Management of the Scaphoid During Four-Corner Fusion - A Cadaveric Study, *The Journal of Hand Surgery*, 2003: 904-9.
6. Ozyurekoglu T., Turker T. Results of a Method of 4-Corner Arthrodesis Using Headless Compression Screws, *Elsevier*, 2012: 486-492.
7. Dutly-Guinand M., von Schroeder H.P. Three-Corner Midcarpal Arthrodesis and Scaphoidectomy, *Techniques in Hand and Upper Extremity Surgery*, 2009: 54-8.
8. Pogliacomi F., Concaro G., Vaienti E. Hahn-Steinthal Fracture: Report of Two Cases, 2005; 76; 178-184.
9. Orr J.D., Glisson R.R., Nunley J.A. Jones Fracture Fixation, *The American Journal of Sports Medicine*, 2012; 40; 691-8.
10. Grewal R., Assini J., Sauder D., Ferreira L., Johnson J., Faber K. A Comparison of Two Headless Compression Screws for Operative Treatment of Scaphoid Fractures, *Journal of Orthopaedic Surgery and Research*, 2011; 6; 27-33.
11. Kilic A., Sokucu S., Parmaksizoglu A.S., Glu M., Kabukcuoglu Y.S. Comparative Evaluation of Radiographic and Functional Outcomes in the Surgical Treatment of Scaphoid Non-Unions, *Turkish Association of Orthopaedics and Traumatology*, 2011: 399-405.
12. Panchal A., Kubiak E.N., Keshner M., Fulkerson E., Paksima N. Comparison of Fixation Methods for Scaphoid Nonunions, *Bulletin of the NYU Hospital for Joint Disease*, 2007; 65(4); 271-5.
13. Morin P., Reindl R., Berry G.K., Harvey E.J. Incorrect Radiographic Evaluation of Vascularized Bone Grafting for Scaphoid Fracture or Nonunion, *Canadian Journal of Plastic Surgery*, 2011: 7-9.
14. Meermans G., Verstreken F. Influence of Screw Design, Sex, and the Approach in Scaphoid Fracture Fixation, *Clinical Orthopaedics and Related Research*, 2011.
15. Beadel G.P., Ferreira L., Johnson J.A., King G.J.W. Interfragmentary Compression Across a Simulated Scaphoid Fracture – Analysis of 3 Screws, *The Journal of Hand Surgery*, 2004: 273-8.
16. Crawford L.A., Powell E.S., Trail I.A., The Fixation Strength of Scaphoid Bone Screws: An In Vitro Investigation Using Polyurethane Foam, *Elsevier*, 2012: 255-260.
17. Amadio P.C. What's New in Hand Surgery, *The Journal of Bone and Joint Surgery*, 2003; 85; 389-393.
18. Bond C.D., Shin A.Y., McBride M.T., Doa K.D. Percutaneous Screw Fixation or Cast Immobilization for Nondisplaced Scaphoid Fractures, *The Journal of Bone and Joint Surgery*, 2001; 83; 483-8.
19. Brutus J-P, Palmer A.K., Mosher J.F., Harley B.J., Loftus J.B. Use of a Headless Compressive Screw for Distal Interphalangeal Joint Arthrodesis in Digits: Clinical Outcome and Review of Complications, *The Journal of Hand Surgery*, 2008; 31; 85-9.

Bibliography

20. Gruska D.S., Burkhart K.J., Nowak T.E., Achenbach T., Rommens P.M., Müller L.P. The Durability of the Intrascaphoid Compression of Headless Compression Screws: In Vitro Study, *The Journal of Hand Surgery*, 2012; 37; 1142-1150.
21. Lee S.K., Kubiak E.N., Liporace F.A., Parisi D.M., Iesaka K, Posner M.A. Fixation of Tendon Grafts for Collateral Ligament Reconstructions: A Cadaveric Biomechanical Study, *The Journal of Hand Surgery*, 2005; 30; 1051-5.
22. Jeon I-H, Micic I.D., Oh C-W, Park B-C, Kim P-T. Percutaneous Screw Fixation for Scaphoid Fracture: A Comparison Between the Dorsal and Volar Approaches, *The Journal of Hand Surgery*, 2009; 34; 228-236.
23. Bedi A, Jebson P.J.L., Hayden R.J., Jacobson J.A., Martus J.E. Internal Fixation of Acute, Nondisplaced Scaphoid Waist Fractures Via a Limited Dorsal Approach: An Assessment of Radiographic and Functional Outcomes, *The Journal of Hand Surgery*, 2007; 32; 326.e1-326.e9.
24. Ibrahim T, Qureshi A, Sutton A.J., Dias J.J. Surgical Versus Nonsurgical Treatment of Acute Minimally Displaced and Undisplaced Scaphoid Waist Fractures: Pairwise and Network Meta-Analyses of Randomized Controlled Trails, *The Journal of Hand Surgery*, 2011; 36; 1759-1768.
25. Budoff J.E., Meyers N, Ambrose C.G. The Comparative Stability of Screw Versus Plate Versus Screw and Plate Coronoid Fixation, *The Journal of Hand Surgery*, 2011; 36; 238-245.
26. Jarrett P, Kinzel V, Stoffel K. A Biomechanical Comparison of Scaphoid Fixation With Bone Grafting Using Iliac Bone or Distal Radius Bone, *The Journal of Hand Surgery*, 2007; 32; 1367-1373.
27. Dodds S.D., Panjabi M.M., Slade III J.F. Screw Fixation of Scaphoid Fractures: A biomechanical Assessment of Screw Length and Screw Augmentation, *The Journal of Hand Surgery*, 2006; 31; 405-413.
28. Walsh E, Crisco J.J., Wolfe S.W. Computer-Assisted Navigation of Volar Percutaneous Scaphoid Placement, *The Journal of Hand Surgery*, 2009; 34; 1722-8.
29. Leung Y.F., Ip S.P.S., Cheuk C, Sheung K.T., Wai Y.L. Trepine Bond Grafting Technique for the Treatment of Scaphoid Nonunion, *The Journal of Hand Surgery*, 2001; 26; 893-900.
30. Chan K.W., McAdams T.R. Central Screw Placement in Percutaneous Scaphoid Screw Fixation: A Cadaveric Comparison of Proximal and Distal Techniques, *The Journal of Hand Surgery*, 2004; 29; 74-9.
31. Villani F, Uribe-Echevarria B, Vaienti L, Distal Interphalangeal Joint Arthrodesis for Degenerative Osteoarthritis With Compression Screw: Results in 102 Digits, *The Journal of Hand Surgery*, 2012; 37; 1330-4.
32. Yassae F, Yang S. Mini-Incision Fixation of Nondisplaced Scaphoid Fracture Nonunions, *The Journal of Hand Surgery*, 2008; 33; 1116-1120.
33. Adamany D.C., Mikola E.A., Fraser B.J. Percutaneous Fixation of the Scaphoid Through a Dorsal Approach: An Anatomic Study, *The Journal of Hand Surgery*, 2008; 33; 327-331.
34. Knoll V.D., Allan C, Trumble T.E. Trans-Scaphoid Perilunate Fracture Dislocations: Results of Screw Fixation of the Scaphoid and Lunotriquetral Repair With a Dorsal Approach, *The Journal of Hand Surgery*, 2005; 30; 1145e.1-1145e.11.
35. Nalbantoglu U, Gereli A, Kocaoğlu B, Aktas S, Turkmen M. Capitellar Cartilage Injuries Concomitant With Radial Head Fractures, *The Journal of Hand Surgery*, 2008; 33; 1602-7.
36. Hankins C.L., Budoff J.E. Analysis of Wrist Motion Following Vascularized Bone Graft to the Proximal Scaphoid, *The Journal of Hand Surgery*, 2011; 36; 583-6.
37. Zlotolow D.A., Knutsen E, Yao J. Optimization of Volar Percutaneous Screw Fixation for Scaphoid Waist Fractures Using Traction, Positioning, Imaging, and an Angiocatheter Guide, *The Journal of Hand Surgery*, 2011; 36; 916-921.

38. Slutsky D.J., Nagle D.J. Wrist Arthroscopy: Current Concepts, *The Journal of Hand Surgery*, 2008; 33; 1228-1244.
39. Wong A.S., Baratz M.E. Elbow Fractures: Distal Humerous, *The Journal of Hand Surgery*, 2009; 34; 176-190.
40. Waitayawinyu T, McCallister W.V., Katolik L.I., Schlenker J.D., Trumble T.E. Outcome After Vascularized Bone Grafting of Scaphoid Nonunions With Avascular Necrosis, *The Journal of Hand Surgery*, 2009; 34; 387-394.
41. Leventhal E.L., Wolfe S.W., Walsh E.F., Crisco J.J. A Computational Approach to the "Optimal" Screw Axis Location and Orientation in the Scaphoid Bone, *The Journal of Hand Surgery*, 2009; 34; 677-684.
42. Wang M.L., Bednar J.M. Lunatocapitate and Triquetrohamate Arthrodeses for Degenerative Arthritis of the Wrist, *The Journal of Hand Surgery*, 2012; 37; 1136-1141.
43. Martineau P.A., Bergeron S, Beckman L, Steffen T, Harvey E.J. Reconstructive Procedure for Unstable Radial-Sided Triangular Fibrocartilage Complex Avulsions, *The Journal of Hand Surgery*, 2005; 30; 727-732.
44. Sekiya J.K., Jolly J, Debski R.E. The Effect of a Hill-Sachs Defect on Glenohumeral Translations , In Situ Capsular Forces, and Bony Contact Forces, *The American Journal of Sports Medicine*, 2012; 40; 388-394.
45. Rettig A.C. Athletic Injuries of the Wrist and Hand, *The American Journal of Sports Medicine*, 2003, 31; 1038-1048.
46. Oduwole K.O., Cichy B, Dillon J.P., Wilson J, O'Beirne J. Acutrak Versus Herbert Screw Fixation for Scaphoid Nonunion and Delayed Union, *The Journal of Orthopaedic Surgery*, 2012; 20(1):61-5.
47. Unal A.M., Baran O, Uzun B, Turan A.C., Comparison of Screw-Fixation Stabilities of First Metatarsal Shaft Osteotomies: A Biomechanical Study, *Turkish Association of Orthopaedics and Traumatology*, 2010;44(1): 70-75.
48. Ahmad Z, McGuinness C.N. Thumbs Up! A Novel Use of the Acutrak Screw Fixation System for the Management of Triphalangeal Thumb, *Journal of Plastic Surgery*, 2010;10: 217-222.
49. Lee YS, Hsu TL, Huang CR. Lateral Fixation of AO Type-B2 Ankle Fractures: The Acutrak Plus Compression Screw Technique, *International Orthopaedics*, 2010;34: 903-7.
50. Weisler E.R., Chloros G.D., Kuzma G.R. Arthroscopy in the Treatment of Fracture of the Trapezium, *The Journal of Arthroscopic and Related Surgery*, 2007;23(11): 1248e.1-1248e.4.
51. Cekin T, Tukenmez M, Tezeren G, Comparison of Three Fixation Methods in Transverse Fractures of the Patella in a Calf Model, *Acta Orthop Traumatol Turc*, 2006;40(3): 248-51.
52. Adolfsson L, Lindau T, Arner M. Acutrak Screw Fixation Versus Cast Immobilisation for Undisplaced Scaphoid Waist Fractures, *The Journal of Hand Surgery (European)*, 2001; 26(3): 192-5.
53. Gregory J.J., Mohil R.S., Ng A.B., Warner J.G., Hodgson S.P. Comparison of Herbert and Acutrak Screws in the Treatment of Scaphoid Non-Union and Delayed Union, *Acta Orthopaedica Belgica*, 2008;74(6): 761-5.
54. Moholkar K, O'Sullivan T. Acutrak Fixation of Comminuted Distal Radial Fractures, *Acta Orthopaedica Belgica*, 2004;70: 478-481.
55. Ganapathi M, Joseph G, Savage R, Jones A.R., Timms B, Lyons K. MRI Susceptibility Artefacts Related to Scaphoid Screws: The Effect of Screw Type, Screw Orientation and Imaging Parameters, *Journal of Hand Surgery (Br.)*, 2002;27: 165-170.
56. Loving V.A., Richardson M.L. Scaphoid Fracture Fixation with an Acutrak Screw, *Radiology Case Reports*, 2006;1(2): 58-60.

Bibliography


57. Fadel G.E., Hussain S.M., Sripatha S, Jain A.S. Fixation of First Metatarsal Basal Osteotomy Using Acutrak Screw, *Journal of Foot and Ankle Surgery*, 2008;14: 21-25.
58. Martus J.E., Bedi A, Jebson P.J.L. Cannulated Variable Pitch Compression Screw Fixation of Scaphoid Fractures Using a Limited Dorsal Approach, *Techniques in Hand and Upper Extremity Surgery*, 2005;9(4): 202-6.
59. Yip T.H., Wu W.C. Symposium on Advances in the Management of Scaphoid Problems, *Hong Kong Journal of Orthopaedic Surgery*, 2002;6(2): 76-81.
60. Assari S, Darvish K, Ilyas A.M. Comparison of Five Headless Screws for Fixation of Small Bones.
61. Lenehan B, Street J, Murphy B, Brennan L, McHugh P, Curtin W. A Biomechanical Study Comparing the Compressive Forces Generated by a Conventional 4.5 AO/ASIF Cortical Lag Screw with a Differentially Pitched Cortical Compression Screw, *Acta of Bioengineering and Biomechanics*, 2009;11(1): 31-5.
62. Geissler W.B. Bone Graft Substitutes in the Upper Extremity, *The Hand Clinics*, 2006;22: 329-339.
63. Geissler W.B. Arthroscopic Management of Scaphoid Fractures in Athletes, *The Hand Clinics*, 2009;25: 359-369.
64. Geissler W.B. Operative Fixation of Metacarpal and Phalangeal Fractures in Athletes, *The Hand Clinics*, 2009;25: 409-421.
65. Geissler W.B. Cannulated Percutaneous Fixation of Intra-articular Hand Fractures, *The Hand Clinics*, 2006;22: 297-305.
66. Geissler W.B. Arthroscopic Assisted Fixation of Fractures of the Scaphoid, *The Hand Clinics*, 2003;8: 37-56.
67. Slade III J.F., Moore A.E. Dorsal Percutaneous Fixation of Stable, Unstable, and Displaced Scaphoid Fractures and Selected Nonunions, *Atlas of the Hand Clinics*, 2003;8: 1-18.
68. Slade III J.F., Moore A.E. Percutaneous Treatment of Transscaphoid, Transcapitate Fracture-dislocations with Arthroscopic Assistance, *Atlas of the Hand Clinics*, 2003;8: 77-94.
69. Slade III J.F., Bombardieri D.A. Percutaneous Capitate Arthrodesis Using Arthroscopic or Limited Approach, *Atlas of the Hand Clinics*, 2003;8: 149-162.
70. Slade III J.F., Taksali S, Safanda J. Combined Fractures of the Scaphoid and Distal Radius: A Revised Treatment Rationale Using Percutaneous and Arthroscopic Techniques, *Atlas of the Hand Clinics*, 2005;21: 427-441.
71. Elkowitz S.J., Polatsch D.B., Egol K.A., Kummer F.J., Koval K.J. Capitate Fractures: A Biomechanical Evaluation of Three Fixation Methods, *Journal of Orthopaedic Trauma*, 2002;16(7): 503-6.
72. Elkowitz S.J., Kubiak E.N., Polatsch D, Cooper J, Kummer F.J., Koval K.J. Comparison of Two Headless Screw Designs for Fixation of Capitate Fractures, *Bulletin – Hospital for Joint Disease*, 2003;61(3-4): 123-6.
73. Kocak E, Carruthers K.H., Kobus R.J. Distal Interphalangeal Joint Arthrodesis with the Herbert Headless Compression Screw: Outcomes and Complications in 64 Consecutively Treated Joints, *American Association for Hand Surgery*, 2011;6: 56-9.
74. Leibovic S.J. Arthrodesis of the Interphalangeal Joints With Headless Compression Screws, *The Journal of Hand Surgery*, 2007;32A: 1113-9.
75. Gutow A.P. Percutaneous Fixation of Scaphoid Fractures, *Journal of the American Academy of Orthopaedic Surgeons*, 2007;15(8): 474-485.
76. Carro L.P., Golano P, Vega J. Arthroscopic Subtalar Arthrodesis: The Posterior Approach in the Prone Position, *The Journal of Arthroscopic and Related Surgery*, 2007;23(4): 445.e1-445.e4.

77. Bain G.I., McLean J.M., Turner P.C., Sood A, Pourgiezis N. Translunate Fracture With Associated Perilunate Injury: 3 Case Reports With Introduction of the Translunate Arc Concept, *The Journal of Hand Surgery*, 2008;33(10): 1770-6.
78. Cobb T.K. Arthroscopic Distal Interphalangeal Joint Arthrodesis, *Techniques in Hand & Upper Extremity Surgery*, 2008;12(4): 266-9.
79. Park M.J., Ahn J.H. Arthroscopically Assisted Reduction and Percutaneous Fixation of Dorsal Perilunate Dislocations and Fracture-Dislocations, *The Journal of Arthroscopic and Related Surgery*, 2005;21(9): 1153.e1-1153.e9.
80. Slade J.F., Lozano-Calderon S, Merrell G, Ring D. Arthroscopic-Assisted Percutaneous Reduction and Screw Fixation of Displaced Scaphoid Fractures, *The Journal of Hand Surgery (European Volume)*, 2008;33(3): 350-4.
81. Ruchelsman D.E., Tejwani N.C., Kwon Y.W., Egol K.A. Coronal Plane Partial Articular Fractures of the Distal Humerus: Current Concepts in Management, *Journal of the American Academy of Orthopaedic Surgeons*, 2008;16(12): 716-728.
82. Kohanzadeh S, Gaon M, Ota K, Lefrennierre S, Kulber D. Use of a Headless Compressive Screw (Acutrak) for Radioulnar Joint Fusion in the Sauve-Kapandji Procedure, *The Journal of the American Society of Plastic Surgeons*, 2012;129(4): 759e-761e.
83. Kang S-H, Kim H-M, Jeong C. Modified Volar Percutaneous Screw Fixation for the Scaphoid Fractures: Trans-trapezial Approach, *Journal of the Korean Society for Surgery of the Hand*, 2010;15(4): 175-183. *Article not in English
84. Singiseti K, Aldlyami E, Middleton A. Early Results of a New Implant: 3.0 mm Headless Compression Screw for Scaphoid Fracture Fixation, *Journal of Hand Surgery (European Volume)*, 2012.
85. Johnson D.L., Durbin T.C. Physeal-sparing Tibial Eminence Fracture Fixation With a Headless Compression Screw, *Orthopedics*, 2012;35(7): 604-8.
86. Kim J.K., Kim J.O., Lee S.Y. Volar Percutaneous Screw Fixation for Scaphoid Waist Delayed Union, *Clinical Orthopaedics and Related Research*, 2010;468: 1066-1071.
87. Shin A.Y. Four-Corner Arthrodesis, *Journal of the American Society for Surgery of the Hand*, 2001;1(2): 93-111.
88. De Paula E.J.L., Mattar Jr R, Rezende M.R., Franca E.N. Complex Wrist Injuries: Fractures of the Lunate, *The Journal of Hand Surgery (European Volume)*, 2011;36E(1): S66-S73.
89. Kim J.P., Lee J.S., Park M.J. Arthroscopic Reduction and Percutaneous Fixation of Perilunate Dislocations and Fracture-Dislocations, *The Journal of Arthroscopic and Related Surgery*, 2012;28(2): 196-203.
90. Kilic A, Kabukcuoglu Y, Sokucu S. The Treatment of Talar Body Fractures with Compression Screws: A Case Series, *Cases Journal*, 2009;2: 7953.
91. Mantovani G, Fukushima W.Y., Cho A.B., Aita M.A., Lino Jr W., Faria F.N. Alternative to the Distal Interphalangeal Joint Arthrodesis: Lateral Approach and Plate Fixation, *The Journal of Hand Surgery*, 2008;33(1): 31-34.
92. Dutly-Guinand M, P. von Schroeder H. Three-Corner Midcarpal Arthrodesis and Scaphoidectomy, *Techniques in Hand & Upper Extremity Surgery*, 2009;13(1): 54-8.
93. Bendre A.A., Hartigan B.J., Kalainov D.M. Mallet Finger, *Journal of the American Academy of Orthopaedic Surgeons*, 2005;13(5): 336-344.
94. Sommerkamp T.G. Scaphoid AARIF: Arthroscopically Assisted Reduction and Internal Fixation of Scaphoid Fractures, *Journal of the American Society for Surgery of the Hand*, 2001;1(3): 192-210.
95. Chung K.C. A simplified approach for unstable scaphoid fracture fixation using the Acutrak screw. *Plastic and Reconstructive Surgery* 2002;110(7): 1697-1703.

Bibliography

96. Newport M.L., Williams C.D., Bradley W.D. Mechanical strength of scaphoid fixation. *Journal of Hand Surgery (Br.)* 1996;21B(1): 99-102.
97. Sides,S.,Fetter,N.,Glisson,R.,Nunley,J. Bending Stiffness and Pull Out Strength of Tapered, Variable Pitch Screws, and 6.5 mm. Cancellous Screws in Acute Jones Fractures, *The American Orthopaedic Foot and Ankle Society*,2006: 821-825.
98. Collman,D., Kaas,M.,Schuberth,J. Arthroscopic Ankle Arthrodesis: Factors Influencing Union in 39 Consecutive Patients, *The American Orthopaedic Foot and Ankle Society*,2006:1079-1085.
99. Huber,W.,Reishner,R.,Trieb,K.,Wanivenhaus,A.,Beer,R. In-Vitro Comparison of Biomechanical Efficiency of Three Cannulated Screws for Arthrodesis of the Hindfoot, *The American Orthopaedic Foot and Ankle Society*,2008: 225-230.
100. Wheeler DL, McLoughlin Sw Biomechanical assessment of compression screws *Clinical Orthopaedics and Related Research*, 1998 May;(350):237-45
101. Chuckpaiwong,B.,Easley,M.,Glisson,R. Screw Placement in Subtalar Arthrodesis: A Biomechanical Study, *The American Orthopaedic Foot and Ankle Society*,2009: 133-141.
102. Elkowitz,S.,Polatsch,D.,Egol,K.,Kummer,F.,Koval,K.Capitellum Fractures: A Biomechanical Evaluation of Three Fixation Methods, *Journal of Orthopaedic Trauma*2002,16:7;503-506.
103. Nguyen,C.,Singh,D.,Harrison,M.,Blunn,G.,Herman,A.,Dudkiewicz,I. Biomechanical Properties of New Mini Compression Screws, *Foot and Ankle International*,2009:545-550.
104. Huber,W.,Reishner,R.,Trieb,K.,Wanivenhaus,A.,Beer,R. In-Vitro Comparison of Biomechanical Efficiency of Three Cannulated Screws for Arthrodesis of the Hindfoot, *The American Orthopaedic Foot and Ankle Society*,2008: 225-230.
105. Mighell,M.,Harkins,D.,Klein,D.,Schneider,S.,Frankle,M. Technique for Internal Fixation of Capitellum and Lateral Trochlea Fractures, *Journal of Orthopaedic Trauma*,2006;20:10,699-704.
106. Pirela-Cruz,M.,Battista,V.,Burnette,S.,Hansen,T.A Technical Note on Percutaneous Scaphoid Fixation Using a Hybrid Technique, *Journal of Orthopaedic Trauma*, 2005;19:8.570-573.
107. Haddad F.S., Goddard N.J. Acute Percutaneous Scaphoid fixation using a Cannulated screw. *Annals of Hand and Upper Limb Surgery (Fr.-Annales de Chirurgie de la Main et du Membre Superieur)* 1998;17(2): 119-126.
108. Caterini,R.,Farsetti,P.,Tarantino,U.,Potenza,V.,Ippolito,E. Arthrodesis of the Toe Joints with an Intramedullary Cannulated Srew For Correction of Hammertoe Deformity,2004: 256-261.
109. Toby,E.,Butler,T.,McCormack,T.,Jayaraman,G. A Comparison of Fixation Screws for the Scaphoid during Application of Cyclical Bending Loads, *The Journal of Bone and Joint Surgery*,1997; 79:1190-7.
110. Ring,D.,Quintero,J.,Jupiter,J. Open Reduction and Internal Fixation, of Fractures of the Radial Head, *The Journal of Bone and Joint Surgery*,2009;84:1811-1815.
111. Slade,J.,Gutow,A.,Geissler,W. Percutaneous Internal Fixation of Scaphoid Fractures via an Arthroscopically Assisted Dorsal Approach, *The Journal of Bone and Joint Surgery*,2009;84:31-36
112. Slade,J.,Gutow,A.,Geissler,W.,Merrel,G. Percutaneous Internal Fixation of Selected Scaphoid Nonunions with an Arthroscopically Assisted Dorsal approach, *The Journal of Bone and Joint Surgery*, 2009;84:21-36.
113. Dubberley,J.,Faber,k.,MacDermid,J.,Patterson,S.,King,G.Outcome after Open Reduction and Internal Fixation of a Capitellar and Trochlear Fractures, *The Journal of Bone and Joint Surgery*, 2009;88:46-54.
114. Haisman,J.,Rohde,R.,Weiland,A.,Acute Fractures of the Scaphoid, *The Journal of Bone and Joint Surgery*,2009;88:2750-2758.

115. Ruschelsman,D.,Tejwani,N.,Kwon,Y.,Egol,K. Open Reduction and Internal Fixation of Capitellar Fractures with Headless Screws, *The Journal of Bone and Joint Surgery*, 2009;90:1321-1329.
116. Cole,B.,Pascual-Garrido,C.,Grumet,R.Surgical Management of Articular Cartilage Defects in the Knee,*The Journal of Bone and Joint Surgery*,2009;91:1778-1790.
117. Bhandari,M.,Hanson,B. Acute Nondisplaced Fractures of the Scaphoid, *Journal of Orthopaedic Trauma*,2004;18:4,253-255.
118. Sides,S.,Fetter,N.,Glisson,R.,Nunley,J. Bending Stiffness and Pull-Out Strength of Tapered, Variable Pitch Screws, and 6.5-mm Cancellous Screws in Acute Jones Fractures, *Foot and Ankle International*,2006:821-825.
119. Haddad F.S., Goddard N.J. Acute percutaneous scaphoid fixation. *The Journal of Bone and Joint Surgery (Br.)* 1998;80-B(1): 95-99. Slide presentation attached. London, England: Royal Free Hospital, Department of Orthopaedics. (Presented as a Scientific Exhibit at the 66th Annual Academy of Orthopaedic Surgeons, Anaheim, CA, February 1999.) (A good clinical paper on the Acutrak. Described is a technique and results of acute scaphoid fixation.)
120. Haddad F.S., Goddard N.J. Acute percutaneous scaphoid fixation using a cannulated screw. *Annals of Hand and Upper Limb Surgery (Fr.-Annales de Chirurgie de la Main et du Membre Superieur)* 1998;17(2): 119-126.



ACUMED®

5885 NW Cornelius Pass Road
Hillsboro, OR 97124
(888) 627-9957
www.acumed.net

Distributed by:

SPF60-00-A
Effective: 3/2013
© 2013 Acumed® LLC

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located. Specific questions physicians may have about the availability and use of the products described on these materials should be directed to their particular local sales representative. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.