OsteoMed® ExtremiLock

Ankle Plating System

EN	Product Information and Instructions for Use	2
FR	Informations sur le produit et instructions d'utilisation	6
DE	Produktinformationen und Gebrauchsanweisung.....	10
IT	Informazioni sul prodotto e istruzioni per l'uso.....	15
PT	Informações sobre o produto e instruções de utilização	19
ES	Información sobre el producto e instrucciones de uso	23

Product Information and Instructions for Use

OsteoMed® ExtremiLock Ankle Plating System

EN These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals in the European Union only.

Description

The **OsteoMed ExtremiLock Ankle Plating System** provides a comprehensive solution for ankle fractures and ankle fusion management. The plates are provided in a variety of shapes and sizes, offering surgeons compression and locking hole designs. The **OsteoMed ExtremiLock Ankle Plating System** includes angulated locking screws and standard non-locking screws as well as cannulated screws. The OsteoMed ExtremiLock Ankle Plating System plates, screws, washers and instrumentation are contained in a single instrument tray. The system instrumentation is used to facilitate plate placement, screw insertion and/or removal of implants.

Material

The plates are made of titanium (ASTM F67) or titanium-alloy (ASTM F136). The screws are made of titanium-alloy (ASTM F136). The instrumentation is made from various grades of stainless steel, anodized aluminum, and/or medical grade polymers. The material composition of Acumed implants is described in detail in the available Metal Sensitivity Statement document "Metal Sensitivity Statement" which can be found at <http://www.acumed.net/ifu>.

Intended Purpose & Indications for use

The **ExtremiLock Ankle Plating System** is intended for fixation of fractures, osteotomies, and non-unions of the tibia and fibula as well as arthrodesis of the ankle. The implants are intended for single use only. The 1/3 tubular plates, hook plates, screws, and washers are also intended for use in trauma, general surgery, and reconstructive procedures of bones appropriate for the size of the device. The intended clinical benefit of **OsteoMed ExtremiLock Ankle Plating System** is to provide fixation and bone fusion for patients with fractures as a result of trauma, Osteotomy or reconstruction of the bones appropriate for the size of the device.

The indirect clinical benefit of the encompassed instrumentation is to facilitate implantation protocols for fracture fixation.

Contraindications

- Active or latent infection or marked inflammation of the treatment area
- Sepsis
- Insufficient quantity or quality of bone, osteoporosis, or in patients with certain metabolic diseases
- Patients with confirmed material sensitivity
- Patients who are unwilling or incapable of following post-operative care instructions and/or the limitations of internal rigid fixation implants
- Percutaneous K-wire placement is contraindicated in cases of displaced fractures and compressed fractures

Warnings

- Treatment or implant may fail, including sudden failure, because of:
 - Loose fixation and/or loosening
 - Stress, include stress from inappropriate bending of the implant during surgery
 - Stress concentrations
 - Stress of weight bearing, load bearing, or excessive activity
- Failure is more likely if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing. Failure is more likely if the patient does not follow post-operative care instructions.
- Nerve or soft tissue damage may result from surgical trauma or the presence of an implant.
- Instrument breakage or damage, as well as tissue damage, may occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use, or unintended use.
- Devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects.
- Fractured implants should be removed from patients during surgery. If unable to remove, notify the patient.
- Implants may cause distortion and/or block the view of anatomic structures on radiographic images.
- In patients with a large intramedullary canal, the diameter/length of the Fusion screw provided may not provide adequate compression of the MCP joint.

Caution

- Implants and instruments are intended only for professional use by a licensed physician.
- Use devices and instruments in accordance with the Instructions for Use.
- Do not use the sterile product past the use-by date. Refer to the device label.
- Do not use or re-sterilize an implant provided in sterile packaging if the package has been damaged. The sterility may be compromised, and the cleanliness of the implant may be uncertain. Report damaged packaging to your distributor or Acumed.
- Use of this system has not been evaluated in children or individuals that are not skeletally mature.
- Inspect all components preoperatively to ensure utility. Do not attempt a surgical procedure with faulty, damaged, or suspect instruments or implants. Alternate fixation methods should be available intraoperatively.
- Reprocessing and/or reuse of single use devices may result in infection/cross contamination, and/or sudden failure due to previous stresses.
- Use of implant components from different manufacturers has not been evaluated.
- Use of chemical disinfection may leave residues that adversely affect steam sterilization.
- Steam penetration and device sterilization may be negatively affected by labeling that block tray steam holes.
- Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling biohazardous material and disposal of sharp materials.

Potential Adverse Effects

- Anesthesia-related problems, problems with positioning of the patient (e.g. nausea, vomiting, neurological impairments, etc.) thrombosis, embolism, infection or injury to other critical structures such as blood vessels, excessive bleeding, etc.
- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.
- Swelling, abnormal scarring, impairment in musculoskeletal function, hardware prominence, migration, loosening, bending, or breakage of the implant, malunion, nonunion or delayed union due to prolonged loading or excessive forces which may lead to implant failure and reoperation.
- Metal sensitivity, histological, allergic, or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu>.
- Injury to user.

Important: Any serious incident that has occurred in relation to the device should be reported to Acumed (through customercomplaints@acumed.net or +1.888.627.9957) and the Competent Authority of the Member State in which the user and/or patient is established.

Target Population and Intended Users

The OsteoMed ExtremiLock Ankle Plating System is intended for use by suitably trained and qualified surgeons in a hospital operating room setting to treat skeletally mature, adult patients with sufficient quantity and quality of bone who require surgical intervention to achieve fracture fixation, fusion or reconstruction of the tibia or fibula bones.

MRI Status

Many Acumed implants have been evaluated for safety in the MR environment. Consult our publication "Acumed Implants in the MR Environment" at <http://www.acumed.net/ifu> for more information.

Device Lifetime

Once installed, implants are expected to provide fixation and physiological support and have an effective life during bone healing. The implants are biocompatible and may remain implanted at the discretion of the surgeon or patient. Multiple-use instruments have a lifetime that is affected by usage, handling, and processing. Assess multiple-use instruments for fitness during pre-sterilization inspection.

Summary of Safety and Clinical Performance (SSCP)

The SSCP for the implants may be obtained from the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Use the following search terms: 08456940BUDIO21NQ.

Maintaining Device Effectiveness

- The surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation products and techniques.
- The surgeon must exercise reasonable judgment when deciding which implant or instrument type to use for specific indication.
- Devices should be used in an operating room environment.
- All Acumed plates, screws, and instrumentation may be required for each surgery. Failure to use dedicated, unique Acumed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure that could cause subsequent patient injury. Failed devices may require re-operation and removal.
- Inspect the Acumed implants prior to use. Do not attempt a surgical procedure with faulty, damaged, or suspected instruments or implants. Alternate fixation methods should be available intraoperatively.
- Instrument lifetime is affected by usage, handling, and processing. Assess multiple-use instruments for fitness before and after each procedure to confirm that they are in proper operating condition. Indications that a reusable instrument should not continue to be used include but are not limited to deformation, dulled sharps, cracks, breakage, missing components, corrosion, excessive wear, and restriction or seizure of moving parts. Instruments which are faulty, damaged, or suspected should not be used. They should be replaced or sent to Acumed for disposition and repair.
- When placing more than one screw, ensure that subsequent screw placement does not interfere with other screws. Insert the second screw on the opposite side of the fracture or osteotomy site, and then all remaining screws, following the outlined procedures.

Instructions for Use

- Acumed offers one or more Surgical Techniques to promote the safe and effective use of this system. Consult our Surgical Techniques at www.acumed.net. Refer to the appropriate OsteoMed ExtremiLock Ankle Plating System Surgical Technique Guide (STG) for detailed guidelines on proper screw and plate placement.
- **Important:** Surgical techniques may contain important safety information. The instruments and implants in this system are intended to be used by suitably trained and qualified surgeons in a hospital operating room setting. Before treatment, the surgeon is advised to read and fully understand all instructions and communicate to the patient any relevant medical information provided therein, including the use, limitations, risks (safety communications), and possible adverse effects of the proposed treatment. Consult the most recent version of the instructions for Use and Surgical Technique Guides as they are subject to change. Contact Acumed or an authorized agent to request any additional information.

Cleaning

1. Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
2. Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
3. Acumed recommends the following cleaning and sterilization instructions for re-usable Instrumentation.

Manual Cleaning

IMPORTANT: All non-sterile instruments and implants are to be cleaned at the healthcare facility prior to use. Implants are subjected to automated cleaning only.

4. Rinse the articles to be cleaned under running cool tap water to remove visible soil until visibly clean.
5. Prepare an enzymatic cleaner, Enzol[®], or equivalent, per manufacturer's recommendations. Fully immerse the articles in the solution and soak for a minimum of 2 minutes. Actuate the articles while immersed in the solution to ensure complete penetration of cleaning solution.
6. Using a soft bristled brush, clean the entire article paying close attention to hard to reach areas until all evidence of soil is removed. A syringe may be used to clean the lumens and other hard to reach areas. Actuate the article while brushing in order to clean mated surfaces and movable parts.
7. Prepare a mild detergent such as Prolystica[®] 2X Concentrate Neutral, or equivalent, per manufacturer's recommendations. Fully immerse the articles in the prepared solution and sonicate the articles for a minimum of 10 minutes.
8. Remove the articles from the detergent and rinse under running reverse osmosis / deionized (RO/DI) water until all evidence of detergent is removed.
9. Dry using a clean, soft cloth and filtered pressurized air at ≤40 psi. Visually examine each article for visible soil. Repeat cleaning procedure if necessary.
10. Steam Autoclave per the recommended sterilization Instructions.

Automated Cleaning

Contaminated instruments require manual cleaning prior to automated cleaning. Follow the washer-disinfector manufacturer's instructions explicitly.

11. Prepare a cleaning neutral enzymatic solution with the nominal detergent concentration per the washer/disinfector manufacturer's instructions.
12. Place implants within the appropriate Acumed organizer. Remove all trays from their cases prior to processing. Place the organizers and implants in a suitable washer/disinfector basket. Place instruments in a suitable washer/disinfector basket. Minimally invasive surgery (MIS) injector or irrigation ports may be affixed to cannulated instruments so cleaning spray forces water through the instrument lumens. Using the cleaning solution, process implants and instruments through a standard washer/disinfector cycle in compliance to EN ISO 15883-1 and EN ISO 15883-2 or equivalent national standards. Dry for the nominal time recommended by the washer/disinfector manufacturer. Do not use a lubrication phase.
13. Visually inspect each article under normal lighting and without magnification for residual water in blind holes, visible deterioration, or loss of function. Pay particular attention to all challenging areas. If the device is visibly deteriorated or is unable to function, return it to Acumed.
14. Return the implants and instruments to their sterilization case and steam sterilize according to Sterilization instructions.

Thermal Disinfection

For thermal disinfection, remove trays from outer case and place separately in the washer/disinfector's rack during automated cleaning and thermal disinfection. Utilize a thermal disinfection program consistent with A0 value > 3000 or perform thermal disinfection at 90°C (maximum temperature 95°C) for at least 5 minutes.

Pre-Sterilization Inspection

- Visually inspect all devices under normal lighting to ensure that cleaning was effective. Pay close attention to all challenging areas. Reprocess instruments that are not clean and replace instruments that cannot be cleaned.
- Inspect the implants and instruments for surface damage, such as nicks, scratches, and cracks. Replace affected devices.
- Assess the instruments for proper use. Operate all parts and connecting mechanisms. Give careful attention to drivers, drill bits and reamers, and instruments used for cutting or implant insertion. Critically assess them for wear, sharpness, straightness, and corrosion. Replace any instrument that does not perform as intended.
- Inspect cutting edges under magnification. Replace worn instruments e.g., dull, chipped, cracked, rolled, or otherwise deformed. Run a cotton cloth over the edge to help detect chipping and cracking.
- Verify the legibility of all markings and reference scales. Replace any device that is unreadable.
- Repair, replace, and/or repeat the cleaning of instruments as needed to ensure proper operation before proceeding with sterilization.

- Lubricate (instrument milk) surgical instruments to increase useful life. Do not use silicone-based lubricants, oil, or grease, as these will interfere with steam sterilization. Use a water-based lubricant intended for use on surgical instruments and with steam sterilization. Use the lubricant as directed by the manufacturer. Use critical water if dilution is required.
- Fully replenish the system trays and caddies.

Sterility

- The OsteoMed ExtremiLock Ankle Plating System is supplied **NON-STERILE** and **STERILE**. **STERILE parts are expressly labeled as STERILE.**
- **DO NOT USE IF THE IMPLANT STERILE PACKAGE IS DAMAGED. DO NOT USE IMPLANTS AFTER THE EXPIRATION DATE.**
- Implants and instruments supplied non-sterile must be sterilized prior to use.
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- The user facility must clean and dry devices prior to sterilization per standard hospital procedures.
- Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017, the Association for the Advancement of Medical Instrumentation's (AAMI's) "Comprehensive guide to steam sterilization and sterility assurance in health care facilities."
- Non-sterile devices are sterilizable by steam sterilization (autoclaving). For sterilization of OsteoMed ExtremiLock Ankle System, the following parameters in the table below should be used. The table shows the minimum parameters validated to achieve a required Sterility Assurance Level (SAL) of 10^{-6} for the system.

Pre-Vacuum Steam Sterilization	ExtremiLock Ankle System Tray(s) with Screw Modules		ExtremiLock Ankle System Tray(s) with Screw Modules – Split Case Configuration		ExtremiLock Ankle System Screw Module (Cannulated and Solid Core)		ExtremiLock Ankle System Fusion Module	
Part No.	320-2920, 320-2930, 320-2931, 320-2940		320-2930, 320-2931, 320-2940		320-2924, 320-2925		320-2933	
Temperature ¹	132°C	134°C	132°C	134°C	132°C	134°C	132°C	134°C
Sterilization Time	4 minutes	3 minutes	4 minutes	3 minutes	4 minutes	3 minutes	4 minutes	3 minutes
Tray Preparation ²	Wrapped Tray		Wrapped Tray		Wrapped Block		Wrapped Block	
Wrapped Dry Time Minimum	60 minutes (Without Fusion Module)	120 minutes (With or Without Fusion Module)	60 minutes (Without Fusion Module)	80 minutes (With or Without Fusion Module)	25 minutes	75 minutes	40 minutes	75 minutes
	80 minutes (With Fusion Module)		80 minutes (With Fusion Module)					
Tray Preparation ³	Rigid Containers		Rigid Containers		Rigid Containers		Rigid Containers	
Rigid Container Dry Time Minimum	30 minutes	45 minutes	30 minutes	45 minutes	30 minutes	45 minutes	30 minutes	45 minutes
1. Do not exceed 135°C to avoid compromising functions of polymeric instrumentation.								
2. Wrapping technique recommendation: Individually wrapped in two layers of 1-ply polypropylene wrap (Kimguard KC600-510(k) K082554) using sequential techniques.								
3. OsteoMed recommends using Aesculap rigid containers & paper filters. The ExtremiLOCK Ankle Plating System (E-APS) must be removed from the outer case(s) and split and loaded into 2 rigid containers. The 2 upper insert trays from the E-APS shipping container are to be stacked up and placed inside an Aesculap full-sized basket, P/N JF224R (4-1/8" deep). Then the basket and contents are to be loaded into an Aesculap, 8-5/8" deep, full-size, rigid container, P/N JN446/JK486 (bottom/lid). The lower 3 modules from the E-APS shipping container, 320-2924, 632-2925, & 320-2933, are to be placed in a shorter Aesculap full-size basket, P/N JF221R (1-1/2" deep). Then this basket and contents are to be placed inside a 2nd Aesculap, 3-5/8" deep, full-size, rigid container, P/N JN441/JK486 (bottom/lid). The 2 rigid containers can be sterilized in the same unit at the same time. Use 4 Aesculap paper filters, P/N US751, 2 in the lid and 2 in the bottom for each rigid container.								
Note: Biological Indicator of <i>G. stearothermophilus</i> was used in the sterilization validations.								

Post-Sterilization Inspection



- Do not store or use sterile devices if they are not dry.
- Inspect the sterile barrier for signs of damage. Do not use the product if the sterile barrier has been compromised.

Storage

- Store sterilized devices under controlled conditions in a manner that minimizes the potential for contamination per ANSI/AAMI ST79:2017. Refer to sterilization wrap or rigid container manufacturer's IFU for limits on sterile product storage time and storage requirements for temperature and humidity.
- Devices should be stored in an area that provides protection from dust, pests, and temperature/humidity extremes.

Safe Disposal

Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling biohazardous material and disposal of sharp materials.

Symbols and Definitions					
Note: This is a general list of symbols. Not all symbols apply to all systems.					
Symbol	Description	ISO 15223-1	Symbol	Description	ISO 15223-1
	Manufacturer	5.1.1		Authorized Representative in the European Community	5.1.2
	Date of Manufacture	5.1.3		Use-by Date	5.1.4
	Batch Code (Lot Number)	5.1.5		Catalogue Number	5.1.6
	Serial number	5.1.7		Sterilized using irradiation	5.2.4
	Do Not Resterilize	5.2.6		Non-Sterile	5.2.7
	Single Sterile Barrier System	5.2.11		Single Sterile Barrier System with Protective packaging inside	5.2.14
	Double Sterile Barrier System	5.2.12		Single Sterile Barrier System with Protective packaging inside	5.2.13
 www.acumed.net/ifu	Consult Instructions for Use	5.4.3	Rx only	Caution: U.S. federal law restricts this device to sale by or on the order of a physician.	U.S. 21 CFR 801.109
	Do not use if package is damaged or opened and consult Instructions for Use	5.2.8		Do not re-use	5.4.2
	Medical Device	5.7.7		MR Conditional (ASTM F2503)	
	CE marking of conformity, Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment.				

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CUSTOMER SERVICE (OUTSIDE USA):

+1.888.627.9957

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

PKGI-014297-Rev C | 07-2025 | © 2025 Acumed® LLC

OsteoMed® is a registered trademark of Acumed, LLC.

Subsidiaries:

Acumed Ltd

Huebner House

The Fairground

Andover

Hampshire UK SP11 0QN

Tel: +44 1264 774450

Acumed Iberica

C/ Procion 1

Edificio Oficor

28023 Madrid-Spain

Tel: +34 913516357

Acumed GmbH

Fuhlsbuttler Stasse 300

22307 Hamburg

Deutschland

Tel: +49-40 947 82 093

EMERGO EUROPE

WESTERVOORTSEDIJK 60

6827 AT ARNHEM

THE NETHERLANDS

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ACUMED, LLC