

Surgical Technique



Acumed<sup>®</sup> is a global leader of innovative orthopaedic and medical solutions.

We are dedicated to developing products, service methods, and approaches that improve patient care.



## Radial Head Plating System

The Acumed Radial Head Plating System provides a method for the treatment of fractures when the radial head is salvageable. This system offers a family of precontoured plates designed for use in the "safe zone" of the radial head.

The Radial Head Plating System may be combined with any of the Acumed anatomic radial head replacement systems and may be found at the bottom of the tray. It may also be brought in as a stand-alone system at the surgeon's request. The Acutrak 2<sup>®</sup> Mini and Micro Headless Compression Screw instrumentation is also found at the bottom of each radial head replacement system tray, expanding the treatment options for the surgeon.

#### Indications for Use

The Radial Head Plating System is indicated for fractures, fusions, and osteotomies of the radius.

# Table of Contents

System Features	. 2
Instrument Overview	. 4
Surgical Technique Overview	. 6
Surgical Technique	. 8
Radial Head Plating System Surgical Technique	. 8
Ordering Information	12

# System Features

#### **Optimized Screw Angles**

Previous generations of radial head plates featured screws that came to a point in the center of the radial head. Acumed's radial head plates now have diverging and converging screw angles to capture bone fragments across the entire radial head. The screws are also strategically angled to prevent entering the articular surface of the radial head or colliding with one another, regardless of screw length selected.







#### Anatomic Plate Contour

The plate contour is designed to fit the anatomic contours of the radial head and neck with little or no intraoperative plate bending needed. The thickness of the plate varies along its length, providing a low-profile proximal portion to allow for closure of the annular ligament. The thicker neck portion of the plate helps provide support if there is a fracture line at the radial neck. The 2.3 mm locking and nonlocking screws are designed to sit flush with the plate's surface, designed to minimize hardware prominence and soft tissue irritation for the patient.

Both the 3-hole and 5-hole plate lengths are offered in two proximal curvatures to accommodate varying radial head anatomies. The standard curvature plates, colored gold, are intended for patients with a 25 mm or larger radial head diameter. The small curvature plates, colored silver, are intended for patients with a radial head diameter of 24 mm or less. By including a variety of lengths and plate curvatures, Acumed's system is designed to treat a wide range of radial head sizes and fracture patterns.

# System Features [continued]

#### Innovative Instrumentation

A radiolucent targeting guide is included in the Radial Head Plating System to assist with threading the locking drill guide into the proximal locking holes. Screw length is determined with either the laser marking on the drill bit or with the depth probe, allowing depth measurement to be read directly off of the drill guide.

K-wire holes are included in the plate for provisional fixation distally, and specially designed Plate Tacks may be used for provisional fixation proximally. The benefit of the short length of the Plate Tack is that it is designed to not interfere with the drilling and placement of the proximal screws. The distal tip of the Plate Tack is threaded, allowing it to maintain its provisional fixation until it is removed later in the procedure.





Two plate lengths and head curvatures accommodate variations in patient anatomy and fracture patterns.

## Instrument Overview



Small Plate Tack (80-0248)

**1.5 mm Hex Driver Tip** (HPC-0015)

1

**2.3 mm Bone Tap** (80-0362)



.035" Depth Probe (80-0357)



.045" x 6" Guide Wire (WS-1106ST) **2.0 mm Quick Release Drill** (80-0318)



Small Pointed Reduction Foreceps (OW-1200) **.054" x 6" Guide Wire** (WS-1406ST)



Locking Radial Head Targeting Guide (80-0246)

# Instrument Overview [continued]





**5-44 Targeting Guide Locking Bolt** (80-0247)

2.3 mm Screw Sleeve (MS-SS23)



Cruciform Driver Handle (MS-2210)



**2.0 mm Locking Drill Guide 4–32 mm** (80-0249)

2.0 mm Nonlocking Drill Guide Assembly 4–32 mm (80-0394)



Small Plate Bender (80-0363)

# Surgical Technique Overview

Articular Reduction and Plate Selection



Radial Head Plating System **Plate Placement** 



Provisional Plate Fixation



Nonlocking Distal Screw Fixation



Insert Locking Screw



Final Screw Placement







# Radial Head Plating System Surgical Technique





#### **Exposure and Dissection**

The most common approach is made with the forearm in neutral rotation, followed by exposing the radial head through the Kaplan interval in a line from the lateral epicondyle toward Lister's tubercle. Split the tendon for 3 to 4 cm distal to the epicondyle, with the distal extent limited by the posterior interosseous nerve.

Proximally release the extensor carpi radialis longus (ECRL) origin with the anterior capsule to permit direct access to the front of the radial head.

Figure 2



#### Articular Reduction and Plate Selection

Articular reduction is obtained with the Small Pointed Reduction Forceps (OW-1200) and provisionally held with .045" x 6" K-wires (WS-1106ST) placed to avoid interference with plate and screw placement in subsequent steps.

There are four plates in the system to choose from: 3-hole and 5-hole lengths, each with two different curvatures. The standard curvature plates are colored gold. The plates with the small curvature are for patients with a smaller radial head and are colored silver. Select the plate with the best fit for the patient and that best addresses the fracture pattern.

Small Pointed Reduction Foreceps (OW-1200)

Locking Radial **Head Plate** 3-Hole, Standard Curvature (31 mm) (70-0097)





(70-0099)



Locking Radial Head Plate 5-Hole, Standard Curvature (46 mm) (70-0098)



Locking Radial Head Plate 5-Hole, Small Curvature (46 mm) (70-0100)

# Radial Head Plating System Surgical Technique [continued]

### **Plate Placement**

Attach the Targeting Guide (80-0246) to the plate with the Locking Bolt (80-0247). Place the plate in the "safe zone" of the radial head. With the forearm in neutral rotation, the safe zone includes a 90° angle centered laterally, plus 20° anterior to this. For simplification, place the plate directly lateral with the forearm in 10° supination. The plate is designed to sit approximately half way onto the radial head's annular surface. This may avoid screw penetration of the head that would occur with more proximal placement. Small Plate Benders (80-0363) are included in the system if an adjustment to the plate's contour is needed.



Figure 3 A

Safe Zone

#### **Provisional Plate Fixation**

After placing the plate on the radius, it is provisionally held with the Small Plate Tack (80-0248) by inserting it through one of the proximal holes, or by placing a .054" x 6" K-wire (WS-1406ST) through the K-wire hole in the plate. The 5–44 Targeting Guide Locking Bolt (80-0247) is also cannulated to accept a .054" x 6" K-wire for provisional fixation.





Targeting Guide (80-0246) Small Plate Bender

Locking Radial

Head Plate

(80-0363)



Small Plate Tack (80-0248)

5–44 Targeting

Guide Locking



Figure 3 B

(WS-1406ST) Also used as a K-wire

# Radial Head Plating System Surgical Technique [continued]

Figure 5



#### Nonlocking Distal Screw Fixation

Attach the plate to the shaft with a screw through the slotted hole that allows for a slight adjustment of the plate position prior to final tightening of the screw at this step. Using the 2.0 mm Nonlocking Drill Guide Assembly 4-32 mm (80-0394) and the 2.0 mm Quick Release Drill (80-0318), drill bicortically through the slotted hole. Read the drill depth from the laser line on the drill or by inserting the .035" Depth Probe (80-0357) and hooking the far cortex. Insert a nonlocking screw (CO-N23XX) with the 1.5 mm Hex Driver (HPC-0015). A Bone Tap (80-0362) is available for use in patients with dense bone.

Figure 6



**Insert Locking Screws** Insert the 2.0 mm Locking Drill Guide 4-32 mm (80-0249) into one of the proximal holes and thread it into the plate. Threading the drill guide into the plate will ensure the hole is drilled in the proper trajectory so the locking screw (CO-T23XX) fully seats and locks into the plate. Drill with the 2.0 mm Quick Release Drill (80-0318). Drill depth can be read from the laser line on the drill or determined with the .035" Depth Probe (80-0357).

Repeat this process for inserting all proximal screws through the targeting guide, removing the Plate Tack as necessary.

2.0 mm Nonlocking Drill Guide Assembly 4-32 mm

(80-0394) 1.5 mm Hex Driver Tip (HPC-0015)



2.3 mm Bone Tap (80-0362)



2.0 mm Locking **Drill Guide** 4–32 mm (80-0249)

.035" Depth Probe



2.3 mm Nonlocking Screw (CO-N23XX)

2.3 mm Locking Cortical Screw (CO-T23XX)



# Radial Head Plating System Surgical Technique [continued]

#### Final Screw Placement

After inserting all proximal screws, remove the targeting guide. Then thread the 2.0 mm Locking Drill Guide 4–32 mm (80-0249) into the distal locking hole or holes to drill, measure depth, and insert the locking screw or screws (CO-T23XX).

Remove the K-wires that were used for provisional reduction of the fracture. Follow standard techniques for internal fixation procedures to close the wound/repair soft tissue, ensure proper hardware positioning and inspection of the joint.



#### Postoperative treatment typically is achieved by the overall management of the elbow and limb, more so than specifically the radial head. Early motion usually occurs within one to two days of surgery.

**Postoperative Protocol** 



To extract a Locking Radial Head Plate, use the 1.5 mm Hex Driver Tip (HPC-0015) and Cruciform Driver Handle (MS-2210) to remove all the screws in the plate. Referencing the Screw Removal Brochure (SPF10-00) may aid in implant extraction if difficulty is experienced.





**2.0 mm Locking Drill Guide 4–32 mm** (80-0249)



2.3 mm Locking Cortical Screw (CO-T23XX) Cruciform Driver Handle (MS-2210) 1.5 mm Hex Driver Tip (HPC-0015)

# Ordering Information

#### Implants Locking Radial Head Plate: 5-Hole, Locking Radial Head Plate: 3-Hole, 70-0098 70-0097 Standard Curvature (46 mm) Standard Curvature (31 mm) Locking Radial Head Plate: 5-Hole, Locking Radial Head Plate: 3-Hole, 70-0100 70-0099 Small Curvature (46 mm) Small Curvature (31 mm) Instrumentation Locking Radial Head Plate 5 Small Plate Tack 80-0248 80-0246 Targeting Guide 6) 1.5 mm Hex Driver Tip HPC-0015 5–44 Targeting Guide Locking Bolt 80-0247 2.3 mm Bone Tap 2.3 mm Screw Sleeve 80-0362 MS-SS23 .035" Depth Probe 80-0357 **Cruciform Driver Handle** MS-2210 2.0 mm Quick Release Drill 80-0318 2.0 mm Locking Drill Guide 4-32 mm 80-0249 2.0 mm Nonlocking Drill Guide 10 .054" x 6" Guide Wire\* WS-1406ST 80-0394 Assembly 4-32 mm .045" x 6" ST Guide Wire\* Small Plate Bender 11 WS-1106ST 80-0363 Small Pointed Reduction Foreceps OW-1200

\*Also used as a K-wire



# Ordering Information [continued]

2.3 mm Locking Cortical Screw	NS		
2.3 mm x 12 mm Locking Cortical Screw	CO-T2312	2.3 mm x 22 mm Locking Cortical Screw	CO-T2322
2.3 mm x 14 mm Locking Cortical Screw	CO-T2314	2.3 mm x 24 mm Locking Cortical Screw	CO-T2324
2.3 mm x 16 mm Locking Cortical Screw	CO-T2316	2.3 mm x 26 mm Locking Cortical Screw	CO-T2326
2.3 mm x 18 mm Locking Cortical Screw	CO-T2318	2.3 mm x 28 mm Locking Cortical Screw	CO-T2328
2.3 mm x 20 mm Locking Cortical Screw	CO-T2320		

#### 2.3 mm Nonlocking Screws

2.3 mm x 12 mm Nonlocking Screw	CO-N2312
2.3 mm x 14 mm Nonlocking Screw	CO-N2314
2.3 mm x 16 mm Nonlocking Screw	CO-N2316
2.3 mm x 18 mm Nonlocking Screw	CO-N2318
2.3 mm x 20 mm Nonlocking Screw	CO-N2320

2
4
6
8
2

Notes:	
	-
	-

Notes:

Notes.	

Notes:	



Acumed Headquarters 5885 NW Cornelius Pass Road Hillsboro, OR 97124 Office: +1.888.627.9957 Office: +1.503.627.9957 Fax: +1.503.520.9618 www.acumed.net These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located. Specific questions physicians may have about the availability and use of the products described on these materials should be directed to their particular authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

#### **ELB00-02-B** | Effective: 2017/04 | © 2017 Acumed<sup>®</sup> LLC