
Acutrak 3[®] Screw System

US	Instructions for Use	2
DE	Gebrauchsanweisung	14
EN	Instructions for use	27
ES	Instrucciones de uso	38
FR	Guide d'utilisation	50
IT	Istruzioni per l'uso	62
PT	Instruções de utilização.....	74

Instructions for Use

Acutrak 3® Screw System

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The US instructions are intended for users in the United States and its territories.

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POTENTIAL ADVERSE EFFECTS

- Anesthesia-related problems, problems with positioning of the patient (e.g., nausea, vomiting, neurological impairments, etc.) thrombosis, embolism, infection, or injury to other critical structures such as blood vessels, excessive bleeding, etc.
- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.
- Swelling, abnormal scarring, impairment in musculoskeletal function, hardware prominence, migration, loosening, bending, or breakage of the implant, malunion, nonunion or delayed union due to prolonged loading or excessive forces which may lead to implant failure and reoperation.
- Metal sensitivity, histological, allergic, or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu>.
- Injury to user.

Important: Any adverse event that has occurred in relation to the device should be reported to Acumed through customercomplaints@acumed.net or +1.888.627.9957

SURGICAL TECHNIQUE GUIDE (including Installation and Calibration)

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Consult the most recent versions of the Instructions for Use and Surgical Technique Guides as they are subject to change. Contact Acumed or an authorized agent for additional information.

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STERILITY

- Products are supplied NON-STERILE unless expressly labeled STERILE.
- **DO NOT USE IF THE STERILE PACKAGE IS DAMAGED. DO NOT USE AFTER EXPIRATION DATE.**
- Non-sterile devices are intended to be sterilized before use.
- The user facility must clean and disinfect devices prior to sterilization per standard hospital procedures.
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
Non-sterile devices are sterilizable by steam sterilization (autoclaving).

IMPLANTS

MATERIALS

- The implants are made of titanium alloy Ti-6Al-4V ELI per ASTM F136.
- Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu> for a detailed description of the material composition of Acumed metal implants.

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Single Use Instruments

- Do not reuse single use instruments. Reuse or reprocessing of devices labeled as "single use" can result in decreased mechanical and clinical performance of the devices which may result in patient harm.
- Reuse or reprocessing of single use instruments may create a risk of contamination (e.g. due to the transmission of infection material from one patient to another), which may result in patient harm.
- Dispose of single use instruments after use on a single patient, during a single procedure.

Multiple Use Instruments

- Multiple use instruments are intended for use on a single patient and a single procedure before requiring processing.
- Multiple (limited) use instruments have a limited lifespan. Carefully inspect the instruments before and after each procedure to ensure they are in proper operating condition. Instruments which are faulty, damaged, or suspect should not be used.

- Multiple use instruments potentially contaminated with transmissible spongiform encephalopathy (TSE) agents shall not be processed for reuse.

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IMPORTANT

- Promptly perform the processing steps to limit microbial growth and maximize the effectiveness of sterilization.
- Inspect implants for contamination by blood or tissue and dispose of them when found. Do not process contaminated implants.
- Prevent instrument corrosion by minimizing contact with solutions containing iodine, chlorine, and saline or other metal salts.
- Prevent damage to the protective anodization layer on aluminum instruments by avoiding contact with solutions < 4 pH and > 9 pH, especially if they contain sodium carbonate or sodium hydroxide.
- Repeated processing of anodized metals may cause colors to fade but this does not affect the function of the device.
- Promptly perform the processing steps to limit microbial growth and maximize the effectiveness of sterilization.
- Inspect implants for contamination by blood or tissue and dispose of them when found. Do not process contaminated implants.
- Prevent instrument corrosion by minimizing contact with solutions containing iodine, chlorine, and saline or other metal salts.
- Prevent damage to the protective anodization layer on aluminum instruments by avoiding contact with solutions < 4 pH and > 9 pH, especially if they contain sodium carbonate or sodium hydroxide.
- Repeated processing of anodized metals may cause colors to fade but this does not affect the function of the device.
- Avoid cleaning agents containing aldehydes since they can denature and coagulate proteins (fixation).
- Enzymatic detergents are well suited for loosening protein-based contamination.
 - Use a neutral pH enzymatic detergent.
 - Use a low foaming solution to allow visibility of the device during cleaning.

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1. Rinse the contaminated instruments under running cold utility water to reduce heavy surface contamination.
 2. Dispose of any used instruments intended for single use only.
 3. Place the contaminated instruments in enzymatic solution* until completely submerged to minimize the spraying of solution.
 4. Actuate all moveable parts to allow detergent to contact all surfaces.
 5. Soak for a minimum of ten (10) minutes.
 6. Scrub the instruments using a soft-bristled brush to remove all visible debris. Do not use stainless steel or other abrasives as these may damage the surface.
 - When possible, scrub the instruments when totally submerged to minimize the spraying of fluid.
 7. Some instruments may require special consideration:
 - Clean the instruments with all parts loosened. Clean the instruments disassembled if they are designed to be taken apart.
 - Use a water jet to flood cleaning solution into challenging areas, such as mating surfaces, springs, coils, cannulations, blind holes, flutes, cutting teeth, and flexible parts to flush out any trapped soil.
 - Operate movable parts and rotate (as necessary) while scrubbing to ensure that all crevices are accessible.
 - Carefully clean cannulated parts and challenging areas using an appropriately sized brush.
 - Optionally, sonicate for 10 to 15 minutes using a fresh, neutral pH ultrasonic cleaning solution. Follow the ultrasonic cleaner and detergent manufacturer's instructions.
- Important:** Any previous surface damage may increase due to ultrasonic cleaning.
8. Perform an initial rinse for at least 3 minutes using clean, soft, utility water in the temperature range of 25°C to 35°C (77°F to 95°F) to remove all signs of contamination and cleaning agent.
 - Actuate all movable parts.
 - Flush out cannulations and complex mechanisms.
 9. Repeat the previous processing steps if visible residue remains present.

10. Perform a final rinse for at least 1 minute using critical water to displace minerals and other impurities found in utility water. Do not use saline solutions for final rinsing because they may interfere with disinfection and sterilization.
 - Actuate all moving parts.
 - Pay particular attention to cannulations and blind holes as well as hinges and joints between mating surfaces.
 - Rinse cannulations at least three times with a syringe (volume 1-50ml).
 11. Remove excess moisture from the instruments using a clean, absorbent, non-shedding wipe.
 12. Allow the instruments to thoroughly dry. Any moisture may affect sterilization, and devices may remain wet after the drying period.
- * Manual cleaning was validated using STERIS Prolystica 2X Concentrate Enzymatic Presoak and Cleaner.

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Important: Contaminated instruments require manual cleaning prior to automated cleaning. Follow the washer-disinfector manufacturer's instructions explicitly.

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PRE-STERILIZATION INSPECTION

- Visually inspect all devices under normal lighting to ensure that cleaning was effective. Pay close attention to all challenging areas.
 - Re-process instruments that are not clean.
 - Replace instruments that cannot be cleaned.
- Inspect the implants and instruments for surface damage, such as nicks, scratches, and cracks. Replace affected devices.
- Assess the instruments for proper use. Operate all parts and connecting mechanisms. Give careful attention to drivers, drill bits and reamers, and instruments used for cutting or implant insertion. Critically assess them for wear, sharpness, straightness, and corrosion. Replace any instrument that does not perform as intended.
- Inspect all cutting edges under magnification. Replace worn instruments e.g., dull, chipped, cracked, rolled, or otherwise deformed. Run a cotton cloth over the edge to help detect chipping and cracking.
- Verify the legibility of all markings and reference scales. Replace any device that is unreadable.
- Repair, replace, and/or repeat the cleaning of instruments as needed to ensure proper operation before proceeding with sterilization.

- Lubricate (instrument milk) surgical instruments to increase useful life. Do not use silicone-based lubricants, oil, or grease, as these will interfere with steam sterilization. Use a water-based lubricant intended for use on surgical instruments and with steam sterilization. Use the lubricant as directed by the manufacturer. Use critical water if dilution is required.
- Fully replenish the system trays and caddies.
- Disassemble the Acutrak 3 Screw Sizer into the Piston (80-4166) and Body (80-4165) and store within the instrument tray (80-4178 or 80-4179) as indicated on the tray for sterilization.

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* Association for the Advancement of Medical Instrumentation (AAMI). Comprehensive guide to steam sterilization and sterility assurance in health care facilities. AAMI ST79:2017. Arlington, VA.

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Acutrak 3 System:

	The <i>Acutrak 3 System</i> is validated for steam sterilization using one of the below sets of parameters. ³	
Pre-Vacuum Steam Sterilization	Storage Case Part Numbers: Base: (80-4173, 80-4175) Case Lid: (80-4174, 80-4176, 80-4177)	
Exposure Temperature:	270°F (132°C)	273°F (134°C) ²
Exposure Time²	4 minutes	3 minutes
Condition:¹	Wrapped tray ¹	
Wrapped Dry Time Minimum:	30 minutes	30 minutes
Wrapped Configuration Details:	Wrap in two layers of 1- ply polypropylene wrap using sequential envelope techniques.	
¹ US customers must use FDA-cleared sterilization packaging/wrap and other accessories appropriate for the cycle parameters recommended in these instructions. Refer to PKGI-76 at www.acumed.net/ifu for sterilization in Aesculap® rigid sterilization containers.		
² The devices are compatible with exposure for 18 minutes at 134°C.		
³ Note: Biological indicator, <i>G. stearothermophilus</i> , was used in sterilization validation.		

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Symbols Glossary

Symbol	Description	ISO 15223-1
 www.acumed.net/ifu	Consult the electronic instructions for use (eIFU) at www.acumed.net/ifu	5.4.3
STERILE R	Sterilized using irradiation	5.2.4
	Double sterile barrier system	5.2.12
	Non-sterile	5.2.7
	Use-by date	5.1.4
REF	Catalogue number	5.1.6
LOT	Batch code	5.1.5
	Manufacturer	5.1.1
	Date of manufacture	5.1.3
	Do not re-sterilize	5.2.6
	Do not re-use	5.4.2
	Do not use if package is damaged and consult instructions for use / do not use if the product sterile barrier system or its packaging is compromised.	5.2.8
R ONLY	Caution: U.S. federal law restricts this device to sale by or on the order of a physician.	U.S. 21 CFR 801.109
	The reticle is a registered trademark of Acumed. It may appear alone or with the Acumed name.	
	MR Conditional: An item with demonstrated safety in the MR environment within defined conditions.	

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Instructions for use

Acutrak 3® Screw System

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The implants and instruments are intended only for use by a licensed and trained health care professional.

DEVICE DESCRIPTION

The Acutrak 3 System consists of headless compression screws intended to provide fixation for fractures and osteotomies during the healing process. The system includes instrumentation designed to assist with implant insertion. Screws are available in various diameters and lengths to accommodate anatomical location, patient size, and fracture pattern. The system is available in four sizes: Nano, Micro, Mini, and Standard.

INTENDED PURPOSE

The Acutrak 3 Headless Compression Screw System is intended as a fixation device designed to provide fixation of various fractures and osteotomies while they heal.

INDICATIONS AND CLINICAL BENEFIT

The Acutrak 3 Screw System is indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, fracture repair, and fracture fixation of bones appropriate for the size of the device. The screws are not intended for interference or soft tissue fixation.

The intended clinical benefit of the Acutrak 3 Screw System is to provide fixation for patients with fractures as a result of trauma, osteotomy or reconstruction. The indirect clinical benefit of the encompassed instrumentation is to facilitate implantation protocols for fracture fixation or bone fusion.

TARGET POPULATION AND INTENDED USERS

The Acutrak 3 Screw System is intended to be used on skeletally mature, adult patients with sufficient quality and quantity of bone who require surgical intervention to achieve fracture fixation, fusion or reconstruction of bones of the appropriate size of the device.

CONTRAINDICATIONS

- Active or latent infection or marked inflammation of the treatment area.
- Sepsis.
- Insufficient quantity or quality of bone, osteoporosis, or in patients with certain metabolic diseases.
- Patients with confirmed material sensitivity.
- Patients who are unwilling or incapable of following post-operative care instructions.
- Soft tissue fixation.

These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS and PRECAUTIONS

Warnings

- Treatment or implant may fail, including sudden failure, because of:
 - Loose fixation and/or loosening
 - Stress, including stress from inappropriate bending of the implant during surgery
 - Stress concentrations
 - Stress of weight bearing, load bearing, or excessive activity

Failure is more likely if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing. Failure is more likely if the patient does not follow post-operative care instructions.

- Nerve or soft tissue damage may result from surgical trauma or the presence of an implant.
- Instrument breakage or damage, as well as tissue damage, may occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use, or unintended use.
- Devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects.
- Fractured implants should be removed from the patient during surgery. If unable to remove, notify the patient.
- Implants may cause distortion and/or block the view of anatomic structures on radiographic images.

Precautions

- Implants and instruments are intended only for professional use by a licensed physician.
- Use devices and instruments in accordance with the Instructions for Use.
- Do not use the sterile product past the use-by date. Refer to the device label.
- Do not use or re-sterilize an implant provided in sterile packaging if the package has been damaged. The sterility may be compromised, and the cleanliness of the implant may be uncertain. Report damaged packaging to your distributor or Acumed.
- Use of this system has not been evaluated in children or individuals that are not skeletally mature.
- Inspect all components preoperatively to ensure utility. Do not attempt a surgical procedure with faulty, damaged, or suspect instruments or implants. Alternate fixation methods should be available intraoperatively.
- Reprocessing and/or reuse of single use devices may result in infection/cross contamination, and/or sudden failure due to previous stresses.
- Use of implant components from different manufacturers has not been evaluated.
- Use of chemical disinfection may leave residues that adversely affect steam sterilization.
- Steam penetration and device sterilization may be negatively affected by labeling that block tray steam holes.
- Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.
- Do not resharpen drill bits or reamers as these devices have critical dimensions and geometries that cannot be restored once the instrument has been consumed.

POTENTIAL ADVERSE EFFECTS

- Anesthesia-related problems, problems with positioning of the patient (e.g., nausea, vomiting, neurological impairments, etc.) thrombosis, embolism, infection, or injury to other critical structures such as blood vessels, excessive bleeding, etc.

- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.
- Swelling, abnormal scarring, impairment in musculoskeletal function, hardware prominence, migration, loosening, bending, or breakage of the implant, malunion, nonunion or delayed union due to prolonged loading or excessive forces which may lead to implant failure and reoperation.
- Metal sensitivity, histological, allergic, or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu>.
- Injury to user.

Important: Any serious incident that has occurred in relation to the device should be reported to Acumed (through customercomplaints@acumed.net or +1.888.627.9957) and the Competent Authority of the Member State in which the user and/or patient is established.

SURGICAL TECHNIQUE GUIDE

Acumed offers one or more Surgical Technique Guides to promote the safe and effective use of this system. Consult our Surgical Technique Guides at <https://www.acumed.net/resources/documents/>. The Surgical Technique Guides also include information about the installation and calibration of the device.

Important: Surgical Technique Guides may contain important safety information.

Important: The instruments and implants in this system are intended to be used by suitably trained and qualified surgeons in a hospital operating room setting. Before treatment, the surgeon is advised to read and fully understand all instructions and communicate to the patient any relevant medical information provided therein, including the use, limitations, risks (safety communications), and possible adverse effects of the proposed treatment.

Consult the most recent versions of the Instructions for Use and Surgical Techniques as they are subject to change. Contact Acumed or an authorized agent to request any additional information.

MRI SAFETY INFORMATION

Many Acumed implants have been evaluated for safety in the MR environment. Consult our publication "Acumed Implants in the MR Environment" at <https://www.acumed.net/resources/documents/>.

LIFETIME

- Once installed, implants are expected to provide fixation, fusion and physiological support and have an effective life during bone healing. The implants are biocompatible and may remain implanted at the discretion of the surgeon or patient.
- Multiple-use instruments have a lifetime that is affected by usage, handling, and processing. Assess multiple-use instruments for fitness during pre-sterilization inspection.

STERILITY

- Products are supplied NON-STERILE unless expressly labeled STERILE.
- **DO NOT USE IF THE STERILE PACKAGE IS DAMAGED. DO NOT USE AFTER EXPIRATION DATE.**
- Non-sterile devices are intended to be sterilized before use.
- The user facility must clean and disinfect devices prior to sterilization per standard hospital procedures.
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.

- Non-sterile devices are sterilizable by steam sterilization (autoclaving).

IMPLANTS

MATERIALS

- The implants are made of titanium alloy Ti-6Al-4V ELI per ASTM F136
- Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu> for a detailed description of the material composition of Acumed metal implants.

SINGLE USE

- Implants are intended for single use only, as indicated on the label. This product has not been validated for multiple uses. Reusing single use implants may increase the risk of failure and cross- contamination.
- Dispose of any unused implant that is contaminated with human blood or tissue in accordance with hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material. Do not process a contaminated implant.

IMPORTANT

- Physiological dimensions limit implant sizes. Select the type and size of implant that best meets the patient's requirements for close adaptation and firm seating with adequate support.
- Implants are designed to hold fractured bones in place during bone healing at the site of fracture. Implants are not designed for excessive loadbearing. Improper selection or improper implantation of the device may increase the possibility of loosening or migration.
- Do not bend the implant except as indicated in the corresponding Surgical Technique Guide. Repeated or excessive bending may weaken the implant and cause failure at a later time.
- Only combine implants when they are intended for that purpose.
- Protect implants against scratching and nicking to prevent stress concentrations, which can result in failure.
- Prevent unused implants from becoming soiled.
- The color of anodized implants may change over time due to processing. This color change does not affect the mechanical properties of the implants.

INSTRUMENTS

MATERIALS

The instruments are manufactured from various grades of titanium, stainless steel, aluminum, silicone, and/or medical grade plastic.

Single Use Instruments

- Do not reuse single use instruments. Reuse or reprocessing of devices labeled as “single use” can result in decreased mechanical and clinical performance of the devices which may result in patient harm.
- Reuse or reprocessing of single use instruments may create a risk of contamination (e.g. due to the transmission of infection material from one patient to another), which may result in patient harm.
- Dispose of single use instruments after use on a single patient, during a single procedure

Multiple Use Instruments

- Multiple use instruments are intended for use on a single patient and a single procedure before requiring processing.
- Multiple (limited) use instruments have a limited lifespan. Carefully inspect the instruments before and after each procedure to ensure they are in proper operating condition. Instruments which are faulty, damaged, or suspect should not be used.
- Multiple use instruments potentially contaminated with transmissible spongiform encephalopathy (TSE) agents shall not be processed for reuse.

IMPORTANT

- Protect instruments against scratching and nicking to prevent stress concentrations, which can lead to instrument failure.
- Near the point of use: Wipe excess contamination from instruments and prevent any soil from drying. Instruments with substantial or dried soil are particularly difficult to reliably process. Transport contaminated instruments for processing as soon as possible after use.
- Avoid prolonged instrument contact with iodine and saline.
- Handle and transport soiled instruments in a manner that avoids contamination of any unused implants.

PROCESSING

Important: Processing personnel must be qualified with suitable training and experience. Use proper personal protective equipment (PPE) when working with contaminated devices.

IMPORTANT

- Promptly perform the processing steps to limit microbial growth and maximize the effectiveness of sterilization.
- Inspect implants for contamination by blood or tissue and dispose of them when found. Do not process contaminated implants.
- Prevent instrument corrosion by minimizing contact with solutions containing iodine, chlorine, and saline or other metal salts.
- Prevent damage to the protective anodization layer on aluminum instruments by avoiding contact with solutions < 4 pH and > 9 pH, especially if they contain sodium carbonate or sodium hydroxide.
- Repeated processing of anodized metals may cause colors to fade but this does not affect the function of the device.
- Avoid cleaning agents containing aldehydes since they can denature and coagulate proteins (fixation).
- Enzymatic detergents are well suited for loosening protein-based contamination.
 - Use a neutral pH enzymatic detergent.
 - Use a low foaming solution to allow visibility of the device during cleaning.

- Closely follow the manufacturer's instructions for the safety, storage, mixing, water quality, exposure time, temperature, replacement, and disposal of cleaning agents.
- Devices potentially contaminated with transmissible spongiform encephalopathy (TSE) agents shall not be processed or reused. These processing instructions are not suitable for the inactivation of TSE agents. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling and disposal of devices potentially contaminated with TSE agents.
- **Utility water:** Refer to AAMI TIR34* when instructed to use utility water. Utility water is typically municipal or tap water but may require additional treatment to be suitable for use.
- **Critical water:** Refer to AAMI TIR34* when instructed to use critical water. Critical water is highly treated and has very low organic and inorganic content with an endotoxin level under 10 EU/mL. Suitable water may also be specified in national pharmacopeias, national standards, and hospital protocols.

* Association for the Advancement of Medical Instrumentation (AAMI). Water for the reprocessing of medical devices. AAMI TIR34:2014/(R)2017. Arlington, VA.

MANUAL CLEANING

Important: Implants are subjected to automated cleaning only. Contaminated instruments require manual cleaning prior to automated cleaning.

1. Rinse the contaminated instruments under running cold utility water to reduce heavy surface contamination.
 2. Dispose of any used instruments intended for single use only.
 3. Place the contaminated instruments in enzymatic solution* until completely submerged to minimize the spraying of solution.
 4. Actuate all moveable parts to allow detergent to contact all surfaces.
 5. Soak for a minimum of ten (10) minutes.
 6. Scrub the instruments using a soft-bristled brush to remove all visible debris. Do not use stainless steel or other abrasives as these may damage the surface.
 - When possible, scrub the instruments when totally submerged to minimize the spraying of fluid.
 7. Some instruments may require special consideration:
 - Clean the instruments with all parts loosened. Clean the instruments disassembled if they are designed to be taken apart.
 - Use a water jet to flood cleaning solution into challenging areas, such as mating surfaces, springs, coils, cannulations, blind holes, flutes, cutting teeth, and flexible parts to flush out any trapped soil.
 - Operate movable parts and rotate (as necessary) while scrubbing to ensure that all crevices are accessible.
 - Carefully clean cannulated parts and challenging areas using an appropriately sized brush.
 - Optionally, sonicate for 10 to 15 minutes using a fresh, neutral pH ultrasonic cleaning solution. Follow the ultrasonic cleaner and detergent manufacturer's instructions.
- Important:** Any previous surface damage may increase due to ultrasonic cleaning.
8. Perform an initial rinse for at least 3 minutes using clean, soft, utility water in the temperature range of 25°C to 35°C (77°F to 95°F) to remove all signs of contamination and cleaning agent.
 - Actuate all movable parts.

- Flush out cannulations and complex mechanisms.
9. Repeat the previous processing steps if visible residue remains present.
 10. Perform a final rinse for at least 1 minute using critical water to displace minerals and other impurities found in utility water. Do not use saline solutions for final rinsing because they may interfere with disinfection and sterilization.
 - Actuate all moving parts.
 - Pay particular attention to cannulations and blind holes as well as hinges and joints between mating surfaces.
 - Rinse cannulations at least three times with a syringe (volume 1-50ml).
 11. Remove excess moisture from the instruments using a clean, absorbent, non-shedding wipe.
 12. Allow the instruments to thoroughly dry. Any moisture may affect sterilization, and devices may remain wet after the drying period.
- * Manual cleaning was validated using STERIS Prolystica 2X Concentrate Enzymatic Presoak and Cleaner.

AUTOMATED CLEANING – THERMAL DISINFECTION

Important: Contaminated instruments require manual cleaning prior to automated cleaning.
Follow the washer-disinfector manufacturer's instructions explicitly.

- Process the devices using a standard washing and thermal disinfection cycle in a washer-disinfector compliant with EN ISO 15883-1 and EN ISO 15883-2 or equivalent national standards.
 - Remove all implants if a lubrication phase will be used.
 - Process all trays removed from the case.
 - Attach the following instruments to the minimally invasive surgery (MIS) injector or irrigation ports:
 - 80-0663 Medium Ratcheting Driver Handle
 - 80-4071 - Small Ratcheting Driver Handle
 - Thermal disinfection has been validated* for an $A_0 \geq 3000$ (at least 5 minutes at 90°C).
- * Washer-disinfector processing was validated for fully loaded trays with all parts placed appropriately, with recommended evaluations as listed in EN ISO 15883-1 and EN ISO 15883-2 using a STERIS Reliance® Genfore Washer-Disinfector and STERIS Prolystica 2X Enzymatic Detergent.

PRE-STERILIZATION INSPECTION

- Visually inspect all devices under normal lighting to ensure that cleaning was effective. Pay close attention to all challenging areas.
 - Re-process instruments that are not clean.
 - Replace instruments that cannot be cleaned.
- Inspect the implants and instruments for surface damage, such as nicks, scratches, and cracks. Replace affected devices.
- Assess the instruments for proper use. Operate all parts and connecting mechanisms. Give careful attention to drivers, drill bits and reamers, and instruments used for cutting or implant insertion. Critically assess them for wear, sharpness, straightness, and corrosion. Replace any instrument that does not perform as intended.
- Inspect all cutting edges under magnification. Replace worn instruments e.g., dull, chipped, cracked, rolled, or otherwise deformed. Run a cotton cloth over the edge to help detect chipping and cracking.

- Verify the legibility of all markings and reference scales. Replace any device that is unreadable.
- Repair, replace, and/or repeat the cleaning of instruments as needed to ensure proper operation before proceeding with sterilization.
- Lubricate (instrument milk) surgical instruments to increase useful life. Do not use silicone-based lubricants, oil, or grease, as these will interfere with steam sterilization. Use a water-based lubricant intended for use on surgical instruments and with steam sterilization. Use the lubricant as directed by the manufacturer. Use critical water if dilution is required.
- Fully replenish the system trays and caddies.
- Disassemble the Acutrak 3 Screw Sizer into the Piston (80-4166) and Body (80-4165) and store within the instrument tray (80-4178 or 80-4179) as indicated on the tray for sterilization.

STERILIZATION

- Perform sterilization using a dynamic-air-removal (pre-vacuum) autoclave.
 - Gravity displacement sterilization is not recommended.
 - Immediate use (flash) sterilization is not recommended.
- Ensure the sterilizer's maximum load limit is not exceeded when sterilizing multiple sets or devices.
- Do not stack containers as this might prevent the penetration of steam and inhibit drying.
- Refer to the sterilizer manufacturer's instructions and ensure proper installation, calibration, use, and ongoing maintenance.
- The sterilized items should be allowed to cool to room temperature before handling. This allows for safe handling and prevents condensation.
- Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017*.

* Association for the Advancement of Medical Instrumentation (AAMI). Comprehensive guide to steam sterilization and sterility assurance in health care facilities. AAMI ST79:2017. Arlington, VA

The following table shows the minimum parameters validated* to achieve a required Sterility Assurance Level (SAL) of 10^{-6} for the system.

Important:

- Sterilization parameters are only valid for devices that have been cleaned per these instructions and are thoroughly dry.
- Sterilization parameters are only valid when the devices are properly housed in the Acumed storage case part numbers identified in the table.

Acutrak 3 System:

The <i>Acutrak 3 System</i> is validated for steam sterilization using one of the below sets of parameters. ³		
Pre-Vacuum Steam Sterilization	Storage Case Part Numbers: Base: (80-4173, 80-4175) Case Lid: (80-4174, 80-4176, 80-4177)	
Exposure Temperature:	270°F (132°C)	273°F (134°C) ²
Exposure Time²	4 minutes	3 minutes
Condition:¹	Wrapped tray ¹	
Wrapped Dry Time Minimum:	30 minutes	30 minutes
Wrapped Configuration Details:	Wrap in two layers of 1- ply polypropylene wrap using sequential envelope techniques.	

¹ US customers must use FDA-cleared sterilization packaging/wrap and other accessories appropriate for the cycle parameters recommended in these instructions. Refer to PKGI-76 at www.acumed.net/ifu for sterilization in Aesculap® rigid sterilization containers.

² The devices are compatible with exposure for 18 minutes at 134°C.

³ Note: Biological indicator, G. steartophilus, was used in sterilization validation.

*Sterilization was validated using a STERIS Amsco 3023 Vacumatic Prevacuum sterilizer and KimGuard KC600 One-Step wrap.

POST-STERILIZATION INSPECTION

- Do not store or use sterile devices if they are not dry.
 - Moisture supports the survival of microorganisms.
 - Moisture remaining on wrapped or contained products after sterilization could compromise the sterile barrier.
 - Moisture can corrode metal and dull sharp edges.
- Inspect the sterile barrier for signs of damage. Do not use the product if the sterile barrier has been compromised.

STORAGE CONDITIONS

STORAGE OF SETS AFTER STEAM STERILIZATION

- Store sterilized devices under controlled conditions in a manner that minimizes the potential for contamination per ANSI/AAMI ST79:2017.
- Refer to sterilization wrap or rigid container manufacturer's IFU for limits on sterile product storage time and storage requirements for temperature and humidity.

STORAGE OF PACKAGED NON-STERILE AND STERILE PRODUCT

- Devices should be stored in an area that provides protection from dust, pests, and temperature/humidity extremes.

SAFE DISPOSAL

Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance Report for the implants may be obtained from the European dataset on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Use the following search terms: 0806378BUDI76QK

Important: Any serious incident that has occurred in relation to the device should be reported to Acumed and the competent authority of the Member State in which the user and/or patient is established. To report a complaint about this product, call Customer Service at +1.888.627.9957 or send an email to customercomplaints@acumed.net.

Symbols Glossary

Symbol	Description	ISO 15223-1
 www.acumed.net/ifu	Consult the electronic instructions for use (eIFU) at www.acumed.net/ifu	5.4.3
	Sterilized using irradiation	5.2.4
	Double sterile barrier system	5.2.12
	Non-sterile	5.2.7
	Use-by date	5.1.4
	Catalogue number	5.1.6
	Batch code	5.1.5
	Authorized representative in the European Community / European Union	5.1.2
	Medical device	5.7.7
	Manufacturer	5.1.1
	Date of manufacture	5.1.3
	Do not resterilize	5.2.6
	Do not re-use	5.4.2
	Unique Device Identifier: Indicates a carrier that contains Unique Device Identifier information.	5.7.10
	Do not use if package is damaged and consult instructions for use / do not use if the product sterile barrier system or its packaging is compromised.	5.2.8
	MR Conditional: An item with demonstrated safety in the MR environment within defined conditions. (ASTM F2503 7.3.2).	
	The reticle is a registered trademark of Acumed. It may appear alone or with the Acumed name.	
	CE marking of conformity, Article 17 of EU Directive 93/42/EEC or Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment.	

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These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

PKGI-012778-Rev D | Effective 06-2025 | © 2023 Acumed® LLC

Acutrak 3 ® is a registered trademark of Acumed, LLC.

Subsidiaries:**Acumed Ltd**

Huebner House
The Fairground
Andover
Hampshire UK SP11 0QN
Tel: +44 1264 774450

Acumed Iberica

C/ Proción 1
Edificio Oficor
28023 Madrid - Spain
Tel: +34 913516357

Acumed Beijing

Room A1206, Horizon International
Tower
No. 6, Zhichun Road
Haidian District
100088 Beijing, China
Tel: +86 10 82001303

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Acumed GmbH