## Acu-Sinch® Repair System

US	Instructions for use .....	2
DA	Brugsanvisning.....	8
DE	Gebrauchsanweisung.....	15
EL	Οδηγίες χρήσης.....	22
EN	Instructions for use .....	30
ES	Instrucciones de uso .....	37
FI	Käyttöohjeet .....	44
FR	Guide d'utilisation .....	51
IT	Istruzioni per l'uso .....	58
NB	Bruksanvisning .....	65
NL	Gebruiksaanwijzing.....	72
PT	Instruções de utilização.....	79
SV	Bruksanvisning .....	86
TR	Kullanım talimatları.....	93

# Instructions for use

## ACU-SINCH® REPAIR SYSTEM

**These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The US instructions are intended for users in the United States and its territories.**

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### INTENDED USE/ INDICATIONS FOR USE

The Acu-Sinch® Repair System is indicated to be used in conjunction with a clavicle plate of the Congruent Bone Plate System to provide fixation during the healing of clavicle fractures.

### CONTRAINDICATIONS

- Active or latent infection or marked inflammation of the treatment area.
- Sepsis.
- Osteoporosis, insufficient quantity or quality of bone and/or soft tissue.
- Patients with material sensitivity. If sensitivity is suspected, tests are to be performed prior to implantation.
- Patients who are unwilling or incapable of following post-operative care instructions.
- The device is not indicated as a sole treatment for chronic ligament and tendon injuries.
- The system may not be suitable for skeletally immature patients and must not disturb the growth plate.

### WARNINGS and PRECAUTIONS

#### Warnings

- For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device.
- Treatment or implant may fail, including sudden failure, because of:
  - Loose fixation and/or loosening
  - Stress, including stress from inappropriate bending of the implant during surgery
  - Stress concentrations
  - Stress of weight bearing, load bearing, or excessive activity
- Device breakage (including suture breakage) is more likely if the implant experiences increased loads due to delayed union, nonunion, or incomplete soft tissue healing.

- The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant, and the possibility of soft tissue damage related to either surgical trauma (i.e. damage to neurovascular structures during preparation such as drilling or implant insertion) or purely due to the presence of an implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.
- The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity.
- Nerve or soft tissue damage may result from surgical trauma or the presence of an implant.
- Improper insertion of the device during implantation can increase the possibility of implant loosening and migration, neurovascular injury, or bony fracture.
- Devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects.
- Implants may cause distortion and/or block the view of anatomic structures on radiographic images.
- Drilling into or through bone may result in immediate or delayed bone fracture.

### **Precaution**

- Implants and instruments are intended only for professional use by a licensed physician.
- Use devices and instruments in accordance with the Instructions for Use.
- Do not use the sterile product past the use-by date. Refer to the device label.
- Do not use or re-sterilize an implant provided in sterile packaging if the package has been damaged. The sterility may be compromised, and the cleanliness of the implant may be uncertain. Report damaged packaging to your distributor or Acumed.
- Use of this system has not been evaluated in children or individuals that are not skeletally mature.
- Inspect all components preoperatively to ensure utility. Do not attempt a surgical procedure with faulty, damaged, or suspect instruments or implants. Alternate fixation methods should be available intraoperatively.
- Reprocessing and/or reuse of single use devices may result in infection/cross contamination, and/or sudden failure due to previous stresses.
- Use of implant components from different manufacturers has not been evaluated.
- Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.
- Inspect suture for signs of damage or fraying.
- Prior to use, inspect the product package for signs of tampering, water contamination, and seal integrity.
- The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures.
- Revision surgeries with implants are not uncommon.

### **POTENTIAL ADVERSE EFFECTS**

- Anesthesia-related problems, problems with positioning of the patient (e.g. nausea, vomiting, neurological impairments, etc.) thrombosis, embolism, infection or injury to other critical structures such as blood vessels, excessive bleeding, etc.
- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.

- Swelling, abnormal scarring, impairment in musculoskeletal function, hardware prominence, migration, loosening, bending, or breakage of the implant or suture, malunion, nonunion or delayed union due to prolonged loading or excessive forces during insertion and during the healing period as well as inadequate postoperative protocol =may lead to implant failure and reoperation.
- Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu>.
- Injury to user.

**Important:** Any adverse event that has occurred in relation to the device should be reported to Acumed through [customercomplaints@acumed.net](mailto:customercomplaints@acumed.net) or +1.888.627.9957

### **SURGICAL TECHNIQUE GUIDE (including Installation and Calibration)**

Acumed offers one or more Surgical Technique Guides to promote the safe and effective use of this system. Consult our Surgical Technique Guides at <https://www.acumed.net/resources/documents/>. The Surgical Technique Guides also include information about the installation and calibration of the device.

**Important:** Surgical techniques may contain important safety information.

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Consult the most recent versions of the Instructions for Use and Surgical Technique Guides as they are subject to change. Contact Acumed or an authorized agent for additional information.

### **MRI SAFETY INFORMATION**

Many Acumed implants have been evaluated for safety in the MR environment. Consult our publication "Acumed Implants in the MR Environment" at [www.acumed.net/ifu](http://www.acumed.net/ifu) for more information.

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### **STERILITY**

- Implants and instruments are provided sterile.
- Sterility was achieved by using the ethylene oxide gas sterilization method to Sterility Assurance Level (SAL) of 10<sup>-6</sup>.
- Acumed does not recommend resterilization of sterile-packaged product.
- If sterile packaging is damaged, the incident must be reported to Acumed.
- The product must not be used and must be returned to Acumed.

## IMPLANTS

### MATERIALS

- **Implants:** The metal implants are made out of implantable grade titanium.
- **Suture:** The suture is composed of Ultra High Molecular Weight Polyethylene (UHMWPE).
- Consult our document “Metal Sensitivity Statement” at [www.acumed.net/ifu](http://www.acumed.net/ifu) for a detailed description of the material composition of Acumed metal implants.

### SINGLE USE

- Implants are intended for single use only, as indicated on the label. This product has not been validated for multiple uses. Reuse of single use implants may increase the risk of failure and cross-contamination.
- Dispose of any unused implant that is contaminated with human blood or tissue in accordance with hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material. Do not process a contaminated implant.

### IMPORTANT

- This product is provided sterile and should not be re-cleaned.
- Physiological dimensions limit implant sizes. Select the type and size of implant that best meets the patient’s requirements for close adaptation and firm seating with adequate support.
- Implants are designed to hold fractured bones in place during bone healing at the site of fracture. Implants are not designed for excessive loadbearing. Improper selection or improper implantation of the device may increase the possibility of loosening, migration, neurovascular injury, or bony fracture.
- Do not bend the implant except as indicated in the corresponding Surgical Technique Guide. Repeated or excessive bending may weaken the implant and cause failure at a later time.
- Only combine implants when they are intended for that purpose.
- Protect implants against scratching and nicking to prevent stress concentrations, which can result in failure.
- Prevent unused implants from becoming soiled.
- The color of anodized implants may change over time due to processing. This color change does not affect the mechanical properties of the implants.

## INSTRUMENTS

### MATERIALS

The instruments are manufactured from various grades of titanium, stainless steel, anodized aluminum, medical grade plastic, and polymers evaluated for biocompatibility.

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#### Single Use Instruments

- Do not reuse single use instruments. Reuse or reprocessing of devices labeled as “single use” can result in decreased mechanical and clinical performance of the devices which may result in patient harm.

- Reuse or reprocessing of single use instruments may create a risk of contamination (e.g. due to the transmission of infection material from one patient to another), which may result in patient harm.
- Dispose of single use instruments after use on a single patient, during a single procedure.

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**STORAGE CONDITIONS****STORAGE OF PACKAGED STERILE PRODUCT**

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## Symbols Glossary

Symbol	Description	ISO 15223-1
 <small>www.acumed.net/ifu</small>	Consult the electronic instructions for use (eIFU) at <a href="http://www.acumed.net/ifu">www.acumed.net/ifu</a>	5.4.3
	Sterilized using ethylene oxide	5.2.3
	Double sterile barrier system	5.2.12
	Use-by date	5.1.4
	Catalogue number	5.1.6
	Batch code	5.1.5
	Manufacturer	5.1.1
	Date of manufacture	5.1.3
	Do not re-sterilize	5.2.6
	Do not re-use	5.4.2
	Do not use if package is damaged and consult instructions for use / do not use if the product sterile barrier system or its packaging is compromised.	5.2.8
	Caution: U.S. federal law restricts this device to sale by or on the order of a physician.	U.S. 21 CFR 801.109
	MR Conditional: An item with demonstrated safety in the MR environment within defined conditions. (ASTM F2503 7.3.2)	
	The reticle is a registered trademark of Acumed. It may appear alone or with the Acumed name.	

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# Instructions for use

## ACU-SINCH® REPAIR SYSTEM

**These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The implants and instruments are intended only for use by a licensed and trained health care professional.**

### DEVICE DESCRIPTION

The Acu-Sinch® Repair System implants and instruments are designed as an adjunct in clavicle fracture repair to provide additional fixation during the healing process when used in conjunction with a clavicle plate and screws from the Congruent Bone Plate System.

### INTENDED PURPOSE

The Acumed Acu-Sinch® Repair System implants and instruments are designed as an adjunct in clavicle fracture repair to provide additional fixation during the healing process when used in conjunction with a clavicle plate and screws from the Congruent Bone Plate System.

### INDICATIONS AND CLINICAL BENEFIT

The Acumed Acu-Sinch® Repair System is indicated to be used in conjunction with a clavicle plate of the Congruent Bone Plate System to provide fixation during the healing of clavicle fractures.

The clinical benefit of the Acu-Sinch® Repair System is providing unicortical fixation, which does not require drilling of the inferior cortex of the coracoid. This allows the Acu-Sinch® Repair System to address CC ligament disruptions associated with the clavicle fractures using a suture anchor in the coracoid tethered by suture to the clavicle plate.

### TARGET POPULATION AND INTENDED USERS

For both the Acumed Clavicle Plating System and the Acu-Sinch® Repair System, the target populations and users are:

Intended user: The system is intended to be used by suitably qualified surgeons knowledgeable about the device, with training and knowledge of the indications and techniques.

Population: The system is intended to be used on adults with sufficient quantity and quality of bone.

### CONTRAINDICATIONS

- Active or latent infection or marked inflammation of the treatment area.
- Sepsis.
- Osteoporosis, insufficient quantity or quality of bone and/or soft tissue.
- Patients with material sensitivity. If sensitivity is suspected, tests are to be performed prior to implantation.
- Patients who are unwilling or incapable of following post-operative care instructions.
- The device is not indicated as a sole treatment for chronic ligament and tendon injuries.

- The system may not be suitable for skeletally immature patients and must not disturb the growth plate.

## **WARNINGS and CAUTIONS**

### **Warnings**

- For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device.
- Treatment or implant may fail, including sudden failure, because of:
  - Loose fixation and/or loosening
  - Stress, including stress from inappropriate bending of the implant during surgery
  - Stress concentration
  - Stress of weight bearing, load bearing, or excessive activity
- Device breakage (including suture breakage) is more likely if the implant experiences increased loads due to delayed union, nonunion, or incomplete soft tissue healing.
- The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant, and the possibility soft tissue damage related to either surgical trauma (i.e. damage to neurovascular structures during preparation such as drilling or implant insertion) or purely due to the presence of an implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.
- The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity.
- Nerve or soft tissue damage may result from surgical trauma or the presence of an implant.
- Improper insertion of the device during implantation can increase the possibility of implant loosening and migration, neurovascular injury, or bony fracture.
- Devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects.
- Implants may cause distortion and/or block the view of anatomic structures on radiographic images.
- Drilling into or through bone may result in immediate or delayed bone fracture.

### **Precautions**

- Implants and instruments are intended only for professional use by a licensed physician.
- Use devices and instruments in accordance with the Instructions for Use.
- Do not use the sterile product past the use-by date. Refer to the device label.
- Do not use or re-sterilize an implant provided in sterile packaging if the package has been damaged. The sterility may be compromised, and the cleanliness of the implant may be uncertain. Report damaged packaging to your distributor or Acumed.
- Use of this system has not been evaluated in children or individuals that are not skeletally mature.
- Inspect all components preoperatively to ensure utility. Do not attempt a surgical procedure with faulty, damaged, or suspect instruments or implants. Alternate fixation methods should be available intraoperatively.
- Reprocessing and/or reuse of single use devices may result in infection/cross contamination, and/or sudden failure due to previous stresses.
- Use of implant components from different manufacturers has not been evaluated.
- Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.
- Inspect suture for signs of damage or fraying.

- Prior to use, inspect the product package for signs of tampering, water contamination, and seal integrity.
- The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures.
- Revision surgeries with implants are not uncommon.

## POTENTIAL ADVERSE EFFECTS

- Anesthesia-related problems, problems with positioning of the patient (e.g., nausea, vomiting, neurological impairments, etc.) thrombosis, embolism, infection, or injury to other critical structures such as blood vessels, excessive bleeding, etc.
- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.
- Swelling, abnormal scarring, impairment in musculoskeletal function, hardware prominence, migration, loosening, bending, or breakage of the implant or suture, malunion, nonunion or delayed union due to prolonged loading or excessive forces which may lead to implant failure and reoperation.
- Metal sensitivity, histological, allergic, or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu>.
- Injury to user.

**Important:** Any serious incident that has occurred in relation to the device should be reported to Acumed (through [customercomplaints@acumed.net](mailto:customercomplaints@acumed.net) or +1.888.627.9957) and the Competent Authority of the Member State in which the user and/or patient is established.

## SURGICAL TECHNIQUE GUIDE

Acumed offers one or more Surgical Technique Guides to promote the safe and effective use of this system. Consult our Surgical Technique Guides at <https://www.acumed.net/resources/documents/>. The Surgical Technique Guides also includes information about the installation and calibration of the device.

**Important:** Surgical Technique Guides may contain important safety information.

**Important:** The instruments and implants in this system are intended to be used by suitably trained and qualified surgeons in a hospital operating room setting. Before treatment, the surgeon is advised to read and fully understand all instructions and communicate to the patient any relevant medical information provided therein, including the use, limitations, risks (safety communications), and possible adverse effects of the proposed treatment.

Consult the most recent versions of the Instructions for Use and Surgical Technique Guides as they are subject to change. Contact Acumed or an authorized agent for any additional information.

## MRI SAFETY INFORMATION

Many Acumed implants have been evaluated for safety in the MR environment. Consult our publication "Acumed Implants in the MR Environment" at <http://www.acumed.net/ifu>.

## LIFETIME

- Once installed, implants are expected to provide fixation, fusion and physiological support and have an effective life during bone healing. The implants are biocompatible and may remain implanted at the discretion of the surgeon or patient.
- Multiple-use instruments have a lifetime that is affected by usage, handling, and processing. Assess multiple-use instruments for fitness during pre-sterilization inspection.

## STERILITY

- Implants and instruments are provided sterile.
- Sterility was achieved by using the ethylene oxide gas sterilization method to Sterility Assurance Level (SAL) of 10<sup>-6</sup>.
- Acumed does not recommend resterilization of sterile-packaged product.
- If sterile packaging is damaged, the incident must be reported to Acumed.
- The product must not be used, and must be returned to Acumed.

## IMPLANTS

### MATERIALS

- Implants: The metal implants are made out of implantable grade titanium.
- Suture: The suture is composed of Ultra High Molecular Weight Polyethylene (UHMWPE).
- Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu> for a detailed description of the material composition of Acumed metal implants.

### SINGLE USE

- Implants are intended for single use only, as indicated on the label. This product has not been validated for multiple uses. Reusing single use implants may increase the risk of failure and cross-contamination.
- Dispose of any unused implant that is contaminated with human blood or tissue in accordance with hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material. Do not process a contaminated implant.

### IMPORTANT

- This product is provided sterile, and should not be re-cleaned.
- Physiological dimensions limit implant sizes. Select the type and size of implant that best meets the patient's requirements for close adaptation and firm seating with adequate support.
- Implants are designed to hold fractured bones in place during bone healing at the site of fracture. Implants are not designed for excessive loadbearing. Improper selection or improper implantation of the device may increase the possibility of loosening, migration, neurovascular injury, or bony fracture.
- Do not bend the implant except as indicated in the corresponding Surgical Technique Guide. Repeated or excessive bending may weaken the implant and cause failure at a later time.
- Only combine implants when they are intended for that purpose.
- Protect implants against scratching and nicking to prevent stress concentrations, which can result in failure.
- Prevent unused implants from becoming soiled.

- The color of anodized implants may change over time due to processing. This color change does not affect the mechanical properties of the implants.

## INSTRUMENTS

### MATERIALS

The instruments are manufactured from various grades of titanium, stainless steel, anodized aluminum, and polymers evaluated for biocompatibility.

### SINGLE USE

Acu-sinch repair system only consists of sterile single-use instruments.

#### Single Use Instruments

- Do not reuse single use instruments. Reuse or reprocessing of devices labeled as "single use" can result in decreased mechanical and clinical performance of the devices which may result in patient harm.
- Reuse or reprocessing of single use instruments may create a risk of contamination (e.g. due to the transmission of infection material from one patient to another), which may result in patient harm.
- Dispose of single use instruments after use on a single patient, during a single procedure.

### IMPORTANT

- Protect instruments against scratching and nicking to prevent stress concentrations, which can lead to instrument failure.

## PROCESSING

- Not Applicable as Acu-sinch Repair System is provided sterile for single use only.

## STORAGE CONDITIONS

### STORAGE OF PACKAGED STERILE PRODUCT

Devices should be stored in an area that provides protection from dust, pests, and temperature/ humidity extremes.

### SAFE DISPOSAL

Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance Report for the implants may be obtained from the European dataset on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Use the following search terms: 0806378BUDI26Q4

**Important:** Any serious incident that has occurred in relation to the device should be reported to Acumed and the competent authority of the Member State in which the user and/or patient is established. To report a complaint about this product, call Customer Service at +1.888.627.9957 or send an email to [customercomplaints@acumed.net](mailto:customercomplaints@acumed.net).

## Symbols Glossary

Symbol	Description	ISO 15223-1
	Consult the electronic instructions for use (eIFU) at <a href="http://www.acumed.net/ifu">www.acumed.net/ifu</a>	5.4.3
	Sterilized using ethylene oxide	5.2.3
	Double sterile barrier system	5.2.12
	Use-by date	5.1.4
	Catalogue number	5.1.6
	Batch code	5.1.5
	Authorized representative in the European Community / European Union	5.1.2
	Medical device	5.7.7
	Manufacturer	5.1.1
	Date of manufacture	5.1.3
	Do not re-sterilize	5.2.6
	Do not re-use	5.4.2
	Unique Device Identifier: Indicates a carrier that contains Unique Device Identifier information.	5.7.10
	Do not use if package is damaged and consult instructions for use / do not use if the product sterile barrier system or its packaging is compromised.	5.2.8
	MR Conditional: An item with demonstrated safety in the MR environment within defined conditions. (ASTM F2503 7.3.2)	
	The reticle is a registered trademark of Acumed. It may appear alone or with the Acumed name.	
	CE marking of conformity, Article 17 of EU Directive 93/42/EEC or Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment.	

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Fax: +1.503.520.9618

[www.acumed.net](http://www.acumed.net)

*These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.*

PKGI-60-rev T | Effective 06-2025 | © 2025 Acumed® LLC

Acu-sinch® Repair System is a registered trademark of Acumed, LLC.

### Subsidiaries:

#### Acumed Ltd

Huebner House  
The Fairground  
Andover  
Hampshire UK SP11 0QN  
Tel: +44 1264 774450

#### Acumed Iberica

C/ Proción 1  
Edificio Oficor  
28023 Madrid – Spain  
Tel: +34 913516357

#### Acumed Beijing

Room A1206, Horizon International Tower  
No. 6, Zhichun Road  
Haidian District  
100088 Beijing, China  
Tel: +86 10 82001303

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#### Acumed GmbH