

Case Study

The Use of Closed Reduction Intramedullary Fixation for an Oblique Fifth Metacarpal Neck Fracture



Kathleen McKeon, MD

Andrews Sports Medicine & Orthopaedic Center
Birmingham, AL

Dr. McKeon is a board certified, fellowship trained orthopedic surgeon who specializes in hand and wrist surgery. Dr. McKeon received the Subspecialty Certificate in Orthopedic Surgery of the Hand after earning her medical degree from Washington University School of Medicine in St. Louis where she also completed her residency. She was fellowship trained in hand surgery at Vanderbilt University. She treats carpal tunnel syndrome, trigger finger, arthritis, finger/hand/wrist fractures just to name a few.

Case Presentation

A 56-year-old, right-hand dominant male patient presented to clinic the day after injuring his left hand. X-rays showed an oblique, displaced fracture of the left fifth metacarpal neck. Nonoperative and operative treatment options were discussed with the patient. Given the amount of shortening and radial translation present at the fracture site, reduction and operative fixation were recommended.

Pre-op Plan

Intramedullary (IM) fixation, without exposed K-wires, allows for earlier return of motion and decreased risk of infection. Fixation with a rigid intramedullary nail allows for fewer restrictions in the immediate postoperative period. Dr. McKeon had successfully used the INnate™ nail for midshaft metacarpal fractures, and wanted to extend those benefits to this patient as well. Having never used the INnate nail for this fracture pattern, she had other fixation options available as backup.

Operative Findings and Approach

Dr. McKeon made two 2 mm incisions to allow placement of a percutaneous bone reduction clamp to achieve and maintain fracture reduction. The guide wire was placed percutaneously from the metacarpal head into the metacarpal shaft. Once the guide wire was in good position, she made a 2 mm incision around the wire. She then used the cannulated drill to drill over the guide wire and threaded the cannulated INnate nail until the head was just beneath the articular cartilage, to achieve distal purchase just below subchondral bone. Proximal purchase was achieved at the isthmic level within the IM canal. The reduction clamp was removed, and she verified that the fracture was stable. A soft dressing was placed, and the patient was allowed to return to activities as tolerated.

Preoperative



Postoperative



Follow-up

On the eighth day post-op, the patient had minimal pain. He had full, active extension of the small finger. He could actively flex the small finger to within 1 cm of the distal palmar crease. By four weeks post-op the patient had returned to full activity, including push-ups. He never required formal therapy.

Discussion

Dr. McKeon has been pleased with the INnate nail. It is her first choice for treating many different patterns of metacarpal fractures, including distal oblique. It is minimally invasive, does not compress the fracture, and does not require casting, allowing immediate mobilization and rapid recovery.

**SPF70-30-A**

Effective: 2025/07

© 2025 Acumed® LLC

www.acumed.net

Acumed Oregon Campus
5885 NE Cornelius Pass Road
Hillsboro, OR 97124
+1.888.627.9957

Acumed Texas Campus
3885 Arapaho Road
Addison, TX 75001
+1.800.456.7779

Acumed Iberica Campus
C. Proción, 1
Edificio Oficor
28023 Madrid, Spain
+34.913.51.63.57

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located. Specific questions physicians may have about the availability and use of the products described on these materials should be directed to their particular authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

ExsoMed™ is a wholly owned subsidiary of Acumed LLC.

Acumed®, Exsomed™, and INnate™ are registered trademarks of Acumed LLC.