Description

The **OSTEOEMD Temporary Condylar Attachment System** is a solid condylar head which attaches with fastening screws to an OsteoMed Locking Fracture Reconstruction Plate. The instruments include screwdrivers to facilitate the placement of the temporary condylar attachment.

Material

The **OSTEOEMD Temporary Condylar Attachment (TCA)** is made from Titanium Alloy (ASTM-F-136). The OsteoMed Locking Fracture Reconstruction Plates are made from Titanium (ASTM-F-67) and the fastening screws are made from Titanium Alloy (ASTM-F-136). The instrumentation is made from various grades of stainless steel, anodized aluminum, and/or medical grade plastic.

Clinical Indications

The **OSTEOEMD Temporary Condylar Attachment System** is indicated for temporary reconstruction of the mandibular condyle in patients who have undergone resective procedures to remove malignant or benign tumors requiring the removal of the mandibular condyle. **This device is intended for implantation for no more than one (1) year.** This device is not intended for permanent implantation, for patients with TMJ or traumatic injuries, or for treatment of temporomandibular joint disease (TMD). The **OSTEOEMD Temporary Condylar Attachment** implants are intended for single patient use only. OsteoMed single use devices/instruments cannot be reused and/or reprocessed. The product has been labeled as single use only for patient safety. The design of the device and the intricacies of the surfaces may not facilitate cleaning and sterilization after contact with body tissues or fluids, so there is an increased risk of contamination if reused. This may lead to potential risks of cross infection/contamination associated with using inadequately cleaned and sterilized devices. For cutting instruments the cutting efficacy may be reduced requiring the surgeon to use increased force that might cause patient harm and/or increase the risk of thermal necrosis. Because these products have not been validated for multiple uses, OsteoMed cannot guarantee the safety and effectiveness of the device if it is used on more than one patient.

Contraindications

Use of the **OSTEOEMD Temporary Condylar Attachment System** is contraindicated in cases of active or suspected infection or in patients who are immunocompromised; in patients previously sensitized to titanium; or in patients with certain metabolic diseases; in patients exhibiting disorders which would cause the patient to ignore the limitations of rigid fixation.

Warnings

- Use of the **OSTEOEMD Temporary Condylar Attachment** with appropriate reconstruction plate and screws in areas of high functional stresses, greater than 88N (i.e. non-reconstructed patients without a restricted diet), may lead to implant system fracture and/or failure.
- Plates and screws, wires, or other appliances of dissimilar metals should not be used together in or near the implant site.
- OsteoMed recommends the use of OsteoMed products in a sterile environment.

MRI Status:

The **OSTEOEMD Temporary Condylar Attachment** has not been evaluated for safety and compatibility in the MR environment. It has not been tested for radio frequency (RF) heating, migration due to magnetically induced displacement, or image artifact in the MR environment. The safety of the implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Maintaining Device Effectiveness

- The surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation products and techniques.
- The surgeon must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
- The **OSTEOEMD Temporary Condylar Attachment System** plates and screws are not intended to endure excessive abnormal functional stresses (greater than 88N), such as those associated with an unrestricted diet.
- The **OSTEOEMD Temporary Condylar Attachment System** is intended for temporary fixation.
- All OsteoMed plates, screws, and instrumentation may be required for each surgery. Failure to use dedicated, unique OsteoMed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Carefully inspect the OsteoMed implants prior to use. Inspect the instruments before and after each procedure to assure they are in proper operating condition. Instruments which are faulty, damaged, or suspect should not be used. They should be replaced or sent to OsteoMed for disposition and repair.
- OsteoMed recommends the use of OsteoMed products in a sterile environment.

Instructions for Use: This is a TRACKED MEDICAL DEVICE.

- Contour OsteoMed Locking Reconstruction Plate to the patient's anatomy.
- Ablate the desired portion of the mandible.
- Mate the OsteoMed Temporary Condylar Attachment (TCA) Implant without screws to the Locking Reconstruction Plate and check for proper position of the condyle in relationship to the patients Temporal Mandibular Joint.
- Assemble the TCA Implant using the mounting screws provided and an OsteoMed Ø2.4mm driver stem. Check for proper position relative to the TMJ. Use a minimum of 3 screws to attach the TCA Implant to the Locking Reconstruction Plate.
- Attach the Locking Reconstruction Plate and TCA Implant to the patient's mandible using OsteoMed Ø2.4mm locking screws.
- Check Locking Reconstruction Plate/TCA Implant function against joint.
- The TCA Implant is Temporary and is intended for implantation for no more than one (1) year.** Standard OsteoMed Ø2.4mm drivers can be used to remove the mounting screws and TCA Implant from the Locking Reconstruction Plate.

Cleaning

- Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- OsteoMed recommends the following cleaning and sterilization instructions for Instrumentation:
 - Clean all instruments thoroughly using mild detergent, soft brush and warm water. Ensure that dried blood, bone chips and other deposits are removed from the instruments and sterilization tray.
 - Thoroughly rinse all instruments and the sterilization tray with water.
 - Arrange all the instruments in the sterilization case and ensure that the lid is in place and properly closed.
 - Steam Autoclave per the following Sterilization Instructions.

Sterility

- Product is supplied **NON-STERILE**.
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- The user facility must clean and disinfect devices prior to sterilization per standard hospital procedures.
- Non-sterile devices are sterilizable by steam sterilization (autoclaving). For sterilization of OsteoMed implant systems, the following parameters should be used.

Pre-Vacuum Steam Sterilization for **TCA Implant and fastening Screws:**

Temperature: **270°F (132°C)**

Time: **6 minutes**

Dry Time: **0 minutes**

Configuration: **Pouched**

Do not exceed 275°F (135°C), to avoid compromising functions of polymeric instrumentation.

Storage

The Temporary Condylar Attachment System is to be stored at a controlled room temperature between 15°C to 30°C (59°F to 86°F)

Caution

- Federal (United States) law restricts this device for sale by or on the order of a medical practitioner licensed to do so.**
- Do not attempt a surgical procedure with faulty, damaged, or suspect OsteoMed instruments or implants. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.**

OSTEOEMD
3885 Arapaho Road
Addison, Texas 75001 USA
Customer Service: 800/456-7779
Outside USA: 972/677-4600

Symbols and Definitions			
	Single Use Only	REF	Catalogue Number
	Batch Code (Lot Number)	SN	Serial Number
	Date of Manufacture		Manufacturer
	Consult Instructions for Use		Federal Law (U.S.A) Restricts this device to sale by or on the order of a physician.
	Attention, See Instructions for Use Caution, Consult Accompanying Documents		
	Storage Temperature		Printed in U.S.A. OSTEOEMD July 2025

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