Part No. **030-1503**
Rev. **U**
EN

OSTEOMED
1st MPJ Hemi Implant System
ESSENTIAL PRESCRIBING INFORMATION

Device Description

The **OSTEOMED 1st MPJ Hemi Implant System** is a one piece implant to supplement the first metatarsophalangeal joint. The Hemi implant is a phalangeal base implant used in arthroplasty of the first MPJ and is available in four sizes. The instrumentation for the system includes trials, sizers, punch and impactor for placement of device.

Material

The **OSTEOMED 1st MPJ Hemi Implant System** is comprised of implants made from cobalt chromium (ASTM F-799). The instrumentation is made of stainless steel, anodized aluminum and/or medical grade plastic.

Indications for Use

The **OSTEOMED 1st MPJ Hemi Implant System** is intended for use in the treatment of patients with inflammatory arthritis in the first metatarsophalangeal joint in the presence of good bone stock and integrity of the first metatarsal head, along with the following clinical conditions: hallux valgus, hallux rigidus and an unstable or painful metatarsophalangeal joint. **OSTEOMED 1st MPJ Hemi Implant System** implants are intended for single patient use only. OsteoMed single use devices/instruments cannot be reused and/or reprocessed. The product has been labeled as single use only for patient safety. The design of the device and the intricacies of the surfaces may not facilitate cleaning and sterilization after contact with body tissues or fluids, so there is an increased risk of contamination if reused. This may lead to potential risks of cross infection/contamination associated with using inadequately cleaned and sterilized devices. For cutting instruments the cutting efficacy may be reduced requiring the surgeon to use increased force that might cause patient harm and/or increase the risk of thermal necrosis. Because these products have not been validated for multiple uses, OsteoMed cannot guarantee the safety and effectiveness of the device if it is used on more than one patient.

Contraindications

Use of the **OSTEOMED 1st MPJ Hemi Implant System** is contraindicated in cases of active or suspected infection or in patients who are immunocompromised; in patients previously sensitized to the implant material; in patients with certain metabolic diseases; or in patients exhibiting disorders which would cause the patient to ignore the limitations of artificial joint replacement. This system is further contraindicated in serious systemic diseases. The **OSTEOMED 1st MPJ Hemi Implant System** should not be used in cases where the remaining bone is too diminished to provide adequate width or height to surround the implant. Implant failure may occur in cases where there is insufficient available bone or poor bone quality or in patients with blood supply limitations. Use of the implant is further contraindicated in cases of significant angular or biomechanical deformities or in cases which would subject the implants to excessive stress or wear due to patient weight and/or activity level.

Warnings

- The **OSTEOMED 1st MPJ Hemi Implant System** implants are not intended to endure excessive functional stresses. Such stresses may lead to accelerated implant wear and/or failure.
- Protect Implant surfaces from scratching or other surface damage. Such damage may lead to accelerated implant failure. Avoid direct impact on implants with hammers or other instruments.

MRI Status:

The **OSTEOMED 1st MPJ Hemi Implant System** has not been evaluated for safety and compatibility in the MR environment. It has not been tested for radio frequency (RF) heating, migration due to magnetically induced displacement, or image artifact in the MR environment. The safety of the implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Maintaining Device Effectiveness

- The Surgeon should have specific training, experience and thorough familiarity with the use of joint replacement implant devices. Surgeon should familiarize him/herself with the surgical instruction manual.
- All OsteoMed implants and instrumentation may be required for each surgery. Failure to use dedicated, unique OsteoMed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Carefully inspect the OsteoMed implants prior to use. Inspect the instruments before use to assure they are in proper operating condition. Instruments which are faulty, damaged, or suspect should not be used. They should be replaced or sent to OsteoMed for disposition.
- OsteoMed products are to be used in a sterile environment.

Individualization of Treatment

The Surgeon must exercise reasonable judgment when selecting the size of implant to use for specific indications and should consider the quality of bone, bone type, functional loads, soft tissue constraints and biomechanic implications of the procedure as well as the patient's post-operative goals and compliance.

Patient Information

Patients should be informed and cautioned in writing regarding the limitations of artificial joint replacement.

Storage

- Sterile packaged implants should be stored at controlled room temperature out of direct sunlight.
- Product package should be inspected prior to use for signs of damage or tampering.
- Sterile implants have a shelf life of 4 years.

Cleaning (Instruments)

- Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- OSTEOMED** recommends the following cleaning and sterilization instructions for Instrumentation:
 - Clean all instruments thoroughly using mild detergent, soft brush and warm water. Ensure that dried blood, bone chips and other deposits are removed from the instruments and sterilization tray.
 - Thoroughly rinse all instruments and the sterilization tray with water.
 - Arrange all the instruments in the sterilization case and ensure that the lid is in place and properly closed.
 - Steam Autoclave per the following Sterilization Instructions.

Sterility


- OSTEOMED 1st MPJ Hemi Implant System** implants are packed **STERILE** (Gamma Sterilized). **DO NOT USE IF STERILE PACKAGE IS DAMAGED.**
- DO NOT USE AFTER EXPIRATION DATE.**
- Instruments are supplied **NON-STERILE** and **MUST** be sterilized prior to use.
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- The user facility must clean and disinfect devices prior to sterilization per standard hospital procedures.
- Non-sterile instruments are sterilizable by steam sterilization (autoclaving). For sterilization of **OSTEOMED 1st MPJ Hemi Implant system** instruments, the following parameters should be used:
 - Pre-Vacuum Steam Sterilization
 - Temperature: 270°F (132°C)
 - Time: 15 minutes
 - Dry Time: 10 minutes
 - Configuration: Individually wrap the case in two layers of 1-ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques with a surgical towel placed between the wraps and the bottom of the case.

Do not exceed 275°F (135°C), to avoid compromising functions of polymeric instrumentation.

Note: Biochemical indicator of *G. stearothermophilus* was used in sterilization validation.

Caution

- Federal Law (USA) restricts this device to sale by or on the order of a medical practitioner licensed to do so.**
- Do not attempt a surgical procedure with faulty, damaged, or suspect OsteoMed instruments or implants. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.**

Read manual prior to use

OSTEOMED
3885 Arapaho Road
Addison, Texas 75001 USA
Customer Service: 1-800-456-7779
Outside USA: 1-972-677-4600

Symbols and Definitions			
	Single Use Only	REF	Catalogue Number
	Use By (Date)		Do not use if sterile package is damaged
	Batch Code (Lot Number)		Consult Instructions for Use
	Date of Manufacture		Manufacturer
	Attention, See Instructions for Use Caution, Consult Accompanying Documents		
	Federal Law (U.S.A) Restricts this device to sale by or on the order of a physician.		Sterile, Method of Sterilization Using Irradiation
	Keep Away from Sunlight		Keep Dry

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