Part No. 030-1552
Rev. AF
EN

OSTEOMED
ICON
Modular Locking Fixation System
CFx™ Rigid Fixation Systems
M3/M4 /Craniomaxillofacial Trauma Fixation Systems

Product Information and Instructions for Use

Description

The **OSTEOMED** **ICON** **Modular Locking Fixation System** is comprised of plates, screws, IMF/MMF screw and instrumentation for fracture fixation in maxillofacial trauma reconstruction, mandibular reconstruction and orthognathic reconstruction. The **OSTEOMED CFx™ Rigid Fixation Systems** are comprised of plates, screws and instrumentation used for fixation of craniofacial, maxillofacial and mandibular fractures and reconstructions. Notes as **OSTEOMED Fixation Systems** throughout this document

These systems feature 1.2mm, 1.6mm, 2.0mm and 2.4mm plates ranging from 0.25mm to 2.8mm thick, Locking plates ranging from 1.0mm to 2.8mm thick, 1.6mm Low Profile Midface Plates (LPM) 0.5mm thick, 1.6mm LPM Mesh ranging 0.5mm to 0.7mm thick, 1.2 mm to 2.7mm diameter standard screws in lengths from 2.0mm to 22.0mm, 1.6mm Auto-Drive® screws in lengths from 3.5mm to 8.0mm, 2.4 diameter IMF (Intermaxillary Fixation) screws in lengths from 10mm to 20mm and 2.0mm diameter IMF Auto-Drive® screws in 8mm and 11mm lengths, 2.0mm diameter Auto-Drive® Lag in lengths of 11mm and 14mm, 2.0mm diameter Locking Screws in lengths of 6.0mm to 18.0mm, and 2.4mm diameter Locking Screws in lengths of 6.0mm to 22.0mm. The IMF screws have holes on the head to allow wire through to fixate fracture segments in the oral cavity. The system also includes Low Profile 1.6mm Auto-Drive® screws in lengths of 3.5mm to 6.0mm. Angled Locking Plates and screws allow up to 20 degrees of angulations within screw placement. The Mandible 2.0mm Angled Locking System features 2.0mm Angled Locking Standard Screws with lengths 4mm to 18mm, 2.0mm Angled Locking Auto-Drive® Screws with lengths of 5mm to 8mm, and 2.0mm Angled Locking Safety Screws with lengths ranging from 4mm to 18mm. The mandible plates are 1.0mm, 1.5mm, 2.0mm, 2.5, and 2.8mm in thickness. The Mandible 2.4mm Angled Locking System features 2.4mm Angled Locking Standard Screws with lengths of 6mm to 22mm, 2.4mm Angled Locking Auto-Drive® Screws with lengths of 6mm to 8mm, and 2.7mm Angled Locking Safety Screws with lengths ranging from 6mm to 22mm. The fracture plates range in thicknesses of 1.5mm to 2.0mm and the reconstruction plates are 2.5mm to 2.8mm thick. The OsteoMed Low Profile Orbital Floor Plates offer hole sizes which are able to accept 1.2mm and 1.6mm standard and auto-drive bone screws and 1.5mm and 1.9mm safety screws and are offered in thicknesses ranging from .2mm to 5mm.

The instruments include drill bits, plate bending forceps, plate holding forceps, plate cutters, cannulae, taps, countersinks, plate bending pliers, mesh benders, plate cutters, drill guides and screwdrivers to facilitate the placement of screws and modification of plates. The system includes wire cutters and twisters to facilitate the ligation and placement of the wire in the IMF/MMF Screws.

The **OSTEOMED Systems** implants and drills are intended for **single use only**. System instruments are reusable. OsteoMed single use devices/instruments cannot be reused and/or reprocessed. The product has been labeled as single use only for patient safety. The design of the device and the intricacies of the surfaces may not facilitate cleaning and sterilization after contact with body tissues or fluids, so there is an increased risk of contamination if reused. This may lead to potential risks of cross infection/contamination associated with using inadequately cleaned and sterilized devices. For cutting instruments the cutting efficacy may be reduced requiring the surgeon to use increased force that might cause patient harm and/or increase the risk of thermal necrosis. Because these products have not been validated for multiple use, OsteoMed cannot guarantee the safety and effectiveness of the device if it is used on more than one patient.

Material

The screws are made from Titanium Alloy (ASTM F-136). The plates are made from Titanium Alloy (ASTM F-136) or commercially pure Titanium (ASTM F-67). The instrumentation is made from various grades of stainless steel, anodized aluminum, and/or medical grade plastic or acetal copolymer.

Clinical Indications: Maxillofacial Indications Only

- The **OSTEOMED** **ICON** **Modular Locking Fixation System** is comprised of plates, screws and instrumentation for fracture fixation in maxillofacial trauma reconstruction, mandibular reconstruction and orthognathic reconstruction.
- The **OSTEOMED MMF / IMF Screws** are indicated for temporary ligature and wire lock fixation for temporary constriction and stabilization of fractured bone segments in the oral cavity in conjunction with primary fixation devices.
- The **OSTEOMED 2.0 Locking Plate System** is indicated for oral, maxillofacial surgery; trauma; reconstructive surgery; and orthognathic surgery (surgical correction of dentofacial deformities). Screws and Plates are intended for single patient use only.
- The **OSTEOMED** **ICON** **Modular Locking Fixation System (2.8mm thickness plates)** is indicated for fracture fixation in maxillofacial trauma reconstruction, mandibular reconstruction and orthognathic reconstruction.
Note: The ICON 2.8mm thickness plates are not to be placed on craniofacial fractures (PN: 231-2801, 231-2802, 231-2811, 231-2812, 231-2813, 231-2814, 231-2821, 231-2822)
- The **OSTEOMED Low Profile Orbital Floor Plate System** is intended to be used in the reconstruction of the floor and/or medial wall of the orbit. The OsteoMed Low Profile Plate System is indicated for the reconstructive treatment of orbital floor and/or medical wall trauma or bone excision in adult patients.

Clinical Indications: Cranial Maxillofacial Indications Only

- The **OSTEOMED CFx™ Rigid Fixation Systems** are indicated for fracture fixation during cranial and facial trauma reconstructions, orthognathic reconstruction, mandibular reconstruction and surgery involving osteotomies and trauma.
- The **OSTEOMED Angulated Locking Fixation System** is indicated for mandibular fracture fixation in cranio-maxillofacial trauma reconstruction, mandibular reconstruction and orthognathic reconstruction.

Contraindications

- Use of the **OSTEOMED** *Fixation Systems* are contraindicated in cases of active or suspected infection or in patients who are immunocompromised; in patients previously sensitized to titanium or stainless steel; or in patients with certain metabolic diseases.
- It is further contraindicated in patients exhibiting disorders which would cause the patient to ignore the limitations of rigid fixation plate and screw implants and/or in patients where there is insufficient bone or poor bone quality.
- The **OSTEOMED** *MMF/ IMF Screws* are also contraindicated in patients with disorders which could prevent the patient from following the limitations of temporary ligature and wire fixation.
- The use of 2.4mm Angled Locking Auto-Drive® screws are contraindicated in bi-cortical bone placement on the mandible.

Warnings

- Manipulation of maxillofacial bones or tissues can potentially result in cardiac complications such as ectopic beats, atrioventricular block, bradycardia, syncope, vomiting, and asystole. These are known complications during maxillofacial surgery and are not specifically related to any device.
- Use of an undersized plate or screw in areas of high functional stresses may lead to implant fracture and failure.
- Plates and screws, wires, or other appliances of dissimilar metals should not be used together in or near the implant site.
- Multiple bending may weaken the plate and could result in implant fracture and failure.
- Use of screws in high density bone may lead to implant fracture or failure upon insertion.
- It is recommended to remove any fractured implants from patients during surgery. If unable to remove, notify patient.
- Use of excessive torque during insertion of screws may lead to implant failure.
- OSTEOMED** *MMF/ IMF Screws* are not intended to endure excessive abnormal function stresses. Use of an undersized wire or screw in areas of high abnormal functional stresses may lead to wire or screw fracture and failure.
- OSTEOMED** *MMF/ IMF Screws* and wire are intended for temporary fixation and are to be implanted for a maximum of four weeks.
- Selection of the screw size must be carefully considered by the operating surgeon and should take into consideration the quality of bone, bone type functional loads exerted on the bone(s) and post-operative patient compliance. Screw site selection shall also consider the avoidance of crucial anatomic structures, which should be verified through a radiographic technique.
- Cannula is never to be used as a cheek retractor when drilling.
- When placing additional screws, ensure that subsequent screw placement does not interfere with the other screws.
- When using the drill guide, do not apply a side load on the drill. This may result in friction, which may generate a thermal burn. Axial loading should always be used.
- When using the transbuccal approach, ensure cheek retractor is used to protect soft tissue.

MRl Status:

The **OSTEOMED** *Fixation Systems* has not been evaluated for safety and compatibility in the MR environment. It has not been tested for radio frequency (RF) heating, migration due to magnetically induced displacement, or image artifact in the MR environment. The safety of the implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Maintaining Device Effectiveness

- The surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation and IMF/MMF products and techniques.
- The surgeon must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
- The **OSTEOMED** *Fixation Systems* plates and screws are not intended to endure excessive abnormal functional stresses.
- When loading a screw onto driver, apply a perpendicular force to engage screw cruciform with driver. The surgeon should avoid multiple insertions of driver into the same screw to maintain self-retention feature of screw to driver.
- The **OSTEOMED** *Fixation Systems* are intended for temporary fixation only until osteogenesis occurs.
- Dedicated OsteoMed plates, screws, and instrumentation are required for each surgery. Failure to use dedicated, unique OsteoMed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Carefully inspect the **OSTEOMED** implants prior to use. Inspect the instruments before and after each procedure to assure they are in proper operating condition. Instruments which are faulty, damaged, or suspect should not be used. They should be replaced or sent to **OSTEOMED** for disposition and repair.
- When placing more than one screw, ensure that subsequent screw placement does not interfere with other screws. Insert the second screw on the opposite side of the fracture or osteotomy site, and then all remaining screws, following the outlined procedures.
- OSTEOMED** recommends the use of OsteoMed products in a sterile environment.
- Drill using the appropriate pilot drill. Note: Speed and torque parameters must be in accordance to the Power System Instruction for Use. Use irrigation when pilot drilling.
NOTE: Drill Color Bands Denotes Size: 1.2mm drills are Yellow, 1.6mm drills are White, 2.0mm drills are Blue and 2.4mm drills are Green.
- After cutting plate, utilize diamond file on the plate cutter handle (220-0584) to remove sharp edges on the plate.
- When using Plate Benders 220-0548 and 220-0529, note they can only be used with Reconstruction plates.
- Depth gauge marking tolerance: ±0.25mm

Instructions for Use, Auto-Drive® screws

The Auto-Drive® screws are self drilling and can be inserted in one step. Insert the screw in a TaperLock™ screwdriver and drive into the bone at a 90° angle using moderate pressure until the head is flush with the surface of the bone/plate. Higher torque may be required to fully engage the threads than when using a normal screw with a pilot drill.

Note: In high density bone, pilot drilling may be necessary.

Instructions for Use, Locking Screw System for Mandible Applications.

- Expose and reduce fracture: After completing the preoperative plan, expose the fracture or osteotomy site. For trauma, reduce the fracture as required.
- Select and prepare implants: Select the appropriate template / plate depending on the indication.
Note: For 1.6mm, 2.0mm and 2.4mm plates, it is recommended to allow at least two screws per bone segment when possible. For Reconstruction plates, allow for at least 3 screws per bone segment. Orient the plate so the topside is facing out. Determine the appropriate screw type, locking or non-locking.
Note: Angulation and screw selection are contingent on patient anatomy and surgeon preference.
- Place screw inserts into holes selected to receive screws. Contour the template to match the anatomy. An exact match is not required when using locking screws, as plate stability is not dependent on plate-to-bone contact when screws are locked. Cut and contour plate to match template form. Plates can be cut with plate cutters.
- Position the plate: Place the plate over the fracture or osteotomy site. Use the plate holding forceps to secure the plate to the bone, if desired.
- Drill the first hole: Select the threaded drill guide and insert it into the first plate hole nearest the fracture or osteotomy site. Rotate the drill guide clockwise to engage the threads in the plate. Drill using the appropriate drill bit.
Note: When using Auto-Drive™ Locking screws, which are self-drilling, higher torque may be required to seat the screw. Pilot drilling may be necessary in dense bone.
- Measure screw length: Remove guide and use depth gauge to measure hole depth to determine appropriate screw length.
- Insert the screw: Insert the proper length locking or non-locking screw through the plate and tighten until secure. Locking screws can be inserted at a +/- 10 degree angle.
Note: When using Auto-Drive™ Locking screws, which are self-drilling, higher torque may be required to seat the screw. Pilot drilling may be necessary in dense bone. 8. Drill and place the remaining screws: Insert the second screw on the opposite side of the fracture or osteotomy site, and then all remaining screws, following the previously outlined procedure. Securely tighten all screws unless mandible resection is to follow.
- Note:** Speed and torque parameters must be in accordance to the Power System Instructions for Use. Use irrigation when pilot drilling.

Attention: Instructions for Use are in the Surgical Technique Guide (030-1559) or OsteoMed MFx™ Surgical Technique Guide 030-1623 (for Angulated Locking) or PI 030-1567.

Instructions for Use, IMF/MMF Screws (Screws are not to be inserted on a plate)

- Screws are placed through the mucosa without making an incision, taking care to avoid the roots of the teeth. Choose screw sites away from the root apices usually medial or distal to the cuspid tooth.
- Prepare bone with a pilot drill that is longer than the intended screw when using 2.4mm IMF/MMF screws, taking care to avoid critical structures.
Note: 2.0mm AutoDrive IMF/MMF Screws are self-drilling, no pilot drilling is necessary.
- Note:** In high density bone, pilot drilling may be necessary for AutoDrive screws.
- Select appropriate IMF/MMF screw length.
- Engage driver tip lightly into screw cruciform and apply moderate pressure. Vertically retract driver and screw from organizer and verify length with length gauge. Insert the screw into the pilot hole and drive the screw to the appropriate depth, leaving the wire-passing hole exposed. Do not over torque or bottom out the screw.
Note: The Auto-Drive screws are self-drilling and can be inserted in one step. Insert the screw in a TaperLock™ screwdriver and drive into the bone at a 90 angle using moderate pressure until the head is flush with the surface of the bone/plate. Higher torque may be required to fully engage the threads than when using a normal screw with a pilot drill.
- For secondary screw ensure placement in the mandible is 5mm inferior and medial or lateral to the canine tooth roots.
- A minimum of three pairs of IMF/MMF screws are recommended to ensure adequate stability. A pair consists of one screw in the mandible and an opposing screw in the maxilla.
- Wire screws using 24 gauge stainless steel wire (207-0120) through exposed wire passing holes into maxillary and opposing mandibular IMF/MMF screw heads in a vertical and "X" pattern. Tighten only enough to provide provisional fixation.
- Establish occlusion and tighten wire fully.

Instructions for Use, Orbital Floor Plates

- Expose the fracture site and retract soft tissue
- Select the orbital floor plate to be used based upon fracture location and patient anatomy.
- If the orbital floor plate requires modification by cutting, plate/mesh cutter (P/N 220 0077) may be used to trim unneeded material.
- Use deburring file (P/N 220 0584) if desired to remove any sharp edges.
If contouring is required prior to implantation, plate benders (P/Ns 220 0046 & 220 0049) may be used to pre bend the plate as needed. (Note: To prevent weakening of the plate, excessive and repeated bending should be avoided.)
- Insert and position plate with plate holding forceps (P/N 220 0027).

- Fixate plate to bone using OsteoMed's patented Auto Drive™ self drilling screws, stop insertion of the screw, once the screw is fully seated on the plate.

Cleaning

- Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- OSTEOMED** recommends the following cleaning and sterilization instructions for instrumentation:
 - Rinse the articles to be cleaned under running cool tap water (<40°C) to remove visible soil until visibly clean. Prepare an enzymatic detergent such as Enzol® or Klenzyme according to the manufacturer's recommendations at 1oz/gal of lukewarm tap water.2. Using a clean soft cloth that has been soaked in the detergent, wipe the entire article. Actuate the articles while immersed in the solution to ensure complete penetration of cleaning solution. Immerse the article in the detergent. Flush with 50mL all hard to reach areas using a syringe with the prepared detergent. Allow the article to soak for 15 minutes.
 - Using a soft bristled brush, brush the entire surface of the article paying close attention to all hard to clean areas such as cracks, crevices and matted surfaces. Actuate the article while brushing to ensure the reach of all areas. Use a syringe to flush any lumen or matted surface with at least 50mL of prepared detergent.
 - Brush any lumen a minimum of 5 times from both ends using a lumen brush or similar.
 - Rinse the articles in RO/DI (reverse-osmosis/deionized) water until all visible evidence of detergent is removed. Flush any lumens or matted surfaces with the RO/DI water. Once all evidence of detergent is removed continue to rinse for an additional 30 seconds.
 - Drain excess water from the article and dry using a clean soft cloth or filtered pressurized air.
 - Prepare a non-enzymatic detergent such as Renu-Klenz per manufacturer's recommendations at ¼ oz/gal of lukewarm tap water.
 - Fully immerse the articles in the detergent and flush all lumens with at least 50mL of the detergent.
 - Allow the articles to soak for a minimum of 15 minutes.
 - Thoroughly brush the exterior of the articles using a soft bristled brush. Flush all lumens and matted surfaces a minimum of 5 times.
 - Prepare a fresh cleaning detergent such as Renu-Klenz in an ultra sonic cleaner. Fully immerse the articles in the solution and allow the articles to sonicate for a minimum of 10 minutes.
 - Rinse the articles in running reverse osmosis / deionized (RO/DI) water until all visible evidence of detergent is removed. Flush any lumens or matted surfaces with the RO/DI water. Once all evidence of detergent is removed, continue to rinse for an additional 30 seconds.
 - Drain excess water from the article and dry using a clean soft cloth or filtered pressurized air.
 - Visually examine each instrument for cleanliness. If visible soil is not noted on the instruments, repeat cleaning procedure.
 - Steam autoclave per the following Sterilization Instructions.

Sterility

- Product is supplied **NON-STERILE** and **STERILE**.
- Select plates and screws are available sterile packaged (Gamma Sterilized). **DO NOT USE IF STERILE PACKAGE IS DAMAGED. DO NOT USE AFTER EXPIRATION DATE.**
- OsteoForm Mesh is also available sterile packaged (Gamma Sterilized). **DO NOT USE IF STERILE PACKAGE IS DAMAGED. DO NOT USE AFTER EXPIRATION DATE.**
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- The user facility must clean and disinfect devices prior to sterilization per standard hospital procedures.
- Non-sterile devices are sterilizable by steam sterilization (autoclaving). For sterilization of **OSTEOMED** implant systems, the following parameters should be used.
 - If sterilizing the 210-3018 mini 4 hole straight plate, rigid, within the CFx Rigid Fixation Sustum- Aluminum Tray, the plate will be placed in the locations marked for the 210-0030 plate. A physical depiction of the 210-3018 4 hole plate within the system is shown in the figure below.

Pre-Vacuum Steam Sterilization	ICON Modular Fixation System Tray: 220-0801 MFx Tray	CFx™ Rigid Fixation Systems - Plastic Tray	CFx™ Rigid Fixation Systems - Aluminum Tray	ICON 2-Bay Case Assembly 220-0084
Temperature:	270°F (132°C)	273°F (134°C)	270°F (132°C)	270°F (132°C)
Time:	4 minutes	30 minutes	10 minutes	4 minutes
Dry Time:	80 minutes	55 minutes	55 minutes	60 minutes
Configuration:	Wrapped tray	Wrapped tray	Wrapped tray	Wrapped tray
Wrapping Technique:	Wrap in two layers of 1-ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques	Wrapped tray in two layers of 1-ply polypropylene wrap (Kimguard KC600 – 510(k) K082554)	Wrap in two layers of 1-ply polypropylene wrap (Kimguard KC600-510(k) K082554) with towel placed between the wraps and bottom of tray.	Wrap in two layers of 1-ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope folding techniques

Storage

Sterile packaged implants should be stored at controlled room temperature out of direct sunlight. Product package should be inspected prior to use for signs of damage or tampering.

Caution

- Federal (United States) law restricts this device for sale by or on the order of a medical practitioner licensed to do so.**
- Do not attempt a surgical procedure with faulty, damaged, or suspect OsteoMed instruments or implants. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.**

OsteoMed
3885 Arapahoe Road, USA
Addison, Texas 75001 USA
Customer Service: 1-800-456-7779
Outside USA: 1-972-677-4600

Symbols and Definitions

Single Use Only

Catalogue Number

Use by (Date)

Do not use if sterile package is damaged

Batch Code (Lot Number)

Consult Instructions for Use

Date of Manufacture (MFG Date)

Manufacturer (MFR)

Attention, See instructions for Use Caution, Consult Accompanying Documents

Sterile, Method of Sterilization Using Irradiation

Federal Law (U.S.A.) Restricts this device to sale by or on the order of a physician.

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