Device Description
The OsteoMed EnCompass/Hemi Extractor is an instrument designed to facilitate in removing the OsteoMed EnCompass and/or Hemi Implants from the metatarsal heads and/or phalangeal bases during revision surgeries. The device mates with sizes 8mm through 21mm EnCompass implants and with sizes Extra Small through Large Hemi implants.

Intended Use
The OsteoMed EnCompass/Hemi Extractor is designed to aid in removing existing OsteoMed implants in the metatarsophalangeal (MTP) joint.


Warnings and Precautions

- The OsteoMed EnCompass/Hemi Extractor is not designed to cut into the cortical wall. Use other instrumentation as necessary to cut through the cortical wall and begin separating the implant from the bone.
- The level of effort required to remove the OsteoMed EnCompass/Hemi implant, and the amount of residual bone removed with the implant, will vary based on the length of time it has been implanted, bone growth and other anatomical variables that cannot be predicted.

MRI Status:
The OsteoMed EnCompass/Hemi Extractohas not been evaluated for safety and compatibility in the MR environment. It has not been tested for radio frequency (RF) heating, migration due to magnetically induced displacement, or image artifact in the MR environment. The safety of the implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Maintaining Device Effectiveness

- Carefully inspect the OsteoMed EnCompass/Hemi Extractor prior to use to assure they are in proper operating condition. Instruments which are faulty, damaged, or suspect should not be used. They should be replaced or sent to OsteoMed for disposition.
- OsteoMed products are to be used in a sterile environment.

- Thread the handle (P/N 386-1685) into the clamp (P/N 386-1686) and twist clockwise until the threads are fully engaged.
- The metatarsophalangeal (MTP) joint should be properly exposed utilizing a dorsal incision from the distal third of the shaft of the first metatarsal to a point midway along the shaft of the proximal phalanx. The incision is medial to the extensor hallucis/digitorum longus tendon.
- The incision is deepened to expose the capsule of the metatarsophalangeal joint. The capsular tissues are carefully dissected from the base of the proximal phalanx and metatarsal head.
- Fully expose the metatarsal head or proximal phalanx base, as needed, to ensure full access to the implant.
- If needed, place an osteotome between the implant lip and bony tissue, tapping until the osteotome separates the ingrown tissue from the implant.
- Place the jaws of the EnCompass/Hemi Extractor clamp on the medial and lateral sides of the implant. Turn the handle clockwise to tighten. Ensure that the clamp jaws are fully engaged around the implant lip.
- Thread the EnCompass Shaft (P/N 386-1551) into the EnCompass/Hemi Extractor handle, turning clockwise until it is fully seated.
- Place the EnCompass Slaphammer (P/N 386-1560) at the end of the EnCompass Shaft closest to the EnCompass/Hemi Extractor, and then forcefully pull it away from the Extractor to begin the removal process. Repeat as needed.
- It may be necessary to alternate between the EnCompass/Hemi Extractor and an osteotome during the removal procedure.
- Throughout the removal process, ensure that the EnCompass/Hemi Extractor remains tight around the implant.
- Complete the revision procedure as necessary and close using standard closing techniques.

Cleaning

- Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- OsteoMed recommends the following cleaning and sterilization instructions for the EnCompass/Hemi Extractor:
 - Unscrew the handle from the clamp to clean and sterilize the parts separately.
 - Clean the EnCompass/Hemi Extractor, and associated instruments, thoroughly using a mild detergent, soft brush, and warm water. Ensure that dried blood, bone chips and other deposits are removed from the instruments and sterilization tray (if included).
 - Thoroughly rinse all instruments, and the sterilization tray if included, with water.
 - Arrange all of the instruments in the sterilization tray and ensure that the lid is in place and properly closed. If not using a sterilization tray, arrange all of the instruments in the sterilization pouch and ensure the pouch is properly sealed.

Sterility

- The OsteoMed Encompass/Hemi Extractor is supplied NON-STERILE.
- Use of the sterilizer shall comply with the manufacture's user instructions for sterilizers.
- The user facility must clean and disinfect prior to sterilization per standard hospital procedures.
- Non-sterile devices are sterilizable by steam sterilization (autoclaving). For sterilization of the OsteoMed Encompass/Hemi Extractor , the following parameters should be used:

Pre-Vacuum Steam Sterilization in Encompass Autoclave Tray
Preconditioning Pulses: 3
Temperature: 132° C (270 ° F)
Cycle Time: 15 minutes
Dry Time: 55 minutes
Configuration: Wrap system tray in 1-ply polypropylene wrap using sequential wrapping techniques.
Do not exceed 275 ° F (135 ° C), to avoid compromising functions of polymeric instrumentation

Pre-Vacuum Steam Sterilization in Sterilization Pouch
Preconditioning Pulses: 3
Temperature: 132° C (270 ° F)
Cycle Time: 4 minutes
Dry Time: 10 minutes
Configuration: Seal in two layers of sterilization pouch (double pouched).

Storage
The OsteoMed EnCompass/Hemi Extractor System should be stored at controlled room temperature and used in a sterile environment. Prior to each use, inspect the contents for signs of damage and or defects.

Caution

- Federal (United States) law restricts this device for sale by or on the order of a medical practitioner licensed to do so.**
- Do not attempt a surgical procedure with faulty, damaged, or suspect OsteoMed instruments or implants. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.**

OsteoMed
3885 Arapaho Road
Addison, Texas 75001 USA
Customer Service: 1-800-456-7779
Outside USA: 1-972-677-4600

Symbols and Definitions

	Catalogue Number		Consult Instructions for Use
	Batch Code (Lot Number)		Keep dry
	Date of Manufacture (MFG Date)		Manufacturer (MFR)
	Attention, See Instructions for Use Caution, Consult Accompanying Documents		
	Federal Law (U.S.A) Restricts this device to sale by or on the order of a physician.		Keep Away from Sunlight
	Do not use if sterile package is damaged		

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