Description
The **OSTEOMED ExtremiLOCK Ankle Plating System** is a comprehensive ankle fracture system intended to provide solutions for various bony fractures, including simple to complex fractures of the distal tibia and fibula. The plates are provided in a variety of shapes and sizes, offering surgeons compression and locking hole designs. The **OSTEOMED ExtremiLOCK Ankle Plating System** includes angulated locking screws¹ and standard non-locking screws as well as cannulated screws. Surgical instrumentation is provided to facilitate insertion, modification and/or removal of implants.

Material
The plates are made of titanium (ASTM F67) or titanium-alloy (ASTM F136). The screws are made of titanium-alloy (ASTM F136). K-wires are made of stainless steel (ASTM F138 or F139). The instrumentation is made from various grades of stainless steel, anodized aluminum, and/or medical grade polymers.

Clinical Indications
The **OSTEOMED ExtremiLOCK Ankle Plating System** is intended for fixation of fractures, arthrodesis, osteotomies, and non-unions of the tibia and fibula. The ExtremiLOCK Ankle Plating System implants are intended for single use only. The 1/3 tubular plates, hook plates, screws, and washers are also intended for use in trauma, general surgery, and reconstructive procedures of bones appropriate for the size of the device.

The **OSTEOMED ExtremiLOCK Ankle Plating System** can be used for adult and pediatric patients.

Contraindications
Use of the **OSTEOMED ExtremiLOCK Ankle Plating System** is contraindicated in the following cases:

- Active or suspected infection or in patients who are immunocompromised;
- Patients previously sensitized to titanium or stainless steel or with nickel allergies;
- Patients with certain metabolic diseases;
- Patients who have insufficient bone or poor bone quality;
- Patients exhibiting disorders which would cause the patient to ignore the physician's pre- and/or post-operative instructions and/or the limitations of internal rigid fixation implants;
- Percutaneous K-wire placement is contraindicated in cases of displaced fractures and compressed fractures.

Warnings
The **OSTEOMED ExtremiLOCK Ankle Plating System** is recommended for use in patients with sufficient bone quality to sustain effectiveness and benefits of rigid fixation.

- Use of undersized implants in areas of high functional stress may lead to implant fracture and failure.
- Plates, screws, wires or other appliances of dissimilar metals should not be used together in or near the implant site.
- Bending the plate excessively multiple times may weaken the plate and could result in implant failure.
- Use of screws in highly dense bone may lead to implant fracture or failure upon insertion.
- Drill using the appropriate pilot drill. Note: Use irrigation when drilling.
- When placing additional screws, ensure that subsequent screw placement does not interfere with previously placed screws.
- It is recommended to remove any fractured implants from patients during surgery. If unable to remove, notify patient.
- K-wires and Plate Holding TAKs™ must be removed from the bone fragment prior to compression hole fixation as it will impede the function of the compression hole.
- The appropriate drill guides for locking screws must be used every time a locking screw is inserted to ensure that the locking angle is within ±20° from perpendicular.
- Excessive or off-axis torque may compromise the mechanical lock between the screw and the plate.
- Excessive non-locking screw angulation may cause increased screw head protrusion from the plate. Screw head prominence may cause soft tissue irritation.
- Ensure screws are fully seated/locked, depending on the type screw used to secure the plate.
- Multiple insertions of a locking screw into the same hole may compromise the locking ability of the screw with the plate. If additional insertions are desired, the ability to lock the screw to the plate may decrease. A non-locking screw may be selected for that plate hole, or the surgeon may select a new plate hole location if locking capability is desired.
- Evaluation of the safety and compatibility of the device in the MR environment has not been conducted.
- In considering the evaluation for the safety and compatibility of these devices in the MR environment, the following concerns are raised based on the implant material, per MDD 93/42/EEC and ISO14630: magnetically induced displacement force and torque, radio frequency (RF) heating and image artifacts.

Maintaining Device Effectiveness

- The surgeon must have specific training, experience, and thorough familiarity with the use of internal rigid fixation devices, surgical techniques and post-operative care. It is recommended to follow standard AO operative techniques whenever possible.
- The surgeon must exercise reasonable judgment when deciding which plate and screw to use for specific indications.
- Multi-planar fluoroscopy is recommended throughout screw and plating procedures.
- When loading a screw onto a driver, insert the driver straight into the screw head with force to engage the screw head. Multiple engagements of the driver into the screw head may affect the self-retention feature. To remove the driver from the screw, rock the driver gently from side to side and lift.
- Use the appropriate sized drill guide when drilling or tapping.
- Instruments are marked with the appropriate size for easier identification.
- Plate benders should be used between adjacent holes whenever possible. Bending across empty plate holes may deform the screw holes and prevent a locking screw from engaging the plate.
- After cutting a plate with the plate cutter, utilize the plate file to remove sharp edges on the plate.
- Do not use the compression hole if the fracture is already reduced.
- The **OSTEOMED ExtremiLOCK Ankle Plating System** implants are not intended to endure excessive abnormal functional stresses.
- The **OSTEOMED ExtremiLOCK Ankle Plating System** is intended for temporary fixation only until osteogenesis occurs.
- All **OSTEOMED ExtremiLOCK Ankle Plating System** instrumentation may be required for each surgery. Failure to use dedicated, unique OsteoMed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal. All implants are held in the organizer block.
- Carefully inspect the **OSTEOMED ExtremiLOCK Ankle Plating System** implants prior to use. Inspect the instruments before and after each procedure to ensure they are in proper operating condition. Instruments which are faulty or damaged and/or suspected to be faulty or damaged should not be used. These instruments should be replaced and sent back to OsteoMed for disposition and repair.
- Post-operative instructions should be given to the patient by the surgeon, including the potential for secondary injuries to a surgical site if the patient is non-compliant. Patients should be instructed to closely follow the post-operative instructions.
- OsteoMed recommends the use of OsteoMed products in a sterile environment.
- Depth gauge marking tolerance: 10mm to 70mm is ±0.25mm.

Instructions for Use

- Follow standard rigid fixation technique (i.e. typical AO technique) for placement of the **OSTEOMED ExtremiLOCK Ankle Plating System** implants.
- Refer to the appropriate **OSTEOMED ExtremiLOCK Ankle Plating System** Surgical Technique Guide (STG) for detailed guidelines on proper screw and plate placement; Document No. **030-1839 (Fracture STG)** and **030-1906 (Fusion STG)**.
- Surgical techniques are identified in the title with the following:

FRACTURE	Applicable only to fracture approaches
FUSION	Applicable only to fusion approaches
FRACTURE / FUSION	Applicable to both fracture and fusion approaches

Plating (Contouring)
The ExtremiLOCK ankle plates are pre-contoured to match the anatomy of the distal tibia and fibula. However, if additional contouring is necessary, plate benders or bending irons may be used to achieve the desired contour. The bending irons feature various slot widths to accommodate plate thicknesses. Care should be taken to avoid over-bending or bending plates multiple times to prevent any stress risers.

A plate cutter is also available to cut the distal aspect of the anterior tibia, medial tibia, and antero-lateral plates. The 1/3 tubular plates can also be cut to the desired length. A plate file is available to blunt any sharp edges.

Plating (Provisional Fixation)
FRACTURE – Based on the type of fracture, choose the most appropriate ankle fracture plate. Position the plate over the fracture; care should be taken to ensure adequate points of fixation can be achieved on all sides of the fracture line(s). If desirable, temporary fixation can be achieved using k-wires and Plate Holding TAKs™.

FUSION (Anterior) - Position the anterior fusion plate so that the lateral body tab rests just lateral to the talar neck and the proximal end of the plate rests on the lateral surface of the tibia. If the neck component of the plate does not fit flush with the talus, use a rongeur to remove bone from the dorsal surface of the talar neck. Once the alignment of the plate is satisfactory, sterile marker may be used to identify the medial and lateral edges of the plate on the tibia.

The Box Chisel may be used to flatten the anterior aspect of the tibia.
FUSION (Posterior) - It is important when selecting the position of the posterior fusion plate to be aware of the subtalar joint. The plate should not contact the calcaneus or interfere with subtalar motion. In order to finalize plate position, place a guide wire down the medial talar screw hole to ensure both alignment with the talar neck and the subtalar joint is not violated.

FUSION (Anterior & Posterior) – Temporary fixation can be achieved using k-wires and Plate Holding TAKs. Fixation should begin with the talus.

With the talus secure, the tibial screws are then placed. If compression through the plate is desired, place a non-locking screw through the compression hole of the plate.

Transfixation screws² should be placed next. Proper positioning requires the use of fluoroscopy to ensure talar body capture without violation of the subtalar joint. If desired, the first screw can be placed in a lag fashion to achieve joint compression. The remaining transfixation screws can be either locking or non-locking. Care should be taken to avoid the previously placed talar screws. The remaining proximal tibial screws can be placed per surgeon preference.

Screw Preparation and Insertion Technique
FRACTURE - All circular plate holes can accommodate 2.7mm, 3.5mm and 4.0mm locking and non-locking screws as well as 4.0mm cannulated screws. All locking screws can be locked on-axis with the plate threads or up to 20 degrees angled-locking in any direction (40 degree conical).

FUSION - All circular plate holes, except the tibia positioning hole on the anterior fusion plate, can accommodate 3.5mm and 4.0mm locking and non-locking screws as well as 4.0mm cannulated screws. All locking screws can be locked on-axis with the plate threads or up to 20 degrees angled-locking in any direction (40 degree conical).

FRACTURE / FUSION - Once the appropriate screw diameter has been selected, follow the steps below for screw insertion through a plate:

- Drill** – If using **locking screws**, select the appropriate size angled locking/compression drill guide and insert the guide into the desired plate hole. Ensure the guide is fully engaged in the corresponding plate hole. The cone will ensure the drill remains within the 40° angled locking screw range (±20° from center). If using **non-locking screws**, select the appropriate pilot/over drill guide and insert the pilot drill side through the target plate hole. Continue to drill a pilot hole using the appropriately-sized pilot drill.
- Determine Screw Length** – Select the corresponding depth gauge to determine the correct screw length.
- Screw Insertion** – Select the appropriate screw and insert. Screw length can be verified using the length gauge on the screw block. Final tightening should be performed using a two-finger technique to avoid over-tightening. Locking screws can be locked/unlocked up to three times in a single hole prior to final tightening. Repeat as necessary for additional screw placement.

Cannulated Screw Technique – FRACTURE / FUSION
The ExtremiLOCK 4.0mm cannulated screws can be used with any ankle plate. When using a cannulated screw through a plate, the cannulated screw must be implanted prior to any other screw. Only one cannulated screw can be used per plate. Follow the steps below for cannulated screw insertion through a plate:

- Insert K-Wire** – Insert the k-wire to the appropriate depth under fluoroscopy. Care should be taken to not bend the wire when placing it in bone.
- Determine Screw Length** – Place the cannulated depth gauge over the k-wire and measure the correct screw length. Care should be taken to ensure the distal end of the depth gauge does not damage the plate.
- Drill** – Place the cannulated pilot drill over the k-wire and drill to the corresponding depth.
- Screw Insertion** – Select the appropriate cannulated screw and insert over the k-wire. Screw length can be verified using the length gauge on the screw block. Remove and discard the k-wire.

Plating (Compression Hole Technique) – FRACTURE / FUSION
Compression screw holes are used for providing compression across the fracture site or ankle joint and can only accommodate non-locking screws. Compression is created as the screw travels to the distal side of the compression screw hole. Follow the steps below for screw insertion through a compression hole:

- Drill** – Select the appropriate angled/compression guide and insert into the target compression hole. The arrow will be pointing toward the fracture/fusion site to drill eccentrically. To maximize compression, ensure the drill guide is fully engaged with the target compression hole. Drill using the appropriately-sized pilot drill.
- Determine Screw Length** – Select the corresponding depth gauge to determine the correct screw length.
- Screw Insertion** – Select the appropriate screw and insert. Screw length can be verified using the length gauge on the screw block. Final tightening should be performed using a two-finger technique to avoid over-tightening. Repeat as necessary for additional compression screw placement.

Lag Screw Technique (Solid Core Screws) - FRACTURE
The **OSTEOMED ExtremiLOCK Ankle Plating System** incorporates 2.7mm, 3.5mm and 4.0mm fully-threaded screws and

corresponding instrumentation to address various types of distal tibia and fibula fractures. The following steps describe insertion of a solid core screw using the lag screw technique:

- Drill** – Choose the appropriate overdrill and drill the proximal fragment to create a gliding hole. To maximize compression, the hole should be drilled as perpendicular to the fracture line as possible. Use the appropriate drill guide and corresponding pilot drill to create the pilot hole through the distal fragment. Care should be taken to ensure the fragments are properly reduced prior to drilling through the distal fragment.
- Countersink (Optional)** – If desirable, the countersink may be used to reduce the profile of the screw head and to decrease cortical stress. Place the appropriately-sized countersink in the pre-drilled hole and rotate clock-wise to remove the desired amount of bone to accommodate the screw head. Care should be taken to not remove too much bone.
- Determine Screw Length** – Select the appropriate depth gauge to determine the correct screw length.
- Screw Insertion** – Select the appropriate screw and insert. Screw length can be verified using the length gauge on the screw block. Closely observe compression along the fracture line as the screw is being inserted. Repeat as necessary for additional screw placement.

Lag Screw Technique (Cannulated Screws) - FRACTURE
The **OSTEOMED ExtremiLOCK Ankle Plating System** incorporates 4.0mm fully-threaded cannulated screws and corresponding instrumentation. The following steps describe insertion of a cannulated screw using the lag screw technique:

- Insert K-Wire** – Insert the k-wire to the appropriate depth under fluoroscopy. Care should be taken to not bend the wire when placing it in bone. The cannulated drill guide may be used to facilitate k-wire insertion.
- Determine Screw Length** – Place the Cannulated Depth Gauge / Countersink over the k-wire and measure the correct screw length. If desirable, the distal end of the depth gauge can also be used to remove additional bone to accommodate any screw head prominence.
- Drill** – Place the cannulated overdrill over the k-wire and drill the proximal fragment to create a gliding hole. The Cannulated Drill Guide can be used to facilitate drill insertion. To maximize compression, the hole should be drilled as perpendicular to the fracture line as possible. Care should be taken to ensure the fragments are properly reduced prior to drilling through the distal fragment. Choose the corresponding pilot drill and drill the distal fragment.
- Screw Insertion** – Select the appropriate screw and insert over the k-wire. Screw length can be verified using the length gauge on the screw block. Closely observe compression along the fracture line as the screw is being inserted. Remove and discard the k-wire. Repeat as necessary.

Plating (Transfixation Hole Fixation) - FUSION
The ExtremiLOCK fusion plates feature four transfixation screw holes that provide rigid fixation across the tibia into the talus. Each transfixation screw hole can accommodate 3.5mm and 4.0mm locking or non-locking screws. It is recommended that the first transfixation screw be placed in lag fashion to gain posterior tibiotalar compression.

- Drill** – Select the appropriate overdrill and drill the proximal fragment to create a gliding hole. Care should be taken to ensure the joint is properly reduced prior to drilling through the distal fragment.
- Determine Screw Length** – Select the corresponding depth gauge to determine the correct screw length.
- Screw Insertion** – Select the appropriate screw and insert. Final tightening should be performed using a two-finger technique to avoid over-tightening. Verify fixation with fluoroscopy.

Plating (Tibial Positioning Hole Fixation) - FUSION
The ExtremiLOCK anterior fusion plates feature a tibial positioning hole that allows anterior translation of the tibia relative to the talus. The tibial positioning hole can accommodate 3.5mm and 4.0mm non-locking screws. The following steps describe screw placement through the tibial positioning hole:

- Drill** – Select the uniaxial drill guide and drill through the tibial positioning screw hole.
- Determine Screw Length** – Select the corresponding depth gauge to determine the correct screw length.
- Screw Insertion** – Select the appropriate non-locking screw and insert. The screw can be used for either provisional or permanent fixation. If the screw if used for permanent fixation, care should be taken to ensure the positioning screw does not interfere with any subsequent transfixation screws.

Plating (Syndesmosis Fixation) – FRACTURE
The ExtremiLOCK fibula plates feature a syndesmotic slot that can accommodate up to two non-locking screws. The syndesmotic slot has been designed to place fixation in the correct anatomical location and also accommodate anterior screw angulation. Syndesmotic Reduction Forceps are available to assist with fixation. Follow the steps below for syndesmotic fixation through a plate:

- Drill** – Select the appropriate pilot and drill through the syndesmotic slot.
- Determine Screw Length** – Select the corresponding depth gauge to determine the correct screw length.
- Screw Insertion** – Select the appropriate screw and insert. Final tightening should be performed using a two-finger technique to avoid over-tightening. Verify fixation with fluoroscopy.

Plating (Hook Plate Fixation) - FRACTURE
The ExtremiLOCK Hook plates are universal and intended to be used on either the tibia or fibula. The following steps describe implantation using a hook plate:

- Contour** – If applicable, contour the plate as needed to accommodate the patient's anatomy.
- Pre-Drill (Optional)** – If desirable, the hook drill guide, handle and k-wires or 2.0mm drill may be used to pre-drill for the hooks.
- Placement** – Position the plate so that the hooks capture the bony fragment(s). The hook impactor may be used to secure the hook plate position.
- Screw Insertion** – Insert screws using the previously described Screw Preparation and Insertion Technique.

Cleaning

- Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- OsteoMed recommends the following cleaning and sterilization instructions for re-usable Instrumentation:
- Rinse the articles to be cleaned under running cool tap water to remove visible soil until visibly clean.
- Prepare an enzymatic cleaner, Enzol®, or equivalent, per manufacturer's recommendations. Fully immerse the articles in the solution and soak for a minimum of 2 minutes. Actuate the articles while immersed in the solution to ensure complete penetration of cleaning solution.
- Using a soft bristled brush, clean the entire article paying close attention to hard to reach areas until all evidence of soil is removed. A syringe may be used to clean the lumens and other hard to reach areas. Actuate the article while brushing in order to clean mated surfaces and movable parts.
- Prepare a mild detergent such as Polystycal® 2X Concentrate Neutral, or equivalent, per manufacturer's recommendations. Fully immerse the articles in the prepared solution and sonicate the articles for a minimum of 10 minutes.
- Remove the articles from the detergent and rinse under running reverse osmosis / deionized (RO/DI) water until all evidence of detergent is removed.
- Dry using a clean, soft cloth and filtered pressurized air at ≤40 psi. Visually examine each article for visible soil. Repeat cleaning procedure if necessary.
- Steam Autoclave per the recommended sterilization Instructions.

Sterility

- The **OSTEOMED ExtremiLOCK Ankle Plating System** is supplied **NON-STERILE** and **STERILE**. **STERILE parts are expressly labeled as STERILE**.
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- The user facility must clean and disinfect devices prior to sterilization per standard hospital procedures.
- Non-sterile devices are sterilizable by steam sterilization (autoclaving). For sterilization of the **OSTEOMED ExtremiLOCK Ankle Plating System**, the following parameters should be used:

Pre-Vacuum Steam Sterilization	ExtremiLOCK Ankle System Tray(s) with Screw Modules	Screw Module (Cannulated and Solid Core)	Fusion Module
Part No.	320-2920, 320-2930, 320-2931, 320-2940	320-2924, 320-2925	320-2933
Temperature ¹	270°F (132°C)	270°F (132°C)	270°F (132°C)
Sterilization Time	4 minutes	4 minutes	4 minutes
Tray Preparation ²	Wrapped Tray	Wrapped Block	Wrapped Block
Wrapped Dry Time Minimum	60 minutes (Without Fusion Module)	25 minutes	N/A
Wrapped Dry Time Minimum	80 minutes (With Fusion Module)	N/A	40 minutes
Tray Preparation ³	Rigid Containers	Rigid Container	Rigid Container
Rigid Container Dry Time Minimum	30 minutes	30 minutes	30 minutes
1. Do not exceed 275°F (135°C) to avoid compromising functions of polymeric instrumentation.			
2. Wrapping technique recommendation: Individually wrapped in two layers of 1-ply polypropylene wrap (Kimguard KC600-510(k) K082554) using sequential techniques.			
3. OsteoMed recommends using Aesculap rigid containers & paper filters. The ExtremiLOCK Ankle Plating System (E-APS) must be split and loaded into 2 rigid containers. The 2 upper insert trays from the E-APS shipping container are to be stacked up and placed inside an Aesculap full-sized basket, P/N JF224R (4-1/8" deep). Then the basket and contents are to be loaded into an Aesculap, 8-5/8" deep, full-size, rigid container, P/N JN446/JK486 (bottom/lid). The lower 3 modules from the E-APS shipping container, 320-2924, 632-2925, & 320-2933, are to be placed in a shorter Aesculap full-size basket, P/N JF221R (1-1/2" deep). Then this basket and contents are to be placed inside a 2nd Aesculap, 3-5/8" deep, full-size, rigid container, P/N JN441/JK486 (bottom/lid). The 2 rigid containers can be sterilized in the same unit at the same time. Use 4 Aesculap paper filters, P/N US751, 2 in the lid and 2 in the bottom for each rigid container.			

Note: Biological Indicator of *G. stearotherophilus* was used in sterilization validations.

Storage
OSTEOMED ExtremiLOCK Ankle Plating System should be stored at controlled room temperature, away from direct sunlight, and used in a sterile environment. Prior to each use, inspect the contents of **OSTEOMED ExtremiLOCK Ankle Plating System** for signs of damage and/or defects.

Caution

- Federal (United States) law restricts this device for sale by or on the order of a medical practitioner licensed to do so.**
- Do not attempt a surgical procedure with faulty, damaged, or suspect OsteoMed instruments or implants. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.**

OsteoMed
3885 Arapaho Road
Addison, Texas 75001 USA
Customer Service: 800/456-7779
Outside USA: 972/677-4600

Symbols and Definitions

Single Use Only

REF

Catalogue Number

Use By (Date)

Do not use if sterile package is damaged

LOT

Batch Code (Lot Number)

Consult Instructions for Use

Date of Manufacture (MFG Date)

Manufacturer (MFR)

Attention, See Instructions for Use Caution, Consult Accompanying Documents

Keep Away from Sunlight

STERILE R

Sterile, Method of Sterilization Using Irradiation

Keep Dry

Rx ONLY

Federal Law (U.S.A) Restricts this device to sale by or on the order of a physician.

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