OSTEOMED
OSTEO*MED* ExtremiLOCK Wrist Plating System
Product Information and Instructions for Use

Description

The **OSTEO*MED* ExtremiLOCK Wrist Plating System** is a rigid fixation system consisting of plates and screws in various configurations. Plates are provided in a variety of shapes and sizes and offer surgeons positioning and locking hole designs. The Wrist Plating System includes variable angled locking, non-locking, fully threaded, partially threaded, and cannulated screws, washers, as well as K-wire implants. Surgical instrumentation is provided to facilitate modification and insertion of implants. Variable locking angle range for 2.4mm screws is ±15 degrees and for 2.7mm screws is ±10 degrees.

Material

The plates are made from titanium (ASTM F-67); Screws are made from titanium alloy (ASTM F-136); K-wires are made of stainless steel (ASTM F-138). The instrumentation is made from various grades of stainless steel, titanium alloy, anodized aluminum, and/or medical grade polymers.

Clinical Indications

OSTEO*MED* ExtremiLOCK Wrist Plating System is intended for fracture fixation, fusion, and osteotomies of wrist and other bones appropriate for the size of the device. It is intended for use in trauma, general surgery, and reconstructive procedures.

OSTEO*MED* ExtremiLOCK Wrist Plating System implants are intended for single use only.

Contraindications

Use of the **OSTEO*MED* ExtremiLOCK Wrist Plating System** is contraindicated in the following cases:

- Active or suspected infection or in patients who are immunocompromised;
- Patients previously sensitized to titanium or stainless steel or with nickel allergies;
- Patients with certain metabolic diseases;
- Patients where there is insufficient bone;
- Patients exhibiting disorders which would cause the patient to ignore the physician's pre- and/or post-operative instructions and/or the limitations of internal rigid fixation implants;

Warnings

- If material sensitivity is suspected, tests should be performed prior to implantation.
- Use of undersized implants in areas of high functional stress may lead to implant fracture and failure.
- The implants are not intended to endure excessive abnormal functional stresses.
- The **OSTEO*MED* ExtremiLOCK Wrist Plating System** is intended for temporary fixation only until osteogenesis occurs.
- Plates, screws, K-wires or other devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects. Do not mix different metal implants in the same construct.
- Multiple/excessive bending may weaken the plate and could result in implant fracture and failure.
- Use of screws in high density bone may lead to implant fracture or failure upon insertion.
- The **OSTEO*MED* ExtremiLOCK Wrist Plating System** is recommended for use in patients with sufficient bone quality to sustain effectiveness and benefits of rigid fixation.
- It is recommended to remove any fractured implants from patient during surgery. If unable to remove, notify patient.
- Use of excessive torque during insertion of screws may lead to implant failure.
- Excessive screw angulations may cause increased screw protrusion from plate.
- For proper device placement, use of a drill guide is recommended.
- If screw is not fully seated into plate hole, it may lead to non-union, misaligned fusion, or tissue irritation.
- The **OSTEO*MED* Wrist Plating System** is recommended for use with fluoroscopy during screw insertion and drilling to prevent interference with previously placed devices.
- Locking screws are recommended for unicortical fixation.
- Inform patient if he/she does not follow physician post-operative instructions and limitations of the internal rigid fixation implants, this may lead to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Do not use locking screws with oblong or positioning holes. Locking screws are intended for use through locking plate holes only.
- Multiple engagements of screw to screwdriver tip may result in loss of retention.
- Locking screws and plate holes can be used up to three times.
- Use irrigation when drilling.

MRI Safety Information

- The **OSTEO*MED* ExtremiLOCK Wrist Plating System** has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of **OSTEO*MED* ExtremiLOCK Wrist Plating System** in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Maintaining Device Effectiveness

- The surgeon must have specific training, experience, and thorough familiarity with the use of internal rigid fixation devices, surgical techniques, and post-operative care. Patients must closely follow the post-operative instructions from their surgeon.
- The surgeon must exercise reasonable judgment when deciding which implant or instrument type to use for specific indication.
- Failure to use dedicated, unique OsteoMed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure and subsequent patient injury. Failed devices may require re-operation and removal.
- Carefully inspect the **OSTEO*MED* ExtremiLOCK Wrist Plating System** implants and instruments prior to use.
- Verify screw length on the screw block prior to use.
- Inspect the instruments before and after each procedure to assure they are in proper operating condition. Instruments which are faulty or damaged, and/or suspected to be faulty or damaged should not be used. These instruments should be replaced and sent back to OsteoMed.
- Depth gauge marking tolerance: 8mm to 40mm is ±0.25mm.

Instructions for Use


- Use the **OSTEO*MED* ExtremiLOCK Wrist Plating System** in a sterile environment.
- Follow standard rigid fixation technique (i.e. typical AO technique) for placement of the **OSTEO*MED* ExtremiLOCK Wrist Plating System** implants.
- Refer to the **OSTEO*MED* ExtremiLOCK Wrist Plating System** Surgical Technique Guide for detailed instructions on the instrumentation use and/or modification and implantation of the system implants.

Attention, **OSTEO*MED* ExtremiLOCK Wrist Plating System** Surgical Technique Guide P/N 030-1952

Cleaning

- Reusable instrumentation must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- Depth Gauge Disassembly instructions:**

- Remove nose/cap of depth gauge from main body by gently rocking it off the main body. Exercise caution not to use too much force to remove the nose/cap as the probe may become damaged if done so.
- Remove slider/handle from main body

- Depth Gauge Assembly instructions**
 - Insert slider/handle through the main body
 - Place the nose/cap of the depth gage over the probe and snap it into the main body

Cleaning (continuation)OsteoMed recommends the following cleaning instructions for the OsteoMed reusable instrumentation:

- Rinse the articles to be cleaned under running cool tap water to remove gross soil.
- Use a sterile syringe to flush water through any cracks, crevices, lumens, and hard to reach areas. During rinsing, actuate the articles to ensure thorough rinsing.
- Prepare an enzymatic cleaner such as Enzol® per manufacturer's recommendations at 1oz/gal using lukewarm tap water. Using a clean soft cloth that has been soaked in the solution, wipe the entire article. Use a sterile syringe to flush at least 50mL of enzymatic cleaner through any cracks, crevices, lumens, and hard to reach areas. Actuate the articles while immersed to ensure complete penetration of cleaning solution. Allow the articles to soak in the solution for a minimum of 15 minutes.
- After soaking, thoroughly brush the articles beneath the surface of the prepared cleaner using a soft bristled brush, paying close attention to all hard to reach areas until all evidence of soil is removed. Actuate the articles while brushing in order to clean matted surfaces and movable parts. Using a lumen brush or similar type brush, brush each lumen a minimum of 5 times, ensuring that the lumen brush is passed completely through the entire lumen during brushing. Use a syringe to flush any lumen or matted surface with at least 50mL of solution.
- Rinse the articles in RO/DI (reverse-osmosis/deionized) water until all visible evidence of detergent is removed. Flush any lumens or matted surfaces with the RO/DI water. Once all evidence of detergent is removed continue to rinse for an additional 30 seconds.
- Drain excess water from the article and dry using a clean soft cloth or filtered pressurized air.
- Prepare a non-enzymatic detergent such as MetriWash™ per manufacturer's recommendations at ¼ oz/gal of lukewarm tap water.
- Fully immerse the articles and flush all lumens with at least 50mL of the detergent.
- Allow the articles to soak in the solution for a minimum of 15 minutes.
- Thoroughly brush the exterior of the articles beneath the surface of the prepared cleaner using a soft bristled brush. Using a lumen brush or similar type brush, brush each lumen a minimum of 5 times, ensuring that the lumen brush is passed completely through the entire lumen during brushing. Flush all lumens with 50mL of solution after brushing.
- Use a mild detergent such as Metriwash™ and prepare per manufacturer's recommendations of ¼ oz/gal using lukewarm tap water in a sonication unit. Fully immerse the articles in the prepared detergent. Use a sterile syringe to flush at least 50mL of the prepared detergent through any cracks, crevices, lumens, and hard to reach areas. Actuate the articles while immersed to ensure complete penetration of the cleaning solution. Sonicate the articles for a minimum of 10 minutes.
- Remove the articles from the sonicator and rinse under RO/DI water for a minimum of 2 minutes. While rinsing, actuate the articles to ensure a thorough rinse until all visible evidence of detergent is removed. Use a sterile syringe to rinse any lumens or other hard to reach areas with a minimum of 60mL of RO/DI water.
- Drain excess water and dry the articles using a clean soft cloth and filtered pressurized air at ≤ 40 psi. Visually inspect each article without magnification for visible soil, deterioration, or loss of function. If soiled, repeat cleaning process. If device is visibly deteriorated or is unable to function, device shall be returned to OsteoMed.

Sterility

- The **OSTEO*MED* ExtremiLOCK Wrist Plating System** implants are supplied **STERILE** and **NON-STERILE**
- Implants which are provided **sterile** packaged are Gamma Sterilized. **DO NOT USE IF IMPLANT STERILE PACKAGE IS DAMAGED. DO NOT USE IMPLANTS AFTER EXPIRATION DATE.**
- Implants and instruments supplied **non-sterile** must be sterilized prior to use.
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- The user facility must clean and disinfect reusable instrumentation prior to sterilization per standard hospital procedures and OsteoMed cleaning instructions.
- Do not stack fully loaded trays during sterilization.
- Non-sterile** devices and instrumentation are sterilizable by steam sterilization (autoclaving). For sterilization of the **OSTEO*MED* Wrist Plating System**, the following parameters should be used:

Pre-Vacuum Steam Sterilization	ExtremiLOCK Wrist System Tray			
Configuration	Wrapped Tray ¹	Wrapped Tray ¹	Rigid Container ²	Rigid Container ²
Temperature	270°F (132°C)	273°F (134°C)	270°F (132°C)	273°F (134°C)
Sterilization Time	4 minutes	3 minutes	4 minutes	3 minutes
Minimum Dry Time	30 minutes	40 minutes	30 minutes	30 minutes
1. Wrapping Technique	Individually wrap in two layers of 1-ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential wrapping techniques. NOTE: If specified wrap is not available, only FDA cleared wraps should be used to wrap the subject device.			
2. Rigid Container Technique	OsteoMed recommends using Aesculap rigid containers and paper filters. A single level ExtremiLOCK Wrist Plating System tray and contents may be sterilized in an Aesculap container base P/N JN446 with lid P/N JK489. Use 4 Aesculap paper filters P/N US751, 2 in the lid and 2 in the bottom of the rigid container.			
Do not exceed 278°F (137°C) to avoid compromising functions of polymeric instrumentation.				

Note: Biological Indicator of *G. stearothermophilus* was used in sterilization validation.

The OsteoMed ExtremiLOCK Wrist Plating System implants labeled as sterile are non-pyrogenic. This labeling claim has been validated to the endotoxin limit of no more than 20.0 EU per device

Storage

- OSTEO*MED* ExtremiLOCK Wrist Plating System** should be stored at controlled room temperature, away from moisture and direct sunlight.
- Sterile packaged devices must be stored according to the labeled storage instructions. Do not use sterile packaged devices if the package has been damaged or otherwise tampered.
- Prior to each use, inspect the contents of **OSTEO*MED* ExtremiLOCK Wrist Plating System** for signs of damage and/or defects.

Caution

- Federal (United States) law restricts this device for sale by or on the order of a medical practitioner licensed to do so.**
- Do not attempt a surgical procedure with faulty, damaged, or suspect OsteoMed instruments or implants. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.**

OsteoMed
3885 Arapaho Road
Addison, Texas 75001 USA
Customer Service: 1-800-456-7779
Outside USA: 1-972-677-4600

Symbols and Definitions

Single Use Only

Catalogue Number

Use By
(Date)

Do not use if sterile
package is damaged

Batch Code
(Lot Number)

Consult Instructions for
Use

Date of Manufacture
(MFG Date)

Manufacturer
(MFR)

Attention, See
Instructions for Use
Caution, Consult
Accompanying
Documents

Keep Away from
Sunlight

Sterile, Method of
Sterilization Using
Irradiation

Keep Dry

Federal Law (U.S.A)
Restricts this device to
sale by or on the order
of a physician.

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