OsteoMed® ExtremiFix Cannulated Screw System

EN	Product Information and Instructions for Use	2
FR	Informations sur le produit et instructions d'utilisation.....	6
DE	Produktinformationen und Gebrauchsanweisung ...	10
IT	Informazioni sul prodotto e istruzioni per l'uso	14
PT	Informações sobre o produto e instruções de utilização	18
ES	Información sobre el producto e instrucciones de uso.....	22

Product Information and Instructions for Use

OsteoMed® ExtremiFix Cannulated Screw System

EN

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals in the European Union only.

Description

OSTEOMED EXTREMI-FIX CANNULATED SCREW SYSTEM is a rigid fixation system consisting of screws in diameters of 2.0mm (10-42mm length), 2.4mm (10-50mm length), 3.0mm (10-40mm length), and 4.0mm (12-52mm length). The system instruments include guide wires, drills, drivers, handle, depth gauge, bone clamp, screw extractor and preparation instruments to facilitate the placement of screws.

Material

OSTEOMED EXTREMI-FIX CANNULATED SCREW SYSTEM implants are made from titanium alloy (ASTM F-136). The instrumentation is made from various grades of stainless steel, anodized aluminum, and/or medical grade plastic.

Intended Purpose & Indications for Use

OSTEOMED EXTREMI-FIX CANNULATED SCREW SYSTEM is indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, fracture repair, and fracture fixation of bones appropriate for the size of the device. The screws are intended for single use only. The system drills and guide wires are single use instruments.

The intended clinical benefit of the **OSTEOMED EXTREMI-FIX CANNULATED SCREW SYSTEM** is to provide fixation for patients with fractures of bones appropriate for the size of the device as a result of trauma, Osteotomy or reconstruction of the bones. The indirect clinical benefit of the encompassed instrumentation is to facilitate implantation protocols for fracture fixation.

Contraindications

- Active or latent infection or marked inflammation of the treatment area
- Sepsis
- Insufficient quantity or quality of bone, osteoporosis, or in patients with certain metabolic diseases
- Patients with confirmed material sensitivity
- Patients who are unwilling or incapable of following post-operative care instructions

Warnings

- The treatment or implant may fail, including sudden failure, because of:
 - Loose fixation and/or loosening
 - Stress, including stress from inappropriate bending of the implant during surgery
 - Stress concentrations
 - Stress of weight bearing, load bearing, or excessive activity
Failure is more likely if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing. Failure is more likely if the patient does not follow post-operative care instructions.
- Nerve or soft tissue damage may result from surgical trauma or the presence of an implant.
- Instrument breakage or damage, as well as tissue damage, may occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use, or unintended use.
- Devices of dissimilar material should not be used together in or near the implant site. Dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects.
- Fractured implants should be removed from patients during surgery. If unable to remove, notify the patient.
- Implants may cause distortion and/or block the view of anatomic structures on radiographic images.

Cautions

- Implants and instruments are intended only for professional use by a licensed physician.
- Use devices in accordance with the Instructions for Use.
- Do not use sterile products past the use-by date. Refer to the device label.
- Do not use or re-sterilize an implant provided in sterile packaging if the package has been damaged. The sterility may be compromised, and the cleanliness of the implant may be uncertain. Report damaged packaging to your distributor or Acumed.
- Use of this system has not been evaluated in children or individuals that are not skeletally mature.
- Inspect all components preoperatively to ensure utility. Do not attempt a surgical procedure with faulty, damaged, or suspect instruments or implants. Alternate fixation methods should be available intraoperatively.
- Reprocessing and/or reuse of single use devices may result in infection/cross contamination, and/or sudden failure due to previous stresses.
- Use of implant components from different manufacturers has not been evaluated.
- Use of chemical disinfection may leave residues that adversely affect steam sterilization.
- Steam penetration and device sterilization may be negatively affected by labeling that block tray steam holes.
- Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling of biohazardous material and disposal of sharp materials.

Potential Adverse Effects

- Anesthesia-related problems, problems with positioning of the patient (e.g., nausea, vomiting, neurological impairments, etc.) thrombosis, embolism, infection, or injury to other critical structures such as blood vessels, excessive bleeding, etc.
- Pain, discomfort, or abnormal sensations, nerve or soft tissue damage, necrosis of bone or tissue, bone resorption, or inadequate healing from the presence of an implant or due to surgical trauma.
- Swelling, abnormal scarring, impairment in musculoskeletal function, hardware prominence, migration, loosening, bending, or breakage of the implant, malunion, nonunion or delayed union due to prolonged loading or excessive forces which may lead to implant failure and reoperation.
- Metal sensitivity, histological, allergic, or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at <http://www.acumed.net/ifu>.
- Injury to user.

Note: Any serious incident that has occurred in relation to the device should be reported to Acumed (through mailto customercomplaints@acumed.net or +1.888.627.9957) and the Competent Authority of the Member State in which the user and/or patient is established.

Target Population

The system is intended to be used on skeletally mature, adult patients with sufficient quantity and quality of bone who require surgical intervention to achieve fracture fixation, fusion or reconstruction of bones.

MRI Status

Many Acumed implants have been evaluated for safety in the MR environment. Consult our publication "Acumed Implants in the MR Environment" at <http://www.acumed.net/ifu> for more information.

Device Lifetime

Once installed, implants are expected to provide fixation, fusion and physiological support and have an effective life during bone healing. The implants are biocompatible and may remain implanted at the discretion of the surgeon. Multiple-use instruments have a lifetime that is affected by usage, handling, and processing. Assess multiple-use instruments for fitness during pre-sterilization inspection.

Summary of Safety and Clinical Performance (SSCP)

The SSCP for the implants may be obtained from the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Use the following search term(s): 08456940BUDIO23NU

Maintaining Device Effectiveness

- The surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation products and techniques.
- The surgeon must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
- Devices should be used in an operating room environment.
- All Acumed plates, screws, and instrumentation may be required for each surgery. Failure to use dedicated, unique Acumed instruments for every step of the implantation technique may compromise the integrity of the implanted device, leading to premature device failure that could cause subsequent patient injury. Failed devices may require re-operation and removal.
- Inspect the Acumed implants prior to use. Do not attempt a surgical procedure with faulty, damaged, or suspected instruments or implants. Alternate fixation methods should be available intraoperatively.
- Instrument lifetime is affected by usage, handling, and processing. Assess multiple-use instruments for fitness before and after each procedure to confirm that they are in proper operating condition. Indications that a reusable instrument should not continue to be used include but are not limited to deformation, dulled sharps, cracks, breakage, missing components, corrosion, excessive wear, and restriction or seizure of moving parts. Instruments which are faulty, damaged, or suspected should not be used. They should be replaced or sent to Acumed for disposition and repair.
- When placing more than one screw, ensure that subsequent screw placement does not interfere with other screws. Insert the second screw on the opposite side of the fracture or osteotomy site, and then all remaining screws, following the outlined procedures.

Instructions for Use – OsteoMed ExtremiFix Cannulated Screw System

- Acumed offers one or more Surgical Technique Guides to promote the safe and effective use of this system. Consult our Surgical Technique Guides at www.acumed.net. The Surgical Technique Guides also include information about the installation and calibration of the device.
- **Important:** Surgical Technique Guides may contain important safety information. The instruments and implants in this system are intended to be used by suitably trained and qualified surgeons in a hospital operating room setting. Before treatment, the surgeon is advised to read and fully understand all instructions and communicate to the patient any relevant medical information provided therein, including the use, limitations, risks (safety communications), and possible adverse effects of the proposed treatment. Consult the most recent versions of the Instructions for Use and Surgical Technique Guides as they are subject to change. Contact Acumed or an authorized agent for additional information.

Cleaning

- Products must be carefully cleaned prior to sterilization. Trained personnel must perform cleaning and mechanical inspection prior to sterilization.
- Compliance is required with the equipment manufacturer's user instructions (manual and/or machine cleaning, ultrasound treatment, etc.) and recommendations for chemical detergents.
- OsteoMed recommends the following cleaning and sterilization instructions for Instrumentation:
 1. Rinse the articles to be cleaned under running cool tap water to remove gross soil.
 2. Use a syringe to flush water through any cracks, crevices, lumens, and hard to reach areas. During rinsing, actuate the articles to ensure thorough rinsing.
 3. Prepare a neutral enzymatic cleaning solution per the manufacturer's instructions. Pre-soak soiled instruments for a minimum of 1 minute in the enzymatic solution.
 4. After soaking, thoroughly brush the articles beneath the surface of the cleaning solution using a soft bristled brush, paying close attention to all hard to reach areas until all evidence of soil is removed. Using a lumen brush or similar type brush, brush each lumen a minimum of 5 times, ensuring that the lumen brush is passed completely through or to the bottom of the entire lumen during brushing. Use a twisting action where possible. Actuate the articles while brushing in order to clean mated surfaces and movable parts.
 5. Rinse instruments under cool running water for at least 30 seconds to flush debris out of any tight spaces. Use a syringe to flush each lumen or mated surfaces.
 6. Drain excess water from the article and dry using a clean, soft cloth or filtered pressurized air at < 40psi. Visually inspect each article without magnification for visible soil, deterioration, or loss of function. If soiled, repeat cleaning process. If the device is visibly deteriorated or is unable to function, device shall be returned to OsteoMed.
 7. Return the instruments to their sterilization case.

Thermal Disinfection

For thermal disinfection, remove trays from outer case and place separately in the washer/disinfector's rack during automated cleaning and thermal disinfection. Utilize a thermal disinfection program consistent with A0 value > 3000 or perform thermal disinfection at 90°C (maximum temperature 95°C) for at least 5 minutes.

Pre-Sterilization Inspection

- Visually inspect all devices under normal lighting to ensure that cleaning was effective. Pay close attention to all challenging areas. Reprocess instruments that are not clean and replace instruments that cannot be cleaned.
- Inspect the implants and instruments for surface damage, such as nicks, scratches, and cracks. Replace affected devices.
- Assess the instruments for proper use. Operate all parts and connecting mechanisms. Give careful attention to drivers, drill bits and reamers, and instruments used for cutting or implant insertion. Critically assess them for wear, sharpness, straightness, and corrosion. Replace any instrument that does not perform as intended.
- Inspect cutting edges under magnification. Replace worn instruments e.g., dull, chipped, cracked, rolled, or otherwise deformed. Run a cotton cloth over the edge to help detect chipping and cracking.
- Verify the legibility of all markings and reference scales. Replace any device that is unreadable.
- Repair, replace, and/or repeat the cleaning of instruments as needed to ensure proper operation before proceeding with sterilization.
- Lubricate (instrument milk) surgical instruments to increase useful life. Do not use silicone-based lubricants, oil, or grease, as these will interfere with steam sterilization. Use a water-based lubricant intended for use on surgical instruments and with steam sterilization. Use the lubricant as directed by the manufacturer. Use critical water if dilution is required.
- Fully replenish the system trays and caddies.

Sterility

- Product is supplied **NON-STERILE unless expressly labeled as STERILE.**
- **DO NOT USE IF THE STERILE PACKAGE IS DAMAGED. DO NOT USE AFTER EXPIRATION DATE.**
- Use of the sterilizer shall comply with the manufacturer's user instructions for sterilizers.
- Implants and instruments supplied non-sterile must be sterilized prior to use.
- The user facility must clean and dry devices prior to sterilization per standard hospital procedures.
- Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017, the Association for the Advancement of Medical Instrumentation's (AAMI's) "Comprehensive guide to steam sterilization and sterility assurance in health care facilities."
- Non-sterile devices are sterilizable by steam sterilization (autoclaving). For sterilization of Acumed implant systems, the following parameters in the table below should be used. The table shows the minimum parameters validated to achieve a required Sterility Assurance Level (SAL) of 10⁻⁶ for the system.

The ExtremiFix System is validated for steam sterilization using one of the below sets of parameters.		
Pre-Vacuum Steam Sterilization	OsteoMed ExtremiFix Cannulated Screw System (2.0/2.4mm and 3.0/4.0mm)	
Part Numbers:	316-1250, 316-1350	
Configuration 1:	Wrapped tray	
Temperature:	132°C	134°C
Full Cycle Time:	4 minutes	3 minutes
Wrapped Dry Time Minimum:	30 minutes	55 minutes
Wrapped Configuration Details:	Wrap in two layers of 1- ply polypropylene wrap (Kimguard KC600 – 510(k) K082554) using sequential envelope techniques.	
Configuration 2:	Rigid Sterilization Container (RSC)*	
Temperature:	132°C	134°C
Full Cycle Time:	4 minutes	3 minutes
RSC Dry Time Minimum:	30 minutes	30 minutes
<i>Note: Biological indicator, G. stearothermophilus, was used in sterilization validation.</i>		
*Acumed recommends using Aesculap rigid sterilization containers and paper filters. Place the ExtremiFix System into a half-size Aesculap RSC, such as JN341 base and JK38X lid. Use 2 Aesculap paper filters, US751, 1 in the lid and 1 in the base.		

Post-Sterilization Inspection

- Do not store or use sterile devices if they are not dry.
- Inspect the sterile barrier for signs of damage. Do not use the product if the sterile barrier has been compromised.

Storage


- Store sterilized devices under controlled conditions in a manner that minimizes the potential for contamination per ANSI/AAMI ST79:2017. Refer to sterilization wrap or rigid container manufacturer's IFU for limits on sterile product storage time and storage requirements for temperature and humidity.
- Devices should be stored in an area that provides protection from dust, pests, and temperature/ humidity extremes.

Safe Disposal

Devices may be contaminated with biohazardous materials and/or sharp materials. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling biohazardous material and disposal of sharp materials.

Symbols and Definitions

Note: This is a general list of symbols. Not all symbols apply to all systems.

Symbol	Description	ISO 15223-1	Symbol	Description	ISO 15223-1
	Manufacturer	5.1.1		Authorized Representative in the European Community	5.1.2
	Date of Manufacture	5.1.3		Use-by Date	5.1.4
	Batch Code (Lot Number)	5.1.5		Catalogue Number	5.1.6
	Serial number	5.1.7		Sterilized using irradiation	5.2.4
	Do Not Resterilize	5.2.6		Non-Sterile	5.2.7
	Do not use if sterile package is damaged	5.2.8		Single Sterile Barrier System	5.2.11
	Double Sterile Barrier System	5.2.12		Single Sterile Barrier System with Protective packaging inside	5.2.13
	Single Sterile Barrier System with Protective packaging inside	5.2.14		Single Use Only	5.4.2
 www.acumed.net/ifu	Consult Instructions for Use	5.4.3	Rx only	Caution: U.S. federal law restricts this device to sale by or on the order of a physician.	U.S. 21 CFR 801.109
	Medical Device	5.7.7		MR Conditional (ASTM F2503)	
	CE marking of conformity, Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment.				

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CUSTOMER SERVICE (OUTSIDE USA):

+1.888.627.9957

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

PKGI-014293-Rev A | 11-2024 | © 2024 Acumed® LLC

OsteoMed® is a registered trademark of Acumed, LLC.

Subsidiaries:

Acumed Ltd

Huebner House

The Fairground

Andover

Hampshire UK SP11 0QN

Tel: +44 1264 774450

Acumed Iberica

C/ Procion 1

Edificio Oficor

28023 Madrid-Spain

Tel: +34 913516357

Acumed GmbH

Fuhlsbuttler Stasse 300

22307 Hamburg

Deutschland

Tel: +49-40 947 82 093

EMERGO EUROPE

WESTERVOORTSEDIJK 60

6827 AT ARNHEM

THE NETHERLANDS

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ACUMED, LLC