



RIDGE-FORM™
ALVEOLAR RIDGE RECONSTRUCTION MESH



RIDGE-FORM™

Alveolar Ridge Reconstruction Mesh


- Smooth pre-curved design requiring limited contouring
- Rigid titanium porous mesh providing optimal stability and enhanced vascularity to the graft site
- Anterior mesh for reconstruction for the pre-maxilla/mandible
- Left and right Hemi meshes for full or partial ridge augmentation
- Auto-Drive™ self-drilling fixation screws

Ridge-Form™ Meshes

11mm Ridge Form™ Mesh



224-0030 Anterior Ridge Form™ Mesh




224-0031 Hemi, Left Mandible/Right Maxilla
224-0032 Hemi, Right Mandible/Left Maxilla

16mm Ridge Form™ Mesh



224-0033 Anterior Ridge Form™ Mesh



224-0034 Hemi, Left Mandible/Right Maxilla
224-0035 Hemi, Right Mandible/Left Maxilla

Instruments



220-0019 Taperlock™ Screwdriver Body



220-0061 1.6mm Screwdriver Shaft, Manual



220-0067 1.6mm Screwdriver Shaft, Dental-Latch



3885 Arapaho Road
Addison, Texas 75001
Customer Service: 800.456.7779
Main: 972.677.4600
Fax: 800.390.2620
www.osteomedcorp.com