








The following symbols may be used in Acumed labeling:

Symbol	Description	Reference*
	Manufacturer Indicates the medical device manufacturer, as defined by EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	5.1.1
	Authorized representative in the European Community Indicates the authorized representative in the European Community	5.1.2
	Date of manufacture Indicates the date when the medical device was manufactured.	5.1.3
	Use-by date Indicates the date after which the medical device is not to be used.	5.1.4
	Batch code Indicates the manufacturer's batch code so that the batch or lot can be identified.	5.1.5
	Catalogue number Indicates the manufacturer's catalogue number so that the medical device can be identified.	5.1.6
	Sterilized using ethylene oxide Indicates a medical device that has been sterilized using ethylene oxide.	5.2.3
	Sterilized using irradiation Indicates a medical device that has been sterilized using irradiation.	5.2.4

Symbol	Description	Reference*
	Do not re-sterilize Indicates a medical device that is not to be re-sterilized.	5.2.6
	Non-sterile Indicates a medical device that has not been subjected to a sterilization process.	5.2.7
	Do not use if package is damaged Indicates a medical device that should not be used if the package has been damaged or opened.	5.2.8
	Upper limit of temperature Indicates the upper limit of temperature to which the medical device can be safely exposed.	5.3.6
	Do not re-use Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	5.4.2
	Consult instructions for use Indicates the need for the user to consult the instructions for use.	5.4.3**
	Caution Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	5.4.4

* ISO 15223-1 2nd Ed. 2012-07-01, Medical Devices - Symbols to be used with medical devices labels, labeling, and information to be supplied - Part 1: General requirements

** An eIFU indicator below symbol 5.4.3 indicates the website URL where instructions for use are made available in electronic format.