


Large Cannulated Screw System

6.5/7.3mm Screws



- Self-Drilling and Self-Tapping Lag Screws
- Low-Profile Screw Heads
- Percutaneous Implant Delivery
- Surgeon-Driven Instrumentation



Rethinking Possibilities, Reshaping Lives


Large Cannulated Screw System

6.5/7.3mm Screws

The **Large Cannulated Screw System** is the premier cannulated system for extremity fixation. The system offers 6.5mm/7.3mm screws and washers used for bone fixation following trauma or osteotomy. A wide range of screw lengths are available with both 20mm and 32mm thread lengths to accommodate multiple procedures and patient anatomy. The comprehensive instrument offering can address either traditional open or percutaneous surgical approaches.



Features and Benefits



Low Profile Screw Head
Reduced potential for soft tissue irritation

Washers
Provide a larger head surface area for use in osteoporotic bone



Titanium Alloy
Ti 6Al-4v ELI

Wide Range of Screw Lengths
6.5mm x 40mm-120mm
7.3mm x 40mm-120mm

Reverse Cutting Flute
Facilitates screw removal

Thread Choice
Implants available with partial short or long threaded options

Self-Drilling / Tapping Design
Efficient cutting tip design facilitates time-saving surgical technique

Guide Wires
Partially threaded or smooth guide wire options available

OsteoMed Products



ExtremiLOCK Ankle Plating System



ExtremiLOCK Ankle Fusion System




ExtremiLOCK Foot Plating System



ExtremiFIX Cannulated Screws



ExtremiFuse



EnCompass



EnCompass Lessers



Hemi




ReFlexion




InterPhlex




Talar-Fit



Inion



OsteoVationEX



OsteoVation QWIK



OSTEOMED
3885 Arapaho Rd.
Addison, TX 75001
Customer Service: 800.456.7779
Outside the U.S.: 001.972.677.4600
Fax: 800.390.2620
Fax Outside the U.S.: 001.972.677.4709
E-mail: customer.service@osteomed.com
www.osteomed.com

P/N 030-1437 Rev.D