

# M3<sup>TM</sup>-X

## EXTREMITY FIXATION SYSTEM



- Solid core screw fixation system designed for a variety of indications in the upper and lower extremity
- Titanium, self-tapping and low profile self-retaining screw head
- Lag and fully threaded screws offered in multiple thread lengths and diameters
- Color coded instrumentation to assure accuracy




*Rethinking Possibilities, Reshaping Lives*

# M3-X™

## EXTREMITY FIXATION SYSTEM

The **M3-X™ EXTREMITY FIXATION SYSTEM** is a solid core screw system for small bone osteotomies, fractures, arthrodesis and reconstruction. The system features lag screws and fully threaded screws in various sizes and lengths used for the fixation of the hand and foot. OsteoMed's lag screw design eliminates the need to over drill the proximal segment while the low profile head minimizes countersinking. M3-X™ features 1.2mm, 2.0mm, 2.4mm, and 3.0mm titanium lag screws and 1.2mm, 2.0mm, and 2.4mm fully threaded titanium screws.



**Low Profile Head**  
**Maximum Thread Area**  
**Self-Tapping**  
**Titanium Alloy**

### 1.2mm Fully Threaded Screws

204-1204	1.2mm x 4mm Screw
204-1206	1.2mm x 6mm Screw
204-1208	1.2mm x 8mm Screw
204-1210	1.2mm x 10mm Screw
204-1212	1.2mm x 12mm Screw

### 1.2mm Lag Screws

304-1206	1.2mm x 6mm Lag Screw
304-1208	1.2mm x 8mm Lag Screw
304-1210	1.2mm x 10mm Lag Screw
304-1212	1.2mm x 12mm Lag Screw
304-1214	1.2mm x 14mm Lag Screw
304-1216	1.2mm x 16mm Lag Screw

### 2.0mm Fully Threaded Screws

202-2008	2.0mm x 8mm Screw
202-2010	2.0mm x 8mm Screw
202-2012	2.0mm x 8mm Screw
202-2014	2.0mm x 8mm Screw
202-2016	2.0mm x 8mm Screw
202-2018	2.0mm x 8mm Screw
202-2020	2.0mm x 8mm Screw
202-2022	2.0mm x 8mm Screw

### 2.0mm Lag Screws

302-2008	2.0mm x 8mm Lag Screw
302-2010	2.0mm x 8mm Lag Screw
302-2012	2.0mm x 8mm Lag Screw
302-2014	2.0mm x 8mm Lag Screw
302-2016	2.0mm x 8mm Lag Screw
302-2018	2.0mm x 8mm Lag Screw
302-2020	2.0mm x 8mm Lag Screw
302-2022	2.0mm x 8mm Lag Screw

### 2.4mm Fully Threaded Screws

206-2408	2.4 x 8mm Screw
206-2410	2.4 x 10mm Screw
206-2412	2.4 x 12mm Screw
206-2414	2.4 x 14mm Screw
206-2416	2.4 x 16mm Screw
206-2418	2.4 x 18mm Screw
206-2420	2.4 x 20mm Screw
206-2422	2.4 x 22mm Screw

### 2.4mm Lag Screws

306-2408	2.4 x 8mm Lag Screw
306-2410	2.4 x 10mm Lag Screw
306-2412	2.4 x 12mm Lag Screw
306-2414	2.4 x 14mm Lag Screw
306-2416	2.4 x 16mm Lag Screw
306-2418	2.4 x 18mm Lag Screw
306-2420	2.4 x 20mm Lag Screw
306-2422	2.4 x 22mm Lag Screw
306-2424	2.4 x 24mm Lag Screw
306-2426	2.4 x 26mm Lag Screw
306-2428	2.4 x 28mm Lag Screw
306-2430	2.4 x 30mm Lag Screw
306-2432	2.4 x 32mm Lag Screw

### 3.0mm Lag Screws

308-2412	3.0 x 12mm Lag Screw
308-2414	3.0 x 14mm Lag Screw
308-2416	3.0 x 16mm Lag Screw
308-2418	3.0 x 18mm Lag Screw
308-2420	3.0 x 20mm Lag Screw
308-2422	3.0 x 22mm Lag Screw
308-2424	3.0 x 24mm Lag Screw
308-2426	3.0 x 26mm Lag Screw
308-2428	3.0 x 28mm Lag Screw
308-2430	3.0 x 30mm Lag Screw
308-2432	3.0 x 32mm Lag Screw
308-2434	3.0 x 34mm Lag Screw
308-2436	3.0 x 36mm Lag Screw

## OsteoMed Products



ExtremiLOCK Ankle Plating System



ExtremiLOCK Ankle Plating System



ExtremiLOCK Foot Plating System



ExtremiFix Small Cannulated Screws



ExtremiFix LDA Screws



Large Cannulated Screws



EnCompass



ReFlexion



Hemi



ExtremiFuse



InterPhlex



Talar-Fit



KobyGard



Inion



OsteoVation EX



OsteoVation QWIK



OsteoVation B-TCP



Allograft Offering