



OMITM
ORTHODONTIC ANCHOR SYSTEM



OMI™

Orthodontic Anchor System

Cruciform Driver

Pass Through Hole

Flanged Head Design



MTA™ Screws
(Maximum Thread Area)

AutoDrive™ Screw



- Patented Auto-Drive™ self-drilling screws eliminate the need for pre-drilling.
- 1.2mm and 1.6mm screws available in 6mm, 8mm, 10mm and 12mm lengths.
- Aperture and Flange design of the implant head provides for secure anchoring of dentition.

Part Number Description

Implants

202-1206	1.2mm x 6mm OMI™ Orthodontic Anchor Screw
202-1208	1.2mm x 8mm OMI™ Orthodontic Anchor Screw
202-1210	1.2mm x 10mm OMI™ Orthodontic Anchor Screw
202-1212	1.2mm x 12mm OMI™ Orthodontic Anchor Screw
202-1606	1.6mm x 6mm OMI™ Orthodontic Anchor Screw
202-1608	1.6mm x 8mm OMI™ Orthodontic Anchor Screw
202-1610	1.6mm x 10mm OMI™ Orthodontic Anchor Screw
202-1612	1.6mm x 12mm OMI™ Orthodontic Anchor Screw

Delivery System

202-0001	OMI™ Orthodontic Anchor System Organizer
----------	--

Implant Drivers

202-0003	OMI™ Orthodontic Anchor Finger Driver
220-0019	Taperlock™ Screwdriver Body
220-0061	1.6mm Screwdriver Shaft, Manual
220-0067	1.6mm Screwdriver Shaft, Dental Latch

Part Number Description

1.2mm Drills

220-0001	1.0mm Pilot Drill, J-Latch
220-0165	1.0mm Pilot Drill, 12mm Stop, Dental Latch

1.6mm Drills

220-0062	1.3mm Pilot Drill, J-Latch, 16mm
220-0064	1.3mm Pilot Drill, Dental Latch



3885 Arapaho Road
Addison, Texas 75001
Customer Service: 800.456.7779
Main: 972.677.4600
Fax: 800.390.2620
www.osteomedcorp.com