



smartflex
CRANIAL SPRING DISTRACTION SYSTEM



Rethinking Possibilities, Reshaping Lives

SPRING ASSISTED SURGERY FOR THE TREATMENT OF CRANIOSYNOSTOSIS

OsteoMed's SmartFlex offers early minimally invasive surgical intervention to decrease the morbidity associated with an extensive decompression operation.

SMARTFLEX CRANIAL SPRING ► Cranial Spring


- Pre-formed based on an extensive analysis of infant skulls
- Medical grade stainless steel
- Adjustable to accommodate various skull anatomies



Foot Plate

- Beveled to reduce the possibility of dural tear
- Designed to fit securely under the skull bone

TOP VIEW




Cranial Springs

Part Number	Description
218-3140-SP	SmartFlex SPRING NC - 4N
218-3145-SP	SmartFlex SPRING NC - 4.5N
218-3150-SP	SmartFlex SPRING NC - 5N
218-3155-SP	SmartFlex SPRING NC - 5.5N
218-3160-SP	SmartFlex SPRING NC - 6N
218-3165-SP	SmartFlex SPRING NC - 6.5N
218-3170-SP	SmartFlex SPRING NC - 7N
218-3175-SP	SmartFlex SPRING NC - 7.5N
218-3180-SP	SmartFlex SPRING NC - 8N
218-3185-SP	SmartFlex SPRING NC - 8.5N

Instruments & Accessories

Part Number	Description
220-0770	SmartFlex Sterilization Tray
220-0772	SmartFlex 5mm Rongeur
220-0773	SmartFlex Cutter



Rongeur
220-0772



Spring Cutter
220-0773

Rx ONLY
FEDERAL LAW (U.S.A.)
RESTRICTS THE DEVICE TO
SALE BY OR ON THE ORDER OF
A PHYSICIAN