




Acutrak® 3


Headless Compression Screw Solutions

Continuously Variable Fully

Threaded Pitch provides multiple points of fixation and a wide window of compression. The fracture does not have to be in a specific location due to fully threaded design



Headless design allows the screws to be implanted subchondral in and around articular regions with **minimal risk of impingement or soft tissue irritation**




Trusted 2 Million Times by Surgeons Worldwide

Acutrak 3

2+ Trusted Over
2 Million Times
Million

Acutrak Headless Compression Screw technology has been providing surgeons with operational efficiency and excellent patient outcomes through its ease of use and stable fixation. Since 1994, Acutrak has been trusted over two million times and continues to be refined based on surgeon feedback.

Acutrak 3 is the latest innovation and is composed of 93 unique screws, including a new 2.0 mm diameter Nano screw and extensions to the existing lengths of the Micro, Mini, and Standard screws, expanding on the many indications historically treated with the Acutrak 2 System.



Acutrak 3 Range of Sizes

Additional Sizes from Second Generations



4.0 Standard	3.5 Mini
2 mm increments 16-60 mm	2 mm increments 12-60 mm
2.5 Micro	2.0 Nano
1 mm increments 8-14 mm 2 mm increments 14-50 mm	1 mm increments 8-14 mm 2 mm increments 14-40 mm



2.0 Nano – Scaphoid Fracture



2.0 Nano – PIP Fusion



View More Acutrak 3 X-rays
go.acumed.net/AT3-Xrays

System Features



Parallel Wire Guide

Allows a guide wire to be placed parallel between 2.5 mm to 20 mm in distance

Screw Sizer

Provides ease of use to achieve unique angles for implant measurements and for difficult-to-reach areas due to the magnetic piston design



Contact Your Local Sales Rep for More Information



Learn More About Acutrak 3

Trusted 2 Million Times
by Surgeons Worldwide
go.acumed.net/AT3



www.acumed.net

Acumed USA Campus
5885 NE Cornelius Pass Road
Hillsboro, OR 97124
+1.888.627.9957

OsteoMed USA Campus
3885 Arapaho Road
Addison, TX 75001
+1.800.456.7779

Acumed Iberica Campus
C. de Álvaro Caballero, 14,
28023 Madrid, Spain
+34.913.51.63.57

SPF10-27-A | Effective: 2023/07 | © 2023 Acumed® LLC

OsteoMed® LLC is a wholly owned subsidiary of Acumed LLC.
OsteoMed is a registered trademark of OsteoMed LLC.

Acumed® is a registered trademark of Acumed LLC.