






Acutrak® 3

Headless Compression Screw Solutions

Acutrak 3 offers **solutions for all 28 bones** in the foot, featuring **93 unique screws** available in **4 diameters**



Headless design allows subchondral implantation of the screws in and around articular regions with **minimal risk of impingement or soft tissue irritation**




**A Trusted Solution for Lower Extremities
Acutrak 3 Has Every Step Covered**

Acutrak 3

3+ Trusted Over
3 Million Times
Million

The Acutrak Headless Compression Screw, introduced in 1994, was the first fully threaded screw on the market with a continuously variable threaded pitch. Acutrak 3 is the latest innovation, featuring 93 unique screws, including a new 2.0 mm diameter Nano screw and extended lengths for the Micro, Mini, and Standard screws. These advancements expand the system's ability to address many indications historically treated with the Acutrak 2 System. Indicated for fractures, fusions, and osteotomies, the Acutrak family of screws is a versatile system designed to treat conditions in the forefoot, midfoot, and hindfoot.



Acutrak 3 Range of Sizes

Additional Sizes from Second Generations



4.0 Standard	3.5 Mini
2 mm increments 16-60 mm	2 mm increments 12-60 mm
2.5 Micro	2.0 Nano
1 mm increments 8-14 mm 2 mm increments 14-50 mm	1 mm increments 8-14 mm 2 mm increments 14-40 mm




3.5 Mini– Hallux Valgus



4.0 Standard – Tarsometatarsal Joint Fusion Involving Naviculocuneiform and second and third TMT joints



View More Acutrak 3 X-rays
go.acumed.net/AT3-Xrays



Expanded Solutions for Small Bones

Upgraded Nano and Micro lengths offer a broader range of sizes, delivering enhanced options for small bone procedures like hammertoe

Hallux Valgus Fixation

With compression along its entire length, Acutrak 3 brings innovation and reliability to Hallux Valgus procedures

Contact Your Local Sales Rep for More Information



28 Operative Bones
in the Foot & Ankle
Acumed Can Treat Them All
Learn More: go.acumed.net/AT3



www.acumed.net



@acumed



@acumedglobal

United States | United Kingdom | Spain | Australia | New Zealand

SPF10-29-B | Effective: 2026/03 | © 2026 Acumed® LLC

Acumed® and Acutrak® are registered trademarks of Acumed LLC.